

!!! RNase-free working environment is required for the whole procedure!!!

Protein preparation

. Prepare 1 μ M (1pmol/ μ L, 161.1ng/ μ L) CAS9/CPF1 with transduction buffer

- a. Thaw the CAS9/CPF1 at r/t.
- b. Centrifuge the sample in a desktop microcentrifuge at maximum speed for 5min at r/t.
- c. Dilute CAS9/CPF1 10x with its corresponding 5xtransduction buffer.
- d. Check the concentration with Nanodrop-Protein280.
- e. Dilute to 161.1ng/ μ L (1 μ M) with its corresponding 1x transduction buffer.

Experimental procedure of the *in vitro* endonuclease activity assay

1) Assemble the reaction at room temperature, in the following order:

Nuclease free water	12 μ L
NEB buffer 3	2 μ L
1 μ M sgRNA	4 μ L
1 μ M CAS9 / CPF1	1 μ L
Reaction volume	19 μ L

Pre-incubate at room temperature for 20min

100fmol/ μ L substrate DNA fragment	1 μ L
Total reaction volume	20 μ L

- 2) Mix thoroughly and briefly spin down
- 3) Incubate in a thermocycler for 1hour at 37°C
- 4) Add 4 μ L of 6xDNA loading buffer
- 5) Checking on a 1.5% gel