

## 1-Transformation and colony inoculation

### Purpose:

Amplify target DNA sequence by transforming plasmid or ligation mixture into DH5 $\alpha$  competent cell.

### Materials:

- Plasmid or ligation product
- DH5 $\alpha$  competent cell (GeneMark, 100  $\mu$ l per tube)
- SOC solution or LB broth
- LB agar plate with antibiotic resistance such as Amp<sup>+</sup> or CamR<sup>+</sup>
- LB broth with antibiotic resistance

### Procedures:

[On ice]

- Thaw competent cell and pre-cool eppendorf tube.
- Inoculate 1  $\mu$ l DNA into 30-50  $\mu$ l competent cell and incubate on ice for 30 min. (\*The volume of DNA added should not exceed over 1% volume of competent cell.)
- Heat shock at 42°C for 30 sec, and incubate on ice for 2 min.
- Add 250  $\mu$ l SOC solution (or LB broth) and recover at 37°C with shaking 225rpm for 45 min.
- Spread all the mixture on the agar plate, and incubate at 37°C for 18~20 hr.

[The next day]

- Pick single colony and inoculate it into LB broth. Incubate at 37°C with shaking 225rpm for 18~20 hr.

## 2-Miniprep and Maxiprep

### Purpose:

Purify plasmid and digest pBKS and endo-1,4-glucanase with EcoRI and PstI restriction enzyme for vector and insert.

### Materials:

- Bacterial liquid
- 100% glycerol for stock
- Miniprep plus purification kit (BioTools)
- Maxiprep plus purification kit (GeneMark)

### Procedures:

[Miniprep]

- Transfer 0.5 ml bacterial cells and add 0.5 ml 100% glycerol (1:1) as a stock. The rest 1.5 ml is for Miniprep.

- Harvest bacterial cells in an eppendorf tube by centrifugation at 13,400 rpm for 1 min.
- Resuspend the cell pellet with **200 µl P1 buffer** (RNase A added and stored at 4°C)
- Add **200 µl P2 buffer** and mix thoroughly by inverting the tube 10 times. (clarified)
- Add **300 µl P3 buffer** and mix thoroughly by inverting the tube 10 times. (cloudy)
- Centrifuge at 14,000 rpm for 10 min.
- Transfer supernatant to the spin column and centrifuge at 13,400 rpm for 1 min.
- Add **400 µl PD buffer** and centrifuge at 13,400 rpm for 1 min. Discard the flow-through.
- Wash by adding **600 µl PW buffer** and centrifuge at 13,400 rpm for 1 min. Discard the flow-through.
- Centrifuge at top speed for 3 min. Discard the flow-through, open the lid and stay for a while.
- Place the spin column to a new eppendorf tube. Add **50 µl elution buffer**, stay for a while and centrifuge at top speed for 2 min.
- Measure the concentration of purified plasmid. (ng/µl)

[Maxiprep]

✧ Pick Single Colony and incubation for Maxiprep.

**Materials:**

- Plate with pUC57-Gal-Ste12 to Kpn1 in DH5α
- Plate with pUC57-Light fusion to suffix in DH5α
- Plate with pUC57-Gluxyn to suffix in DH5α
- 500 ml flask \*3

**Procedures:**

- **Pick single colony and stain it into LB broth with volume of 2-5 ml. Incubate at 37°C with shaking 225 rpm for 8 hr.**
- Then prepare 250 ml of LB broth (with AmpR) in the flask.
- Transfer small volume of bacterial liquid to large volume of LB broth.

[Note]: Antibiotic need 1000x dilution, so as 400ml LB broth need 0.4 ml antibiotic.

1. Pellet down cells from 250 ml liquid culture, centrifugation at 14000g for 5 min and pour off the supernatant.
2. Completely resuspend the pellet in 10 ml Solution I by pipetting.
3. Add 10 ml Solution II and mix gently. Incubate at RT for 4 min.
4. Add 10 ml Solution III and mix gently. Incubate on ice for 10 min.
5. Centrifuge the lysate at 14000g for 20 min at 4°C. A compact white pellet will be formed at the bottom of the tube.

6. To remove the floating pellet, hold a Filter Net by hand above a sterile 50 ml centrifuge tube and gently pour the lysate through the net.
7. Open the cap of Maxi-Spin Column provided in 50 ml centrifuge tube and add about 12 ml of cleared lysate into column. Close the cap and wait for 2 min for equilibrium with the membrane. Centrifuge the column at 12000 g for 1 min.
8. Discard the filtrate, and add the remaining lysate into the same Maxi-Spin column. Repeat centrifugation and discard the filtrate.
9. Add 10 ml Endotoxin Removal Wash Solution to the Maxi-Spin Column, wait for 2 min for equilibration with the membrane and centrifuge at 12000g for 1 min. Discard the filtrate.
10. Add 10ml Wash Solution to the column, and wait for 2 min for equilibration with the membrane and centrifuge at 12000g for 1 min and discard the filtrate. Repeat this step once more.
11. Reassemble the column to the tube and centrifuge at 12000 g for 10 min to dry the column.
12. For complete evaporation of ethanol, incubate the column at 60°C for 15 min.
13. To elute the DNA, place the Maxi-Spin Column into a sterile 50 ml centrifuge tube and add 2 ml preheated Elution Solution to the membrane. Let stand for 3 min and centrifuge at 12000g for 5 min to elute DNA. Withdraw eluent from 50 ml centrifuge tube and add to the Maxi-Spin Column, repeat centrifuge from this step to increase DNA yield.
14. Collect the DNA solution into four 1.5 ml Eppendorf and store at -30°C.

### 3-Restriction enzyme analysis and digestion

#### Purpose:

Digest DNA with restriction enzyme to analyze the plasmid mapping or cut the restriction site to make vector and insert.

#### Materials:

- DNA plasmid
- NEB restriction enzyme
- NEB restriction enzyme buffer
- ddH<sub>2</sub>O
- Gel electrophoresis

#### Procedures:

	Vector	Insert	Plasmid mapping
<b>DNA</b>	2000 ng	2000 ng	<b>1000 ng</b>
<b>10X buffer (cutsmart)</b>	5	5	<b>2</b>
<b>EcoRI-HF</b>	0.5	0.5	<b>0.2</b>
<b>PstI-HF</b>	0.5	0.5	<b>0.2</b>

ddH <sub>2</sub> O	Up to total	Up to total	Up to total
Total volume	50	50	20

- Incubate at 37 °C water bath for 1-2 hr.

## 4.1-Gel extraction

### Purpose:

Extract DNA fragments from enzymatic reaction or PCR amplification.

### Materials:

- Electrophoresis with 1% agarose gel
- DNA extraction kit (GeneMark)

### Procedures:

- After electrophoresis, cut out the target DNA band from agarose gel.
- Transfer gel slice into an eppendorf tube and add 2 volumes of binding buffer. Incubate at 60°C until the gel dissolved.
- Transfer the solution into the spin column, and centrifuge at 15,000 g for 1 min. Discard the filtrate.
- Add 500 µl of Binding solution and centrifuge at top speed for 1 min. Discard the filtrate.
- Add 700 µl of Wash solution and wait for 1 min. Centrifuge at 15,000 g for 1 min and discard the filtrate. Repeat one more time.
- Centrifuge at top speed for additional 5 min to remove residual ethanol.
- Transfer the spin column into a new Eppendorf tube, add 30 µl of Elution buffer and wait for 2 min.
- Centrifuge at top speed for 2 min to elute DNA.

## 4.2-PCR/DNA clean-up

### Purpose:

Extract DNA fragments from enzymatic reaction or PCR amplification without running electrophoresis on agarose gel.

### Materials:

- DNA extraction kit (GeneMark)

### Procedures:

- Transfer PCR product or other enzymatic reaction mixture into an eppendorf tube and add 3 volumes of binding buffer. Vortex briefly to mix.

- Transfer the solution into the spin column, and centrifuge at 15,000 g for 1 min. Discard the filtrate.
- Add 700 µl of Wash solution to the spin column and wait for 1 min. Centrifuge at 15,000 g for 1 min and discard the filtrate. Repeat one more time.
- Centrifuge at top speed for additional 5 min to remove residual ethanol.
- Transfer the spin column into a new Eppendorf tube, add 30 µl of Elution buffer and wait for 2 min.
- Centrifuge at top speed for 2 min to elute DNA.

### 5-T4 DNA Ligation

**Purpose:**

Ligate target insert into vector which both digested with the same cutting site and make a recombinant plasmid.

**Materials:**

- Digest vector and target insert
- T4 ligase and 10X T4 ligation buffer (NEB)
- ddH<sub>2</sub>O
- Gel electrophoresis equipment for ligation check

**Procedures:**

[Verification of digested vector and insert before ligation]

- Load: marker 5, uncut 5+1, digest vector 5+2, uncut5+1, digest insert5+2
- According to the brightness of the band on the gel, estimate the ratio of concentration between vector and insert. [Calculation]-vector: insert=the ratio of brightness/ basepairs = estimated ratio of concentration.
- Insert : vector= 3 : 1

	µl
Vector	1.0
Insert	1.5
T4 ligase	1.0
10X ligation buffer	2.0
ddH <sub>2</sub> O	14.5
<b>20.0</b>	

Calculation: according to the band brightness of digested vector and insert, 3/1 multiply 1441/2961 = 1.46/1 (insert: vector)

- Incubate at RT for 30-60 min. (More incubation time (2 hr) for blunt end)

[Note]:

**Ligation information (insert → vector)**

1. Light fusion to suffix (4323 bp, 529ng/μl) →pRS423 (HIS3) (change the plasmid to pRS112)  
→Restriction enzyme: Kpn1 & Pst1
2. Gal-Ste12 to Kpn1 (3082 bp, 302 ng/μl) →YEplac195 (URA3. 5233+8 bp)  
→Restriction enzyme: EcoR1 & Kpn1
3. Gluxyn to suffix (2887 bp, 50 ng/μl-condense to 80 ng/μl) →YEplac181 (LEU2. 5214+27 bp)  
→Restriction enzyme: Kpn1 & Pst1