

# Flanking Region PCR Reaction

## Introduction

Previously cut flanking regions will be added to the AddGene HTT sequence to complete exon 1.

## Materials

- › Amplified AddGene HTT sequence
- › Cut flanking regions
- › Primers for flanking regions
- › Polymerase

## Procedure

### Set-up PCR

1. Mix equal molar amounts of flanking region DNA and AddGene HTT DNA at concentrations (1-500 ng) sufficient for initial PCR.
2. Add primers at standard PCR concentrations (0.2-0.5  $\mu$ M).
3. Initially denature dna at 95°C for 10 min (**with no polymerase**)
4. Add Taq Polymerase.
5. Cycle (20x) at:

PCR Amplification Parameters			
	A	B	C
1	Tempature	Time	Cycle name
2	95°C	30-60 s	Melt
3	T <sub>m</sub> - 5°C ( <a href="https://tmcalculator.neb.com/">https://tmcalculator.neb.com/</a> )	30-60 s	Anneal
4	72°C	1 min / kb	Elongate

6. End cycle as follows

End Cycle Parameters			
	A	B	C
1	72°C	10 min	Final Elongation
2	4°C	Indefinitely	Hold Chilled

7. Gel Verify the amplification

8. Store the DNA at 4°C