Conjugation of antigen to gold nanoparticles

- A. Preparation :
 - Measure the concentration of antigen at O.D. 280 nm
 - Know the quantity of the final product
 - Know the pH of the final product (adjusted by K₂CO₃)
- B. Theory:

Gold nanoparticles tend to conjugate on antigen due to the van der Waals force

- C. Material & Method :
 - Prepare 400 µl 5X gold nanoparticle buffer, diluted by 2mM borax buffer and filtered by 0.22 µm filter
 - add _____ µg antibody (diluted <u>10X</u> in 0.01M PBS), add into the solution in step1
 - RT stir <u>5~6</u> rpm [,] <u>60</u> min (particles are shown in <u>bright red)</u>
 - Add 10%BSA to 0.01M PBS [,] filter through 0.45µm filter [,] add solution <u>220</u>µL to the gold solution [,] the final concentration will be 1% BSA. RT stir for 30 minutes
 - Distribute the solution into <u>2</u> tubes , centrifuge at 4°C , 13000rpm for 30 minutes
 - Remove the supernatant [,] restore the pellet by buffer A, <u>90</u> µL each
 - Store the conjugated antigen in a 4ml vial, 4°Cfreezer
 - Final solution are shown in <u>burgundy red</u>, results: <u>I is larger than 20ppb~50ppb</u>