



TARDIGARD
p r o t o c o l

iGEM Team SIAT—SCIE
Email: siatxscie@gmail.com

Colony Picking and Colony PCR

Picking clones for further processing, may be colony PCR to verify if the plasmid is successfully transferred, or to further culture bacteria in shake incubator.

Colony PCR is to verify whether the desired plasmid is successfully transferred

Colony picking and colony PCR

1. Colony picking

1.1 Materials

Tips

Pipettes

Alcohol lamp

Tweezers

Shake incubator or pre-prepared PCR mix (depending on your next experiment plan)

1.2 Set up and protocol

- a. Put the tip of the tweezer into alcohol lamp's flame for a few seconds
- b. Using the tweezer to pick up half/a small portion of the colony
- c. Mix it into PCR mix you have prepared or shake incubator with antibiotics and liquid medium

2. Colony PCR

2.1 Material

PCR machine

1.5ul centrifuge tube

Autoclaved ddH₂O

Ice in foam box

Pipettes and tips

Primer designed earlier

2.2 Setup and protocol

We are using product from takara, detail of product and protocol can be found in this linked website:

http://r.search.yahoo.com/_ylt=A0SO80_TwKxZr3IARTVXNyoA;_ylu=X3oDMTEzOGdwa2xuBGNvbG8DZ3ExBHBvcwMyBHZ0aWQDVUkwMkMOXzEEc2VjA3Ny/RV=2/RE=1504522579/RO=10/RU=http%3a%2f%2fwww.clontech.com%2fxxclt_ibcGetAttachment.jsp%3fcltemId%3d14060%26minisite%3d10020%26secltmId%3d16633/RK=1/RS=P7D1q7uRnp5FwIPxVQAa_X7VeQE-

a. First prepare an 200 ul system in 1.5ul centrifuge tube

b. Prepare pcr mix

template 1ul

primer 1ul

primer 1ul

prime star premix= 100ul

ddh₂o 90ul

add in order of largest volume to smallest

c. Set up the programme

The programme we are using is

| | | |
|-----|-----|-------|
| | 98C | 3min |
| X30 | 98C | 10sec |
| X30 | 55C | 15sec |
| X30 | 72C | 10sec |
| | 72C | 10min |
| | 12C | -- |

d. Run gel electrophoresis

Analyse the result obtained, if target length of sequence is obtained