

VLP Generation

Sf21 Transfection

For the composition of Grace's medium, refer to: <https://www.thermofisher.com/nl/en/home/technical-resources/media-formulation.71.html> (retrieved: 29-09-2017)

1. Seed the Sf21 cells in a 6-wells plate. Use around 8×10^5 cells/well in 2mL Grace's medium (no antibiotics).
2. Give the cells 2 to 3 hours to attach to the well, or incubate overnight.
3. Combine 10 μ L Bacmid DNA and 90 μ L Grace's medium and set aside. Avoid pipetting up and down!
4. Prepare the transfection mixture by combining 5 μ L Express2 TR and 95 μ L Grace's medium. Add this to the 100 μ L DNA mixture and incubate for 20 minutes at room temperature.
5. Add the full Transfection:DNA mix to the medium of the cells, in a dropwise manner.
6. After 4 hours, replace the culture medium with 2mL fresh Grace's medium+antibiotics. Incubate at [...]°C.
7. In 2-3 days, monitor the cells for presence of baculovirus signs.

Bacmid DNA Isolation

1. Pellet the cells in a table top centrifuge at 2000x g for 15 min, at 4°C.
2. Resuspend the cell pellet in 200 μ L chilled GTE and transfer it to an Eppendorf tube.
3. Add 400 μ L NaOH solution and mix by inverting the tube a few times. Incubate for 5 minutes at RT afterwards.
4. Add 300 μ L KAc solution, invert the tube several times, avoid vortexing! Incubate on ice for 5 minutes.
5. Centrifuge 10 min at max speed in the Eppendorf centrifuge at 4°C.
6. Transfer 1ml supernatant to a new Eppendorf tube and add 0.6x volume (600 μ L) ice-cold 100% isopropanol and mix by inverting the tube.
7. Centrifuge for 15 min at max speed, 4°C. Remove the supernatant afterwards.
8. Wash the pellet with 500 μ L 70% ethanol and centrifuge again for 5 min, max speed at 4°C.
9. Remove the supernatant. If there is still a lot of liquid left, spin an additional 5 min to remove all liquid.
10. Air dry the pellet for 15-30 minutes.
11. Add 40 μ L dH₂O to the pellet and leave at RT for 15-30 minutes.
12. Dissolve the pellet by flicking the tube (!do not pipette up and down!) and store the bacmid DNA at 4°C.

RNA Isolation

Sample Lysis

1. Pellet the cells by centrifugation and discard the supernatant.
2. Add TRIzol™ Reagent to the sample in a final 3:1 ratio (TRIzol:sample). Homogenize the mixture by pipetting up and down.
3. Incubate for 5 minutes at room temperature and add 0.2mL chloroform per mL of TRIzol™ Reagent added previously.
4. Incubate for 3 more minutes, then centrifuge the sample for 15 minutes at 12000x g, 4°C.
5. Transfer the upper aqueous phase to a new tube.

RNA Isolation

6. Add 0.5mL isopropanol per mL of TRIzol™ Reagent used for sample lysis and incubate for 10 minutes.
7. Centrifuge for 10 minutes at 12000x g at 4°C and discard the supernatant.
8. Resuspend the pellet in 1mL of 75% ethanol per mL of TRIzol™ Reagent used for lysis.
9. Vortex briefly, then centrifuge for 5 minutes at 7500x g at 4°C. Discard the supernatant.
10. Air dry the pellet for 10 minutes.
11. Resuspend the pellet in 20-50 μ L MQ. Homogenize by pipetting up and down.
12. Use RNA for further applications or freeze in -20°C.



Semi-Dry Blotting

1. Soak per gel 2x 3 Whatmann filter papers in 1x Semi-dry blotting buffer.
2. Place the Immobulon-P membrane in a 96% ethanol solution.
3. Assemble the following construct onto the blotting machine, using this order:

3x Whatmann paper>Immobulon-P membrane>SDS-gel>3 Whatmann filter papers

4. Attach the lids of the machine and blot with a constant 0.05A/gel for 1 hour.
5. Remove the blot from the machine afterwards and place it into a 50mL tube, with 10mL 1% ELK added.
6. Incubate at room temperature for 1 hour, or overnight at 4°C.
7. Discard the ELK and add 5mL of primary antibody in 1% ELK using the correct dilution.
8. Incubate at room temperature for 1 hour, or overnight at 4°C.
9. Discard the solution and wash 3x with 10mL PBS-T for 5 minutes.
10. Add 5mL of secondary antibody, diluted in PBS-T.
11. Incubate at room temperature for 1 hour, or overnight at 4°C.
12. Discard this solution as well and wash with 3x 10mL of PBS-T for 5 minutes.
13. Add 10mL AP-buffer and incubate for 10 minutes at room temperature.
14. Add 5mL NBT/BCIP solution (made by adding 75µl to 5mL AP-buffer).
15. Develop the membrane until bands appear, then wash with 3x water to stop the reaction.

Baculovirus Titration by End Point Dilution Assay (EPDA)

1. Prepare a dilution series of the baculovirus suspension from 10^{-1} to 10^{-9} in culture medium.
2. The 10^{th} sample will be a negative control, only containing 90µl Sf9-ET medium.
3. Dilute the Sf cells to a concentration of $1.5 \cdot 10^6$ cells/mL.
4. Add 90µl Sf9-ET cell suspension to each tube containing virus dilution and the negative control. Mix well.
5. Add 10µl of each tube to a 6 wells of a 60-wells microtiter plate. Incubate 4-7 days at 27°C in a humid box.
6. Count the number of infected wells and calculate the TCID₅₀/mL (Tissue Culture Infected Dose₅₀) using the following formula:

$$\text{TCID}_{50} = 200 \cdot 10^{(a+b)}$$

In which:

a = the highest dilution of which the percentage of AIW is higher than 50%

b = $(c-50\%)/(c-d)$

c = percentage of AIW of dilution a.

d = percentage of AIW of a ten times dilution of dilution a.

Transmission Electron Microscopy (TEM) sample preparation

1. Prepare tungsten TEM grids in a vacuum evaporator.
2. Apply 5µL of sample and incubate for 1 minute. Remove excess liquid with blotting paper.
3. Apply 5µL of phosphotungstic acid staining and incubate for 30 seconds RT. Remove excess liquid here as well.
4. Use sample for microscopy.

