

## Plasmid Assembly Protocol with Golden Gate Assembly

- ◆ Select a 20 nt protospacer of interest. The 3' protospacer adjacent sequence (PAM) must be NGG, where N is any nucleotide. Preference is given to:
  - ◇ Sequences with purines occupying the last four (3') bases of the protospacer.
  - ◇ Sequences on the non-coding strand.
  - ◇ Sequences in which the last 12 nt of protospacer + 3 nt PAM (15 nt total) are unique in the genome (check by BLAST with all four possible NGG sequences).
- ◆ Design two 24 nt oligonucleotides (4 nt 5' sticky end + 20 nt spacer sequence) with the sticky ends ACGC on the forward primer and AAAC on the reverse primer.
- ◆ For single spacers, anneal spacer oligos as follows
  - ◇ Resuspend both oligos to 100 $\mu$ M in water
  - ◇ Mix 5  $\mu$ L FOR + 5  $\mu$ L REV + 90  $\mu$ L 30 mM HEPES, pH 7.8
  - ◇ Heat to 95  $^{\circ}$ C for 5 min, then ramp to 4  $^{\circ}$ C at 0.1  $^{\circ}$ C/sec
- ◆ Insert annealed spacer (or dual+spacer synthetic construct) by Golden Gate assembly.
  - ◇ X  $\mu$ L Backbone (100 ng)
  - ◇ 0.3  $\mu$ L Insert
  - ◇ 2  $\mu$ L T4 Ligase Buffer (NEB)
  - ◇ 1  $\mu$ L T4 ligase (NEB)
  - ◇ 1  $\mu$ L BbsI (NEB)
  - ◇ Fill up to 20  $\mu$ L with H<sub>2</sub>O

- ♦ Golden Gate Program:

	°C	min
9x	37.0	10
	16.0	10
	50.0	5
	65.0	20
	4.0	PAUSE

- ♦ Transform 3  $\mu\text{L}$  of each reaction to E. coli NEB5alpha by heat shock
- ♦ Plate 10% of recovery culture on selective plates with 10  $\mu\text{L}$  of 0.5 M IPTG and 40  $\mu\text{L}$  of 20  $\text{mg mL}^{-1}$  Bluo-gal (in DMSO).
- ♦ Pick white colonies to selective LB and recover plasmid

**Cobb, R.E., Wang, Y., and Zhao, H.** (2014). High-Efficiency Multiplex Genome Editing of *Streptomyces* Species Using an Engineered CRISPR/Cas System.: 1–12.