

Blunt ligation

To finish vector creation, we need to close the ends.

Protocol

1. Amplify the vector
2. Add 1 μL DpnI restriction enzyme directly to 50 μL PCR product.
3. Incubate 37 °C 2 hours
4. Make a gel purification in 50 μL (using the gel purification protocol)
5. Phosphorylate with 0.5 μL PNK4 and 1 μL T4 ligation buffer.
6. Incubate 37 °C 30 minutes.
7. inactivate by placing in 65 °C 20 minutes
8. Add 0.5 μL T4 ligase
9. Store over night in 4 °C or 2 hours in 37 °C.
10. Transform

Component	10 μL reaction
10x NEB buffer	1.0 μL
Vector	4.0 μL
Restriction enzyme 1	0.25 μL
Restriction enzyme 2	0.25 μL
Nuclease-free water	4.5 μL

Incubation

Incubate in PCR machine

1. Incubate: 37 °C 30 minutes

Now the vector is ready to insert genes in, by using the USER ligation protocol, or for PCR amplification.