

iGEM Team SIAT—SCIE Email: siatxscie@gmail.com

### **Protocol for Radiation Resistance Test**

In order to evaluate the effect of the protein Dsup on *E. coli*, we designed this protocol. Dsup is a self-expressed protein that can protect the bacteria from radiation by binding with DNA<sup>[3]</sup>. It is expected that the LD<sub>50</sub> of UV can be obtained by setting up a dose gradient.

### **Equipment and materials:**

Ultraviolet source, microscope, hemocytometer<sup>[1]</sup>, trypan blue stain (0.8mM in polybutylene succinate, PBS)<sup>[1]</sup>.

# 1. Radiation disposal<sup>[4]</sup>

- We expect to use a wide spectrum of UV radiations with a peak of emission at 312nm (UV-B), which is considered the most cytotoxic and mutagenic wavelength among types of solar radiation<sup>[5]</sup>.
- Set up a gradient of UV doses, e.g. 10kJm<sup>-2</sup> 100kJm<sup>-2</sup>.
- Apply newly cultured bacteria in LB broth respectively on small culture dishes.
- Put the dishes under the UV lamp without cover. The UV lamp was positioned 30 cm above the samples. Expose the sample for a period of 180 min.
- The temperature should be kept constant during the whole process of irradiation, and the surroundings should be kept dark.
- For un-irradiated control, another dish is placed in the same temperature without direct UV radiation.

## 2. Survival rate assay<sup>[1]</sup>

- Mix irradiated cells 1:1 with trypan blue solution (for no more than 30 min because of its toxicity).
- Count the number of blue (dead) and unstained (living) cells respectively. For usage of hemocytometer, see <a href="https://en.wikipedia.org/wiki/Hemocytometer">https://en.wikipedia.org/wiki/Hemocytometer</a>.
- Calculate the proportion of two types of cells.

#### **References:**

- 1. Sandy Westerheide, Northwestern University, Protocol for use of tripan blue dye.
- 2. Hemocytometer Wikipedia -.
- 3. Takuma Hashimoto et al. Extremotolerant tardigrade genome and improved radiotolerance of human cultured cells by tardigrade-unique protein.
- 4. Tiziana Altiero et al. Ultraviolet radiation tolerance in hydrated and desiccated eutardigrades.
- 5. Connelly SJ. et al. Temperature effects on survival and DNA repair in four freshwater cladoceran Daphnia species exposed to UV radiation.