

6/05/17

Monday, June 5, 2017 9:05 AM

Who was in lab today: Jeff, Ayesha, Qingxi, Ana, Salma, Martin

Picked some plasmids to use in future Gibson Assembly to ligate our genes into

- pSB3T5 (TetR, p15A ori) http://beta.labgeni.us/registries/parts_registry/?part=pSB1C3
- pSBK3 (KanR, pMB1 ori) http://beta.labgeni.us/registries/parts_registry/?part=pSB1C3

Plasmids picked have different antibiotic resistances and origins of replications so that the E. Coli will not be forced to uptake/choose only one

Who was in lab today: Jeff, Ayesha, Qingxi, Ana, Salma, Martin

Ligation of Insert and Vector:

20uL total

12uL PET28

5uL cj blue

2uL buffer

1uL DNA Ligase

****ADD THE LIGASE LAST!!!**

Making the agarose gel for gel electrophoresis:

1 % agarose gel (essentially 1g/ 100 mL) is the base range for separating DNA.

0.8% agarose is used for separating plasmids.

0.6% -2% agarose gel is used for separating DNA base pairs.

<0.6% too floppy

We used 500mg agarose in 50mL TAE buffer to create 1% agarose buffer.

Protocol:

1. Microwave for 130 seconds at a time.
2. Add 5 uL of Ethidium Bromide to heated mixture.
3. Add 10 uL of 6X loading dye to DNA digest samples.
4. Run Gel at 100mA.
 - Rule of Thumb: 10 mA per inch of gel you run.

Monday, June 5, 2017

1:17 PM

Gel Extraction/Purification Protocol

1. Excise the DNA fragment from the agarose gel with a clean, sharp razor blade.
2. Weigh the gel slice in a colorless tube. Add 3 volumes Solubilizing Buffer to 1 volume gel (100 mg gel ~ 100 μ l).
3. Incubate on 50°C water bath for 10 minutes until gel is completely dissolved in buffer (option to vortex every 2 minutes to mix)
4. Pipet 800 μ L of the gel/buffer mixture into a spin column, centrifuge for 1 min @ 13000 rcf, discard the flow through from the collection tube, repeat until all of the liquid has been spun down
5. Pipet 500 μ L solubilizing buffer to the white membrane of the spin column
6. Add 720 PE μ L to the white membrane of the spin column, let sit for 5 minutes
7. Centrifuge for 1 min @ 13000 rfu, discard flow through
8. Centrifuge for 2 min @ 13000 rfu, discard flow through
9. Add 50 μ L EB buffer to spin column in clean eppendorf tube, then spin down