

Cloning and DNA Manipulation

Restriction Enzyme Digestion

Source: <https://www.neb.com/protocols/2012/12/07/optimizing-restriction-endonuclease-reactions> (Retrieved: 28-08-2017)

1. Look up the enzyme properties and set up the following reaction:

Component:	Volume:
Restriction Enzyme	1 μ L
DNA	... μ L (1 μ g)
NEBuffer: CutSmart/1.1/2.1/3.1	5 μ l (1X)
MQ	... μ L
(CIP/rSAP/other)*	(2x 1 μ L)
Total Reaction Volume	50 μl
Incubation Time	1 hour (add 1 μ L CIP after 30 minutes)
Incubation Temperature	37°C (enzyme-dependent)

2a. Run 5 μ L on gel to check for correct fragment sizes and heat inactivate the enzymes (65°C/80°C for 20 min, depending on the enzyme). PCR purify the samples afterwards.

2b. Run everything on gel and gel purify the band with the correct size.

* = Phosphatases (CIP/rSAP/other) may be added to prevent back-ligation of any of the parts, if applicable. Add 1 μ L to the initial mix and add an additional 1 μ L of CIP/rSAP after 30 minutes of incubation.

Retrieving DNA from iGEM plates

(Adapted from the iGEM website)

1. With a pipette tip, puncture the foil of the desired well.
2. Add 10 μ L MQ and pipette up and down a few times. Incubate for 5 minutes to make sure the DNA is fully resuspended.
3. Transform 1 μ L of this DNA into your desired strain by a transformation method of your choice and plate the culture on plates with the appropriate antibiotic. Incubate O/N at 37 °C.
4. Pick a single colony and inoculate in LB with the appropriate antibiotic O/N at 37°C, 200rpm.
5. Make glycerol stocks of the bacteria and miniprep the rest of the culture for use.

T4 Ligation

Source: <https://www.neb.com/protocols/1/01/01/dna-ligation-with-t4-dna-ligase-m0202> (Retrieved: 28-08-2017)

1. Set up the following reaction (add DNA ligase last, vector:insert ratio of 3:1)*:

Component:	Volume:
T4 DNA Ligase Buffer (10x)	2 μ L
Vector DNA (Backbone)	... μ L (50ng)
Insert DNA (Construct)	... μ L (37.5ng)
MQ	... μ L
T4 DNA Ligase	1 μ L
Total:	20 μL

2. Mix the solution by pipetting up and down. Centrifuge briefly to collect the liquid.

3a. For sticky ends, incubate at RT for 10 minutes/For blunt ends, incubate at RT for 2 hours.

3b. Alternatively, ligate O/N at 16°C for tricky ligations.

3c. Alternatively, ligate O/N in an endless cycle of 30 seconds at 30°C and 30 seconds at 10°C [1].

4. Heat inactivate the ligase at 65°C for 10 minutes.

5. Chill on ice and transform 2 μ L of the reaction into 50 μ L competent cells/Store at -20°C until transformation.

* = Other vector:insert ratios may be preferable, use the NEBioCalculator (<https://nebiocalculator.neb.com/#!/ligation>) to estimate ratios. Multiple inserts can be used in one reaction, subtract the needed volumes from MQ in that case.



cDNA Synthesis using Reverse Transcriptase

Source: <https://www.thermofisher.com/order/catalog/product/18080044> (Retrieved: 28-09-2017)

N.B. A 20 μ L reaction volume can be used for 1ng-5 μ g of total RNA or 1-500ng of mRNA. For further applications, it may be wise to remove the complementary RNA by adding 1 μ L (containing 2 units) of *E. coli* RNase H and incubate at 37°C for 20 minutes.

1. Set up the following reaction:

Component:	Volume:
50-250ng random hexameric primers <i>or</i> 2pmole gene-specific primers	1 μ L
1ng-5 μ g total RNA <i>or</i> 1-500ng of mRNA	... μ L
dNTP mix (10mM)	1 μ L
MQ	Fill to 12 μ L
Total:	12μL

2. Heat the mix to 65°C for 5 min and quickly chill on ice.

3. Add the following components:

Component:	Volume:
5x First-Strand Buffer	4 μ L
0.1M DTT	2 μ L
RNase Inhibitors (40 units/ μ L)	1 μ L
Total (+earlier volume):	19μL

4. Mix the tube contents gently. Incubate at 42°C for 2 minutes.

5. Add 1 μ L (containing 200 units) of Reverse Transcriptase and mix by pipetting up and down.

6. Incubate the mix at 25°C for 10 min.

7. Increase the temperature to 42°C and incubate for 50 more minutes.

8. Heat inactivate the enzyme by heating the mix to 70°C for 15 minutes.

Golden Gate Assembly

N.B. If the T4 Ligase Buffer negatively affects the chosen enzyme, it can be replaced by any NEBuffer (1.1/2.1/3.1 or CutSmart). Take in mind that 1.5 μ L of a 10mM ATP solution should be added to the mix as well in order for the ligase to work.

1. Set up the following reaction (vector:insert ratio of 1:2, use 50-75ng of vector DNA):

Component:	Volume:
Vector DNA	... μ L
Insert DNA*	... μ L
DpnI	0.35 μ L
Type IIS Enzyme	1 μ L
T4 DNA Ligase	1 μ L
T4 DNA Ligase Buffer	1.5 μ L
BSA (1mg/mL)	1.5 μ L
MQ	... μ L (fill to 15 μ L)
Total:	15 μL

Step:	Temperature (°C):	Time (minutes):
Incubation	37	20
25x	37	3
	16	4
Heat inactivation	80	20
Cooling	16	1
Hold	12	∞

2. Incubate the mix according to the table on the right.

3. Use 1-2 μ L of this mixture for transformation, or store at -20°C until further use.

* = Other vector:insert ratios may be preferable, use the NEBioCalculator (<https://nebiocalculator.neb.com/#!/ligation>) to estimate ratios. Multiple inserts can be used in one reaction, subtract the needed volumes from MQ in that case. To estimate background colonies for transformation, make an identical mix, but leave the insert DNA out.





Gateway cloning

N.B. Generally, try to use the 4 μ L reaction volume. However, if the relevant PCR product concentrations are low(ish), resort to the 8 μ L end volume. To reach the maximum number of recombinants, extend the first incubation time to up to 18 hours.

1. Set up the BP mix as described in the table below:

Component:	1x Volume:	2x Volume:
PCR Product	... (150ng)	... (150ng)
pDONR207	... (150ng)	... (150ng)
TE Buffer	Add up to 4 μ L	Add up to 8 μ L
BP Clonase	1 μ L	2 μ L
Total:	4μL	8μL

2. Incubate for 1 hour at room temperature.

3. Add 1 μ L Protein Kinase to inactivate enzymes and incubate at 37°C for 10 min.

4. Electroporate 1-2 μ L into DH10 β and plate on LB+genta agar plates. Incubate overnight at 37°C.

References

- [1] A. H. Lund, M. Duch, and F. Skou Pedersen, "Increased cloning efficiency by temperature-cycle ligation," *Nucleic Acids Res.*, vol. 24, no. 4, 1996.

