

For the amplification of the promoter and the gene we transformed them into Escherichia coli JM-109 and Dh5-Alpha bacteria.

|          | Sample ID | Method          | Competent cell type | Volume | Annotation |
|----------|-----------|-----------------|---------------------|--------|------------|
| Promoter | P1        | thermal shock   | JM 109              | 150 µL |            |
|          | P2        | thermal shock   | JM 109              | 150 µL |            |
|          | P3        | thermal shock   | JM 109              | 150 µL | no colony  |
|          | P4        | electroporation | DH5-Alpha           | 200 µL |            |
|          | P5        | electroporation | DH5-Alpha           | 200 µL |            |
|          | P6        | electroporation | DH5-Alpha           | 200 µL | no colony  |
| Gene     | G1        | thermal shock   | JM 109              | 150 µL |            |
|          | G2        | thermal shock   | JM 109              | 150 µL |            |
|          | G3        | thermal shock   | JM 109              | 150 µL |            |
|          | G4        | electroporation | DH5-Alpha           | 200 µL |            |
|          | G5        | electroporation | DH5-Alpha           | 200 µL |            |
|          | G6        | electroporation | DH5-Alpha           | 200 µL |            |

We used both heat shock and electroporation in order to get safe results.

During thermal shock the reaction mixes were incubated at 37 °C for 10 min. After the incubation, we transformed E.coli JM109 competent cells with the samples of table (number) using thermal shock by the following protocol.

1. On ice, to 50 µL of E. coli JM109 competent cell solution 5 µL of the reaction mix was confused in an Eppendorf-tube, then the solution was left on ice for 30 minutes. During that, SOC-containing soil was placed on ice as well.
2. After 30 minutes, thermal shock was accomplished at 42 °C for 60 seconds.
3. When the thermal shock was ready, the Eppendorf-tubes were placed back on ice, and right after it, 200 µL of the SOC-containing soil was pipetted into each Eppendorf-tube.

The mixtures were incubated at 37 °C for 1 hour.

In the case of electroporation we used 0.1 cm tubes, and we applied 1.8 kV.

In 10 of 12 Petri-dishes containing LB and Ampicillin, the colonies have successfully developed and we were able to select individual ones and grow them in test tubes. We used a growth medium containing LB and Ampicillin as well, therefore only successfully transformed bacteria with were able to produce.

For the amplification of plasmids pMHE5, pMHE7 and pDSK5 we cultivated the bacteria containing them in mediums including antibiotic. In the case if bacteria containing pDSK5 they were grown in LB medium with Kanamycin and bacteria with pMHE5 and pMHE7 were cultivated in LB medium containing Streptomycin. With this method, only the certain plasmid-containing bacteria were able to grow, preventing being infected by other ones. From every Petri-dishes 16 colonies were separated and put in individual test tubes.

| Petri-dishes | Test-tubes | Plasmid | Antibiotic   |
|--------------|------------|---------|--------------|
| 1.           | pM5 01     | pMHE5   | Streptomycin |
|              | pM5 02     | pMHE5   | Streptomycin |
|              | pM5 03     | pMHE5   | Streptomycin |
|              | pM5 04     | pMHE5   | Streptomycin |
|              | pM5 05     | pMHE5   | Streptomycin |
|              | pM5 06     | pMHE5   | Streptomycin |
|              | pM5 07     | pMHE5   | Streptomycin |
|              | pM5 08     | pMHE5   | Streptomycin |
|              | pM5 09     | pMHE5   | Streptomycin |
|              | pM5 10     | pMHE5   | Streptomycin |
|              | pM5 11     | pMHE5   | Streptomycin |
|              | pM5 12     | pMHE5   | Streptomycin |
|              | pM5 13     | pMHE5   | Streptomycin |
|              | pM5 14     | pMHE5   | Streptomycin |
|              | pM5 15     | pMHE5   | Streptomycin |
|              | pM5 16     | pMHE5   | Streptomycin |
| 2.           | pM7 01     | pMHE7   | Streptomycin |
|              | pM7 02     | pMHE7   | Streptomycin |
|              | pM7 03     | pMHE7   | Streptomycin |
|              | pM7 04     | pMHE7   | Streptomycin |
|              | pM7 05     | pMHE7   | Streptomycin |
|              | pM7 06     | pMHE7   | Streptomycin |
|              | pM7 07     | pMHE7   | Streptomycin |
|              | pM7 08     | pMHE7   | Streptomycin |
|              | pM7 09     | pMHE7   | Streptomycin |
|              | pM7 10     | pMHE7   | Streptomycin |
|              | pM7 11     | pMHE7   | Streptomycin |
|              | pM7 12     | pMHE7   | Streptomycin |
|              | pM7 13     | pMHE7   | Streptomycin |
|              | pM7 14     | pMHE7   | Streptomycin |
|              | pM7 15     | pMHE7   | Streptomycin |
|              | pM7 16     | pMHE7   | Streptomycin |
| 3.           | pD5 01     | pDSK5   | Kanamycin    |
|              | pD5 02     | pDSK5   | Kanamycin    |
|              | pD5 03     | pDSK5   | Kanamycin    |
|              | pD5 04     | pDSK5   | Kanamycin    |
|              | pD5 05     | pDSK5   | Kanamycin    |
|              | pD5 06     | pDSK5   | Kanamycin    |
|              | pD5 07     | pDSK5   | Kanamycin    |
|              | pD5 08     | pDSK5   | Kanamycin    |
|              | pD5 09     | pDSK5   | Kanamycin    |
|              | pD5 10     | pDSK5   | Kanamycin    |
|              | pD5 11     | pDSK5   | Kanamycin    |
|              | pD5 12     | pDSK5   | Kanamycin    |
|              | pD5 13     | pDSK5   | Kanamycin    |
|              | pD5 14     | pDSK5   | Kanamycin    |
|              | pD5 15     | pDSK5   | Kanamycin    |
|              | pD5 16     | pDSK5   | Kanamycin    |

At this point we had 18 test tubes with the E coli containing the synthesized gene and 12 test tubes with the E coli containing the synthesized promoter. We vortexed the test tubes and then we united them in groups of three (Figure 1) into Eppendorf tubes and then centrifuged them thus one Eppendorf tube holds bacteria from 3 test tubes and then we discarded the supernatant.

For the miniprep we used the GeneJET Plasmid Miniprep Kit (Thermo Fisher Scientific, Catalog number: K0502) which contains 15 mL Resuspension Solution, 15 mL Lysis Solution, 20 mL Neutralization Solution, 20 mL Wash Solution, 0.15 mL RNase A and 4 mL Elution Buffer (10 mM Tris-HCl, pH 8.5) 50 pieces of GeneJET Spin Columns 50 pieces of Collection Tubes (2 mL). To prepare the buffer we had to add the RNase A solution to the Resuspension Solution and mix, add 35 mL ethanol (96%) to the Wash Solution.

Before starting the protocol of the kit we had 6 Eppendorf tubes with the E.coli containing the synthesized gene and 4 with the E. coli containing the synthesized promoter. (Figure 1).

1. We resuspended the pelleted cells in 250  $\mu$ L of the Resuspension Solution in 2 mL Eppendorf tubes. We had the bacteria resuspended by vortexing until no cell clumps remained.
2. Then we added 250  $\mu$ L of the Lysis Solution and mixed it by inverting the tube 4-6 times until the solution became viscous and slightly clear.
3. We added 350  $\mu$ L of the Neutralization Solution and mixed it immediately by inverting the tube 4-6 times. The neutralized bacterial lysate became cloudy. We centrifuged them for 5 minutes to pellet cell debris and chromosomal DNA.
4. We transferred the supernatant to the GeneJET spin columns by decanting. We centrifuged them for 1 minute, discard the flow-through and placed the column back into the same collection tube.
5. We added 500  $\mu$ L of the Wash Solution to the GeneJET spin column, centrifuged them for 1 minute and discarded the flow-through and placed the column back into the same collection tube.
6. We repeated the wash procedure using 500  $\mu$ L of the Wash Solution, discard the flow-through and centrifuged for an additional 1 min to remove residual Wash Solution.
7. We transferred the GeneJET spin column into a fresh 1.5 mL Eppendorf tube then we added 50  $\mu$ L of the Elution Buffer (which we rewarmed to 70°C) to the center of GeneJET spin column membrane to elute the plasmid DNA.
8. We incubated them for 2 minutes at room temperature and centrifuged them for 2 minutes. Finally we discard the column and stored the purified plasmid DNAs at -20°C.

| test tubes | After protocol Eppendorf tubes | content  |
|------------|--------------------------------|--|
| G11        | G1                             | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G12        |                                |  |
| G13        |                                |  |
| G21        | G2                             | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G22        |                                |  |
| G23        |                                |  |
| G31        | G3                             | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G32        |                                |  |
| G33        |                                |  |
| G41 X      | G4 X                           | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G42 X      |                                |  |
| G43 X      |                                |  |
| G51        | G5                             | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G52        |                                |  |
| G53        |                                |  |
| G61        | G6                             | Plasmid containing<br>LDH gene +50µL<br>Elution Buffer |
| G62        |                                |  |
| G63        |                                |  |
| P11        | P1                             | Plasmid containing<br>promoter +50µL<br>Elution Buffer |
| P12        |                                |  |
| P13        |                                |  |
| P21        | P2                             | Plasmid containing<br>promoter +50µL<br>Elution Buffer |
| P22        |                                |  |
| P23        |                                |  |
| P41        | P4                             | Plasmid containing<br>promoter +50µL<br>Elution Buffer |
| P42        |                                |  |
| P43        |                                |  |
| P51        | P5                             | Plasmid containing<br>promoter +50µL<br>Elution Buffer |
| P52        |                                |  |
| P53        |                                |  |

## Miniprep

We obtained three plausible working plasmids: pMHE5, pMHE7 and pDSK5. We amplified them by letting the E. coli divide which had these plasmids inside them. We had 16 test tubes with bacterial solution for each plasmid. We vortexed the test tubes and then we united them in groups of four (Figure 1) into Eppendorf tubes and then centrifuged them thus one Eppendorf tube holds bacteria from 4 test tubes and then we discarded the supernatant.

For the miniprep we used the GeneJET Plasmid Miniprep Kit (Thermo Fisher Scientific, Catalog number: K0502) which contains 15 mL Resuspension Solution, 15 mL Lysis Solution, 20 mL Neutralization Solution, 20 mL Wash Solution, 0.15 mL RNase A and 4 mL Elution Buffer (10 mM Tris-HCl, pH 8.5) 50 pieces of GeneJET Spin Columns 50 pieces of Collection Tubes (2 mL). To prepare the buffer we had to add the RNase A solution to the Resuspension Solution and mix, add 35 mL ethanol (96%) to the Wash Solution.

Before starting the protocol of the kit we had 4 Eppendorf tubes for each type of plasmid (a total of 12 Eppendorf tubes). (Figure 1).

1. We resuspended the pelleted cells in 250  $\mu$ L of the Resuspension Solution in 2 mL Eppendorf tubes. We had the bacteria resuspended by vortexing until no cell clumps remained.
2. Then we added 250  $\mu$ L of the Lysis Solution and mixed it by inverting the tube 4-6 times until the solution became viscous and slightly clear.
3. We added 350  $\mu$ L of the Neutralization Solution and mixed it immediately by inverting the tube 4-6 times. The neutralized bacterial lysate became cloudy. We centrifuged them for 5 minutes to pellet cell debris and chromosomal DNA.
4. We transferred the supernatant to the GeneJET spin columns by decanting. We centrifuged them for 1 minute, discard the flow-through and placed the column back into the same collection tube.
5. We added 500  $\mu$ L of the Wash Solution to the GeneJET spin column, centrifuged them for 1 minute and discarded the flow-through and placed the column back into the same collection tube.
6. We repeated the wash procedure using 500  $\mu$ L of the Wash Solution, discard the flow-through and centrifuged for an additional 1 min to remove residual Wash Solution.
7. We transferred the GeneJET spin column into a fresh 1.5 mL Eppendorf tube then we added 50  $\mu$ L of the Elution Buffer (which we rewarmed to 70°C) to the center of GeneJET spin column membrane to elute the plasmid DNA.
8. We incubated them for 2 minutes at room temperature and centrifuged them for 2 minutes. Finally we discarded the co-stored the purified plasmid DNAs at -20°C.

| plasmids | test tubes(16pieces/plasmid) | before protocol Eppendorf tubes | after protocol Eppendorf tubes | content  |
|----------|------------------------------|---------------------------------|--------------------------------|--|
| PMHE5    | PM5 01                       | PM5 1-4                         | PMHE5 1-8                      | PMHE5 plasmids from test tubes PM5 01-08 +100µl Elution Buffer |
| PMHE5    | PM5 02                       |                                 |                                |  |
| PMHE5    | PM5 03                       |                                 |                                |  |
| PMHE5    | PM5 04                       |                                 |                                |  |
| PMHE5    | PM5 05                       | PM5 5-8                         | PMHE5 1-8                      | PMHE5 plasmids from test tubes PM5 01-08 +100µl Elution Buffer |
| PMHE5    | PM5 06                       |                                 |                                |  |
| PMHE5    | PM5 07                       |                                 |                                |  |
| PMHE5    | PM5 08                       |                                 |                                |  |
| PMHE5    | PM5 09                       | PM5 9-12                        | PMHE5 9-16                     | PMHE5 plasmids from test tubes PM5 09-16 +100µl Elution Buffer |
| PMHE5    | PM5 10                       |                                 |                                |  |
| PMHE5    | PM5 11                       |                                 |                                |  |
| PMHE5    | PM5 12                       |                                 |                                |  |
| PMHE5    | PM5 13                       | PM5 13-16                       | PMHE5 9-16                     | PMHE5 plasmids from test tubes PM5 09-16 +100µl Elution Buffer |
| PMHE5    | PM5 14                       |                                 |                                |  |
| PMHE5    | PM5 15                       |                                 |                                |  |
| PMHE5    | PM5 16                       |                                 |                                |  |
| PMHE7    | PM7 01                       | PM7 1-4                         | PM7 1-8                        | PMHE7 plasmids from test tubes PM7 01-08 +100µl Elution Buffer |
| PMHE7    | PM7 02                       |                                 |                                |  |
| PMHE7    | PM7 03                       |                                 |                                |  |
| PMHE7    | PM7 04                       |                                 |                                |  |
| PMHE7    | PM7 05                       | PM7 5-8                         | PM7 1-8                        | PMHE7 plasmids from test tubes PM7 01-08 +100µl Elution Buffer |
| PMHE7    | PM7 06                       |                                 |                                |  |
| PMHE7    | PM7 07                       |                                 |                                |  |
| PMHE7    | PM7 08                       |                                 |                                |  |
| PMHE7    | PM7 09                       | PM7 9-12                        | PM7 9-16                       | PMHE7 plasmids from test tubes PM7 09-16 +100µl Elution Buffer |
| PMHE7    | PM7 10                       |                                 |                                |  |
| PMHE7    | PM7 11                       |                                 |                                |  |
| PMHE7    | PM7 12                       |                                 |                                |  |
| PMHE7    | PM7 13                       | PM7 13-16                       | PM7 9-16                       | PMHE7 plasmids from test tubes PM7 09-16 +100µl Elution Buffer |
| PMHE7    | PM7 14                       |                                 |                                |  |
| PMHE7    | PM7 15                       |                                 |                                |  |
| PMHE7    | PM7 16                       |                                 |                                |  |
| PDSK5    | PDS 01                       | PDS 1-4                         | PDSK 1-8                       | PDSK5 plasmids from test tubes PDS 01-08 +100µl Elution Buffer |
| PDSK5    | PDS 02                       |                                 |                                |  |
| PDSK5    | PDS 03                       |                                 |                                |  |
| PDSK5    | PDS 04                       |                                 |                                |  |
| PDSK5    | PDS 05                       | PDS 5-8                         | PDSK 1-8                       | PDSK5 plasmids from test tubes PDS 01-08 +100µl Elution Buffer |
| PDSK5    | PDS 06                       |                                 |                                |  |
| PDSK5    | PDS 07                       |                                 |                                |  |
| PDSK5    | PDS 08                       |                                 |                                |  |
| PDSK5    | PDS 09                       | PDS 9-12                        | PDS 9-16                       | PDSK5 plasmids from test tubes PDS 09-16 +100µl Elution Buffer |
| PDSK5    | PDS 10                       |                                 |                                |  |
| PDSK5    | PDS 11                       |                                 |                                |  |
| PDSK5    | PDS 12                       |                                 |                                |  |
| PDSK5    | PDS 13                       | PDS 13-16                       | PDS 9-16                       | PDSK5 plasmids from test tubes PDS 09-16 +100µl Elution Buffer |
| PDSK5    | PDS 14                       |                                 |                                |  |
| PDSK5    | PDS 15                       |                                 |                                |  |
| PDSK5    | PDS 16                       |                                 |                                |  |

To separate our gene and promoter from the plasmid they were synthesised into first we had to digest them with the right type of restriction enzyme.(figure 1)

| Eppendorf tube      | G1E1  | G1E2  | G1E3  | G1E4  | G1E5  | G1E6  | G1E7  |
|---------------------|-------|-------|-------|-------|-------|-------|-------|
| NcoI                | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| Green Buffer        | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| G1                  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| Nuclease-free water | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL |
| Total               | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL |

| Eppendorf tube      | P1E1  | P1E2  | P1E3  | P1E4  | P1E5  | P1E6  | P1E7  |
|---------------------|-------|-------|-------|-------|-------|-------|-------|
| BglII               | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| Green Buffer        | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| P1                  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| Nuclease-free water | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL | 14 µL |
| Total               | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL |

Then we filled them into the pockets of two gels (figure 2) and run them with gel electrophoresis, we used the settings 110V with 400mA for 40 minutes. After the run we shined UV lights on the gels for just for 1 or 2 seconds in a dark room and quickly cut out the gel parts where the lines were furthest away from the pockets because those were our gene and promoter. (it is important to use UV light just for a short time not to damage the DNA)

| gel pocket  | 1.     | 2.   | 3.   | 4.   | 5.   | 6.   | 7.   | 8.   |
|-------------|--------|------|------|------|------|------|------|------|
| content     | ladder | G1E1 | G1E2 | G1E3 | G1E4 | G1E5 | G1E6 | G1E7 |
| volume (µL) | 5      | 10   | 10   | 10   | 10   | 10   | 10   | 10   |

| gel pocket  | 1.     | 2.   | 3.   | 4.   | 5.   | 6.   | 7.   | 8.   |
|-------------|--------|------|------|------|------|------|------|------|
| content     | ladder | P1E1 | P1E2 | P1E3 | P1E4 | P1E5 | P1E6 | P1E7 |
| volume (µL) | 5      | 10   | 10   | 10   | 10   | 10   | 10   | 10   |

Following this we used the GeneJET Gel Extraction Kit (Catalog number: K0692) on the sliced gel parts.

This kit contains:

- Binding Buffer (150 mL)
- Wash Buffer (concentrated) (45 mL)
- Elution Buffer (10 mM Tris-HCl, pH 8.5) (30 mL)
- GeneJET Purification Columns (preassembled with collection tubes) (250pieces)

Then we followed the protocol:

1. We placed the gel slices into a pre-weighed 1.5 mL tube and recorded the weight of the gel slices.
2. We added 1:1 volume of Binding Buffer to the gel slices (e.g. 100  $\mu$ L of Binding Buffer for every 100 mg of agarose gel)
3. We incubated the gel mixture at 60  $^{\circ}$ C for 10 min or until the gel slice is completely dissolved. (Mixed the tube by inversion every few minutes to facilitate the melting process). We ensured that the gel is completely dissolved, vortexed the gel mixture briefly before loading on the columns.
4. This step was only necessary for the promoter because it was shorter than 500bp. We added 1 gel volume of 100% isopropanol to the solubilized gel solution (e.g. 100  $\mu$ L of isopropanol should be added to 100 mg gel slice solubilized in 100  $\mu$ L of Binding Buffer) and mixed thoroughly.
5. We transferred up to 800  $\mu$ L of the solubilized gel solution to the GeneJET purification columns. We centrifuged for 1 min, then discarded the flow-through and placed the column back into the same collection tube.
6. We added 700  $\mu$ L of Wash Buffer to the GeneJET purification columns then centrifuged them for 1 min, discarded the flow-through and placed the column back into the same collection tube.
7. We centrifuged the empty GeneJET purification columns for an additional 1 min to completely remove residual wash buffer.
8. We transferred the GeneJET purification columns into clean 1.5 mL microcentrifuge tubes. We added 25  $\mu$ L of Elution Buffer to the center of each purification column membrane and centrifuged them for 1 min.
9. Finally we discarded the GeneJET purification columns and stored the purified DNA at -20  $^{\circ}$ C.

| Tube name | mass (g) | Binding Buffer ( $\mu$ L) | Wash Buffer ( $\mu$ L) | Elution Buffer ( $\mu$ L) | pH buffer | isopropanol ( $\mu$ L) | Eppendorf tube after protocol | Final volume ( $\mu$ L) |
|-----------|----------|---------------------------|------------------------|---------------------------|-----------|------------------------|-------------------------------|-------------------------|
| G1E1      | 0,88     | 880                       | 700                    | 25                        | -         | -                      | G1E1T                         | 25                      |
| G1E2      | 0,57     | 570                       | 700                    | 25                        | -         | -                      | G1E2T                         | 25                      |
| G1E3      | 0,66     | 660                       | 700                    | 25                        | -         | -                      | G1E3T                         | 25                      |
| G1E4      | 0,64     | 640                       | 700                    | 25                        | -         | -                      | G1E4T                         | 25                      |
| G1E5      | 0,6      | 600                       | 700                    | 25                        | -         | -                      | G1E5T                         | 25                      |
| G1E6      | 0,7      | 700                       | 700                    | 25                        | -         | -                      | G1E6T                         | 25                      |
| P1E1      | 0,4      | 400                       | 700                    | 20                        | -         | 400                    | P1E1T                         | 20                      |
| P1E2      | 0,54     | 540                       | 700                    | 20                        | -         | 540                    | P1E2T                         | 20                      |
| P1E3      | 0,44     | 440                       | 700                    | 20                        | -         | 440                    | P1E3T                         | 20                      |
| P1E4      | 0,39     | 390                       | 700                    | 20                        | -         | 390                    | P1E4T                         | 20                      |
| P1E5      | 0,38     | 380                       | 700                    | 20                        | -         | 380                    | P1E5T                         | 20                      |
| P1E6      | 0,29     | 290                       | 700                    | 20                        | -         | 290                    | P1E6T                         | 20                      |
| P1E7      | 0,38     | 380                       | 700                    | 20                        | -         | 380                    | P1E7T                         | 20                      |



To separate our gene and promoter from the plasmid they were synthesised into first we had to digest them with the right type of restriction enzyme, this time we cut out the gene with NcoI and EcoRV too.(figure 1)

| Eppendorf tube      | G1G1   | G1G2   |
|---------------------|--------|--------|
| NcoI                | 1,5 µl | 1,5 µl |
| EcoRV               | 1,5 µL | 1,5 µL |
| Green Buffer        | 3 µl   | 3 µl   |
| G1                  | 5 µL   | 5 µL   |
| Nuclease-free water | 19µL   | 19 µL  |
| Total               | 30 µL  | 30 µL  |

| Eppendorf tube      | P1G12 | P1G34 | P1G56 |
|---------------------|-------|-------|-------|
| BglII               | 2 µL  | 2 µL  | 2µL   |
| Green Buffer        | 2 µl  | 2 µl  | 2 µl  |
| P1                  | 8 µl  | 8 µl  | 8 µl  |
| Nuclease-free water | 8 µL  | 8 µL  | 8 µL  |
| Total               | 20 µl | 20 µl | 20µl  |

Then we filled them into the pockets of two gels (figure 2) and run them with gel electrophoresis, we used the settings 110V with 400mA for 40 minutes. After the run we shined UV lights on the gels for just for 1 or 2 seconds in a dark room and quickly cut out the gel parts where the lines were furthest away from the pockets because those were our gene and promoter. (it is important to use UV light just for a short time not to damage the DNA)

| gel pocket  | 1.     | 2.   | 3.   | 4.   | 5.   | 6.   | 7.   | 8. |
|-------------|--------|------|------|------|------|------|------|----|
| content     | ladder | G1G1 | G1G1 | G1G1 | G1G2 | G1G2 | G1G2 | -  |
| volume (µL) | 5      | 10   | 10   | 10   | 10   | 10   | 10   | -  |

| gel pocket  | 1.     | 2.    | 3.    | 4.    | 5.    | 6.    | 7.    | 8. |
|-------------|--------|-------|-------|-------|-------|-------|-------|----|
| content     | ladder | P1G12 | P1G12 | P1G34 | P1G34 | P1G56 | P1G56 | -  |
| volume (µL) | 5      | 10    | 10    | 10    | 10    | 10    | 10    | -  |

Following this we used the GeneJET Gel Extraction Kit (Catalog number: K0692) on the sliced gel parts.

This kit contains:

- Binding Buffer (150 mL)
- Wash Buffer (concentrated) (45 mL)
- Elution Buffer (10 mM Tris-HCl, pH 8.5) (30 mL)
- GeneJET Purification Columns (preassembled with collection tubes) (250pieces)

Then we followed the protocol:

1. We placed the gel slices into a pre-weighed 1.5 mL tube and recorded the weight of the gel slices.

2. We added 1:1 volume of Binding Buffer to the gel slices (e.g. 100  $\mu$ L of Binding Buffer for every 100 mg of agarose gel)
3. We incubated the gel mixture at 60  $^{\circ}$ C for 10 min or until the gel slice is completely dissolved. (Mixed the tube by inversion every few minutes to facilitate the melting process). We ensured that the gel is completely dissolved, vortexed the gel mixture briefly before loading on the columns.
4. This step was only necessary for the promoter because it was shorter than 500bp. We added 1 gel volume of 100% isopropanol to the solubilized gel solution (e.g. 100  $\mu$ L of isopropanol should be added to 100 mg gel slice solubilized in 100  $\mu$ L of Binding Buffer) and mixed thoroughly.
5. We transferred up to 800  $\mu$ L of the solubilized gel solution to the GeneJET purification columns. We centrifuged for 1 min, then discarded the flow-through and placed the column back into the same collection tube.
6. We added 700  $\mu$ L of Wash Buffer to the GeneJET purification columns then centrifuged them for 1 min, discarded the flow-through and placed the column back into the same collection tube.
7. We centrifuged the empty GeneJET purification columns for an additional 1 min to completely remove residual wash buffer.
8. We transferred the GeneJET purification columns into clean 1.5 mL microcentrifuge tubes. We added 20  $\mu$ L of Elution Buffer to the center of each purification column membrane and centrifuged them for 1 min.
9. Finally we discarded the GeneJET purification columns and stored the purified DNA at -20  $^{\circ}$ C.

| Eppendorf tube | mass (g) | Binding Buffer ( $\mu$ L) | Wash Buffer ( $\mu$ L) | Elution Buffer ( $\mu$ L) | isopropanol | Eppendorf tube after protocol | final volume ( $\mu$ L) |
|----------------|----------|---------------------------|------------------------|---------------------------|-------------|-------------------------------|-------------------------|
| G1G1           | 0,30     | 300                       | 700                    | 20                        | -           | G1G1T                         | 40                      |
| G1G2           | 0,26     | 260                       | 700                    | 20                        | -           |                               |                         |
| G1G3           | 0,24     | 240                       | 700                    | 20                        | -           | G1G2T                         | 40                      |
| G1G4           | 0,22     | 220                       | 700                    | 20                        | -           |                               |                         |
| G1G5           | 0,20     | 200                       | 700                    | 20                        | -           | G1G3T                         | 40                      |
| G1G6           | 0,32     | 320                       | 700                    | 20                        | -           |                               |                         |
| P1G1           | 0,93     | 930                       | 700                    | 40                        | 930         | P1G1T                         | 40                      |
| P1G2           | 0,80     | 800                       | 700                    | 40                        | 800         | P1G2T                         | 40                      |
| P1G3           | 0,92     | 920                       | 700                    | 40                        | 920         | P1G3T                         | 40                      |

To provide place for the promoter in the pMHE5 and pMHE7 plasmids first we digested them with BglII.(figure 1)

|                     |       |       |       |       |       |       |
|---------------------|-------|-------|-------|-------|-------|-------|
| Eppendorf tube      | M5RB1 | M5RB2 | M5RB3 | M7RB1 | M7RB2 | M7RB3 |
| sample              | pMHE5 | pMHE5 | pMHE5 | pMHE7 | pMHE7 | pMHE7 |
| sample volume       | 4 µL  | 4 µL  | 4 µL  | 4 µL  | 4 µL  | 4 µL  |
| Green Buffer        | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| BglII               | 2 µL  | 2µL   | 2 µL  | 2 µL  | 2 µL  | 2 µL  |
| nuclease-free water | 12 µL | 12 µL | 12 µL | 12 µL | 12 µL | 12 µL |
| Total               | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL | 20 µL |

Then we separated the DNA fragments from all the other parts of this solution with gel electrophoresis. We filled them into the gel pockets(figure 2) and run them on the settings 110V with 400mA for 40 minutes. After the run we shined UV lights on the gels for just for 1 or 2 seconds in a dark room and quickly cut out the gel parts where the lines were. (it is important to use UV light just for a short time not to damage the DNA)

|             |        |       |       |       |       |       |       |    |
|-------------|--------|-------|-------|-------|-------|-------|-------|----|
| gel pocket  | 1.     | 2.    | 3.    | 4.    | 5.    | 6.    | 7.    | 8. |
| content     | ladder | M5RB1 | M5RB1 | M5RB2 | M5RB2 | M5RB3 | M5RB3 | -  |
| volume (µL) | 5      | 10    | 10    | 10    | 10    | 10    | 10    | -  |
|             |        |       |       |       |       |       |       |    |
| gel pocket  | 1.     | 2.    | 3.    | 4.    | 5.    | 6.    | 7.    | 8. |
| content     | ladder | M7RB1 | M7RB1 | M7RB2 | M7RB2 | M7RB3 | M7RB3 | -  |
| volume (µL) | 5      | 10    | 10    | 10    | 10    | 10    | 10    | -  |

Following this we used the GeneJET Gel Extraction Kit (Catalog number: K0692) on the sliced gel parts.

This kit contains:

- Binding Buffer (150 mL)
- Wash Buffer (concentrated) (45 mL)
- Elution Buffer (10 mM Tris-HCl, pH 8.5) (30 mL)
- GeneJET Purification Columns (preassembled with collection tubes) (250pieces)

Then we followed the protocol:

1. We placed the gel slices into a pre-weighed 1.5 mL tube and recorded the weight of the gel slices.
2. We added 1:1 volume of Binding Buffer to the gel slices (e.g. 100  $\mu$ L of Binding Buffer for every 100 mg of agarose gel)
3. We incubated the gel mixture at 60 °C for 10 min or until the gel slice is completely dissolved. (Mixed the tube by inversion every few minutes to facilitate the melting process). We ensured that the gel is completely dissolved, vortexed the gel mixture briefly before loading on the columns.
4. We transferred up to 800  $\mu$ L of the solubilized gel solution to the GeneJET purification columns. We centrifuged for 1 min, then discarded the flow-through and placed the column back into the same collection tube.
5. We added 700  $\mu$ L of Wash Buffer to the GeneJET purification columns then centrifuged them for 1 min, discarded the flow-through and placed the column back into the same collection tube.
6. We centrifuged the empty GeneJET purification columns for an additional 1 min to completely remove residual wash buffer.
7. We transferred the GeneJET purification columns into clean 1.5 mL microcentrifuge tubes. We added 20  $\mu$ L of Elution Buffer to the center of each purification column membrane and centrifuged them for 1 min.
8. Finally we discarded the GeneJET purification columns and stored the purified DNA at -20 °C.

| Tube name | mass (g) | Binding Buffer ( $\mu$ L) | Wash Buffer ( $\mu$ L) | Elution Buffer ( $\mu$ L) | pH buffer | isopropanol ( $\mu$ L) | Eppendorf tube after protocol | Final volume ( $\mu$ L) |
|-----------|----------|---------------------------|------------------------|---------------------------|-----------|------------------------|-------------------------------|-------------------------|
| M5E1      | 0,49     | 490                       | 700                    | 20                        | -         | -                      | M5E1T                         | 40                      |
| M5E2      | 0,28     | 280                       | 700                    | 20                        | -         | -                      |                               |                         |
| M5E3      | 0,39     | 390                       | 700                    | 20                        | -         | -                      | M5E2T                         | 40                      |
| M5E4      | 0,45     | 450                       | 700                    | 20                        | -         | -                      |                               |                         |
| M5E5      | 0,38     | 380                       | 700                    | 20                        | -         | -                      | M5E3T                         | 40                      |
| M5E6      | 0,36     | 360                       | 700                    | 20                        | -         | -                      |                               |                         |
| M7E1      | 0,39     | 390                       | 700                    | 20                        | -         | -                      | M7E1T                         | 40                      |
| M7E2      | 0,38     | 380                       | 700                    | 20                        | -         | -                      |                               |                         |
| M7E3      | 0,43     | 430                       | 700                    | 20                        | -         | -                      | M7E2T                         | 40                      |
| M7E4      | 0,33     | 330                       | 700                    | 20                        | -         | -                      |                               |                         |
| M7E5      | 0,36     | 360                       | 700                    | 20                        | -         | -                      | M7E3T                         | 40                      |
| M7E6      | 0,38     | 380                       | 700                    | 20                        | -         | -                      |                               |                         |

## Protocol for ligation of the LDH gene and the plasmid

### Plasmid-LDH ligation

The ligation of the lactate-dehydrogenase (LDH) gene and the pMHE5 / pMHE7 plasmid was carried out with the Thermo Scientific T4 DNA ligase. Plasmid DNA concentration in the samples are shown by the table.

| Sample ID | DNA concentration (ng/ $\mu$ L) |
|-----------|---------------------------------|
| G1G1T     | 11.65                           |
| G1G2T     | 6.54                            |
| G1G3T     | 11.06                           |
| P1G1T     | 5.76                            |
| P1G2T     | 3.37                            |
| P1G3T     | 3.92                            |
| M5E1T     | 4.49                            |
| M5E2T     | 3.55                            |
| M5E3T     | 2.57                            |
| M7E1T     | 2.24                            |
| M7E2T     | 2.03                            |
| M7E3T     | 2.46                            |

Table

Since we knew the exact DNA content of the samples, we could carry out the ligation. The ligation mixture was composed as shown by table.

|  | I. reaction mixture | II. reaction mixture | III. reaction mixture | IV. reaction mixture |
|--|---------------------|----------------------|-----------------------|----------------------|
| Plasmid sample   | M5E1T               | M5E2T                | M7E1T                 | M7E3T                |
| Plasmid sample volume ( $\mu\text{L}$ )                  | 10                  | 13                   | 20                    | 20                   |
| Gene sample  | G1G3T               | G1G2T                | G1G1T                 | G1G3T                |
| Gene sample volume ( $\mu\text{L}$ )                     | 3                   | 4                    | 2                     | 2                    |
| Buffer volume ( $\mu\text{L}$ )                          | 4                   | 4                    | 4                     | 4                    |
| Thermo Scientific T4 DNA ligase volume ( $\mu\text{L}$ ) | 1                   | 1                    | 1                     | 1                    |
| Nuclease-free water volume ( $\mu\text{L}$ )             | 2                   | -                    | -                     | -                    |
| Total volume ( $\mu\text{L}$ )                           | 20                  | 22                   | 27                    | 27                   |

Table

The reaction mixtures were incubated at 37 °C for 10 minutes.

After the incubation, the putative construct of the ligated plasmid and gene were transferred into *Escherichia coli* JM109 competent cells with thermal shock by the following protocol:

1. On ice, to 50  $\mu\text{L}$  of *E. coli* JM109 competent cell solution 5  $\mu\text{L}$  of the reaction mix was confused in an Eppendorf-tube, then the solution was left on ice for 30 minutes.  
During that, SOC-containing soil was placed on ice as well.
2. After 30 minutes, thermal shock was accomplished at 42 °C for 60 seconds.
3. When the thermal shock was ready, the Eppendorf-tubes were placed back on ice, and right after it, 200  $\mu\text{L}$  of the SOC-containing soil was pipetted into each Eppendorf-tube.
4. The mixtures were incubated at 37 °C for 1 hour.

#### Scaling up plasmid DNA

Under Laminar Air-Flow 200  $\mu\text{L}$  of each samples was placed on Petri-dishes containing LB and Streptomycin, and as a control, 50  $\mu\text{L}$  of each samples was placed on Petri-dishes containing only LB. The samples were incubated at 37 °C for a day. At the end of the incubation time, the result was 5 colonies containing the putative pMHE5-LDH construct, and 11 colonies containing the putative pMHE7-LDH construct, both on Petri-dishes containing LB and Streptomycin. Each colony on the LB and Streptomycin containing Petri-dishes was transferred into test tubes containing 2 mL of LB and Streptomycin. These bacterial cultures were incubated at 37 °C for 2 days.

### DNA purification

Plasmid DNAs were purified using the Thermo Scientific GeneJET Plasmid Miniprep Kit.

1. We resuspended the pelleted cells in 250  $\mu\text{L}$  of the Resuspension Solution in 2 mL Eppendorf tubes. Bacteria were resuspended by vortexing until no cell clumps remained.
2. 250  $\mu\text{L}$  of the Lysis Solution was added and mixed by inverting the tube 4-6 times until the solution became viscous and slightly clear.
3. 350  $\mu\text{L}$  of the Neutralization Solution was added and mixed immediately by inverting the tube 4-6 times. The neutralized bacterial lysate became cloudy. They were centrifuged for 5 minutes to pellet cell debris and chromosomal DNA.
4. We transferred the supernatant to the GeneJET spin columns by decanting. We centrifuged them for 1 minute to discard the flow-through and placed the column back into the same collection tube.
5. We added 500  $\mu\text{L}$  of the Wash Solution to the GeneJET spin column then centrifuged them for 60 seconds, discarded the flow-through and placed the column back into the same collection tube.
6. We repeated the wash procedure using 500  $\mu\text{L}$  of the Wash Solution, discarded the flow-through and centrifuged for an additional 60 seconds to remove residual Wash Solution.
7. We transferred the GeneJET spin column into a fresh 1.5 mL Eppendorf tube then we added 50  $\mu\text{L}$  of the Elution Buffer (which was rewarmed to 70°C) to the center of GeneJET spin column membrane to elute the plasmid DNA.
8. We incubated them for 2 minute at room temperature and centrifuged them for 2 minutes. Finally we discarded the column and stored the purified plasmid DNAs at -20°C.

### Control of the putative pMHE5-LDH and pMHE7-LDH construct

We had to make sure that the ligation of the pMHE5 / pMHE7 and the LDH was successful. Therefore, we carried out a gel electrophoresis analysis. As figure shows, the (here the exact number) samples contained an approximately 1000 base pairs long DNA sequence, which must be due to the successfully inserted LDH sequence, so the ligation of the pMHE5 / pMHE7 and the LDH was proven successful.

### Restriction of the purified plasmids with BglII

The purified pMHE5-LDH and pMHE7-LDH constructs were restricted with the Thermo Scientific FastDigest BglII restriction kit. The composition of the reaction mixes are shown by table.

|   |    |
|---|----|
| BglII volume ( $\mu\text{L}$ )                  | 3  |
| Buffer volume ( $\mu\text{L}$ )                 | 3  |
| Volume of plasmid sample ( $\mu\text{L}$ )      | 8  |
| Volume of nuclease-free water ( $\mu\text{L}$ ) | 16 |
| Total volume ( $\mu\text{L}$ )                  | 30 |

Table

## Ligation of the pMHE5-LDH / pMHE7-LDH constructs and the putative promoter

The pMHE5-LDH (S5.1M-S5.5M) / pMHE7-LDH (S7.1M-S7.7M) constructs and the promoter were ligated using the Thermo Scientific T4 DNA Ligase kit. The composition of the reaction mixtures are shown by table.

|                             | I. RM | II. RM | III. RM | IV. RM | V. RM | VI. RM | VII. RM | VIII. RM | IX. RM | X. RM |
|-----------------------------|-------|--------|---------|--------|-------|--------|---------|----------|--------|-------|
| Plasmid sample              | S5.1M | S5.1M  | S5.2M   | S5.3M  | S5.4M | S5.5M  | S7.1M   | S7.2M    | S7.6M  | S7.7  |
| Plasmid sample volume (µL)  | 10    | 10     | 10      | 10     | 10    | 10     | 10      | 10       | 10     | 10    |
| Promoter sample             | P1G1T | P1G1T  | P1G1T   | P1G1T  | P1G1T | P1G1T  | P1G3T   | P1G3T    | P1G3T  | P1G3T |
| Promoter sample volume (µL) | 5     | 5      | 5       | 5      | 5     | 5      | 7       | 7        | 7      | 7     |
| Buffer volume (µL)          | 5     | 5      | 5       | 5      | 5     | 5      | 5       | 5        | 5      | 5     |
| T4 DNA ligase volume (µL)   | 2     | 2      | 2       | 2      | 2     | 2      | 2       | 2        | 2      | 2     |
| Total volume (µL)           | 22    | 22     | 22      | 22     | 22    | 22     | 24      | 24       | 24     | 24    |
| Sample ID                   | LP1   | LP2    | LP3     | LP4    | LP5   | LP6    | LP7     | LP8      | LP9    | LP10  |

RM: abbrev. of Reaction Mixture

Table

The reaction mixes were incubated at 37 °C for 10 min. After the incubation, we transformed E.coli JM109 competent cells with the samples of table (number) using thermal shock by the following protocol.

1. On ice, to 50 µL of E. coli JM109 competent cell solution 5 µL of the reaction mix was confused in an Eppendorf-tube, then the solution was left on ice for 30 minutes.  
During that, SOC-containing soil was placed on ice as well.
2. After 30 minutes, thermal shock was accomplished at 42 °C for 60 seconds.
3. When the thermal shock was ready, the Eppendorf-tubes were placed back on ice, and right after it, 200 µL of the SOC-containing soil was pipetted into each Eppendorf-tube.
4. The mixtures were incubated at 37 °C for 1 hour.

Under Laminar Air-Flow 200 µL of each samples was placed on Petri-dishes containing LB and Streptomycin, and as a control, 50 µL of each samples was placed on Petri-dishes containing only LB. The samples were incubated at 37 °C for a day. At the end of the incubation time, the result was 5 colonies containing the putative pMHE5-LDH-promoter construct, and 11 colonies containing the putative pMHE7-LDH-promoter construct, both on Petri-dishes containing LB and Streptomycin. Each colony on the LB and Streptomycin containing Petri-dishes was transferred into test tubes containing 2 mL of LB and Streptomycin. These bacterial cultures were incubated at 37 °C for 2 days. After 2 days there were colonies on the Petri-dishes containing LB and Streptomycin:

| Sample ID          | LP1 | LP2 | LP3 | LP4 | LP5 | LP6 | LP7 | LP8 | LP9 | LP10 |
|--------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|
| Number of colonies | 1   | 1?  | 2   | 7   | 2   | -   | -   | -   | 1   | 4    |

Table

The colonies were transferred into test tubes containing LB and Streptomycin.