

PCR

Goal	To replicate desired DNA so more copies of DNA are yielded
Materials	<ul style="list-style-type: none">• Pipettes• Filter tips• Sterile Water• dNTPs• Template DNA• DNA primers for region of interest• DNA polymerase buffer• DNA polymerase• 0.2 mL PCR tubes• Centrifuge• Vortex• Thermocycler
Procedure	<ol style="list-style-type: none">1. Ensure that a sterile working environment is maintained through the duration of the PCR2. Remove reaction components from -20°C freezer and thaw completely. After, vortex and centrifuge3. To 0.2 mL PCR tubes add: DNA polymerase buffer, dNTP, primers, template DNA and sterile water4. Vortex and centrifuge5. Add and pipet the DNA polymerase into the reaction6. Move reaction tubes to thermocycler7. Enter the appropriate information into thermocycler8. After thermocycling, remove tubes, centrifuge, vortex, and centrifuge a second time9. Store samples in -20°C or run immediately on a gel