

Knock-Out of LysA Gene in E.coli BW25113 by λ Red Recombination System

A. The Preparation of Targeting Vector

1	The Design of the Primers	<p>a. The knockout primers: agcacttatctggagttgttatgccacattcactgtGTGTAGGCTGGAGCTGCTTC gtcacatgcaaccagcgactaaccgcagtaaagcaATGGGAATTAGCCATGGTCC (The part of lowercase letters is the kan resistance gene and the capital one is the LysA homology arms)</p> <p>b. The detection primers: TAGTAGTCCGACGCTGGTACGTCG TTGCATAGACTCGACATAAATCGA</p>																										
2	PCR	<p>Amplify the targeting vector, using the plasmid pKD4 as template. The PCR system:</p> <table style="width: 100%; border: none;"> <tr> <td style="padding-left: 20px;">PCR MIX</td> <td style="text-align: right;">12.5μl</td> </tr> <tr> <td style="padding-left: 20px;">plasmid pKD4</td> <td style="text-align: right;">1μl</td> </tr> <tr> <td style="padding-left: 20px;">upstream primer(10μmol/L)</td> <td style="text-align: right;">1μl</td> </tr> <tr> <td style="padding-left: 20px;">downstream primer(10μmol/L)</td> <td style="text-align: right;">1μl</td> </tr> <tr> <td style="padding-left: 20px;">sterile water</td> <td style="text-align: right;">9.5μl</td> </tr> <tr> <td style="padding-left: 20px;">paraffin oil</td> <td style="text-align: right;">10μl</td> </tr> </table> <p>The PCR processes:</p> <table border="1" style="width: 100%; border-collapse: collapse; margin-left: 20px;"> <tr> <td style="padding: 2px;">Denaturation</td> <td style="padding: 2px;">94$^{\circ}$C</td> <td style="padding: 2px;">20s</td> <td rowspan="3" style="text-align: center; vertical-align: middle;">33 circulations</td> </tr> <tr> <td style="padding: 2px;">Annealing</td> <td style="padding: 2px;">52$^{\circ}$C</td> <td style="padding: 2px;">20s</td> </tr> <tr> <td style="padding: 2px;">Elongation</td> <td style="padding: 2px;">72$^{\circ}$C</td> <td style="padding: 2px;">5min</td> </tr> <tr> <td style="padding: 2px;">Final elongation</td> <td style="padding: 2px;">72$^{\circ}$C</td> <td style="padding: 2px;">5min</td> <td style="text-align: center; vertical-align: middle;">1 circulations</td> </tr> </table>	PCR MIX	12.5 μ l	plasmid pKD4	1 μ l	upstream primer(10 μ mol/L)	1 μ l	downstream primer(10 μ mol/L)	1 μ l	sterile water	9.5 μ l	paraffin oil	10 μ l	Denaturation	94 $^{\circ}$ C	20s	33 circulations	Annealing	52 $^{\circ}$ C	20s	Elongation	72 $^{\circ}$ C	5min	Final elongation	72 $^{\circ}$ C	5min	1 circulations
PCR MIX	12.5 μ l																											
plasmid pKD4	1 μ l																											
upstream primer(10 μ mol/L)	1 μ l																											
downstream primer(10 μ mol/L)	1 μ l																											
sterile water	9.5 μ l																											
paraffin oil	10 μ l																											
Denaturation	94 $^{\circ}$ C	20s	33 circulations																									
Annealing	52 $^{\circ}$ C	20s																										
Elongation	72 $^{\circ}$ C	5min																										
Final elongation	72 $^{\circ}$ C	5min	1 circulations																									
3	AGE	<p>To analyze whether the PCR product is correct.</p> <p>a. The preparation of the agarose gel</p> <table style="width: 100%; border: none;"> <tr> <td style="padding-left: 20px;">1*TAE</td> <td style="text-align: right;">30mL</td> </tr> <tr> <td style="padding-left: 20px;">agarose</td> <td style="text-align: right;">300mg</td> </tr> <tr> <td style="padding-left: 20px;">genecolour ITM</td> <td style="text-align: right;">3μl</td> </tr> <tr> <td style="padding-left: 20px;">heating</td> <td style="text-align: right;">2min</td> </tr> </table> <p>clotting in the gel container</p> <p>b. Electrophoresis 110V 35min</p> <p>c. Analysis Gel-Imaging System</p> <p>d. Gel extraction Gel Extraction Kit</p>	1*TAE	30mL	agarose	300mg	genecolour I TM	3 μ l	heating	2min																		
1*TAE	30mL																											
agarose	300mg																											
genecolour I TM	3 μ l																											
heating	2min																											
4	Purification	<p>To eliminate the false positives of pKD4</p> <p>a. Enzyme digestion System: 17μl extraction product, 1μl Dpn1, 2μl 10*buffer, 37$^{\circ}$C, 1h</p> <p>b. AGE</p>																										

B. The Transformation of Plasmid pKD46

1	The Preparation of Competent Cells(CaCl ₂)	<ul style="list-style-type: none"> a. Shake the E.coli BW25113 overnight at 30°C, then transfer it into a new liquid no resistant LB medium. b. Detect OD600 each hour. c. When OD600 is between 0.3 and 0.4, take 1.5ml bacterium solution into a 1.5ml sterile centrifuge tube, putting it into ice for 15-30min. d. Centrifuge the bacterium solution in the conditions of 4000rpm, 2min, 4°C, abandoning the supernatant. e. Add 800μl pre-cooling 0.1mol/L sterilized CaCl₂ solution to suspend the sediment. f. 30min ice bath. g. Centrifuge the bacterium solution in the conditions of 4000rpm, 2min, 4°C, abandoning the supernatant. h. Add 200μl pre-cooling 0.1mol/L sterilized CaCl₂ solution to suspend the sediment, take 50μl to transform.
2	The Transformation of Plasmid pKD46	<ul style="list-style-type: none"> a. Add 50μl competent cells and 5μl plasmid pKD46 to 1.5ml sterile centrifuge tube, taking 30min ice bath. b. After a 45s heating at 42°C water, put the tube into ice rapidly for 2min, then adding 500μl liquid no resistant LB medium. c. Shake the tube for 1.5h at 30°C. d. Centrifuge the bacterium solution in the conditions of 5000rpm for 1min, abandoning the supernatant. e. Coat the ampicillin resistant solid LB medium with the bacterium solution, than cultivated overnight at 30°C. f. Picking a single bacterial colony into the ampicillin resistant liquid LB medium, cultivated overnight at 30°C.
3	The Verification	<ul style="list-style-type: none"> a. Bacteria keeping: add 700μl bacterium solution and 700μl 50% sterile glycerin into a 1.5ml sterile centrifuge tube, storing at -20°C. b. Extract plasmids: TIANprep Mini Plasmid Kit. c. AGE

C. Electroporation

1	The Preparation of Competent Cells(glycerin)	<ul style="list-style-type: none"> a. Add 100μl pKD46-transformed E.coli BW25113 to 7ml liquid LB medium, along with 7μl ampicillin solution, shaking overnight at 30°C. b. Take 5% of bacterium solution into 30ml no resistant liquid LB medium within 30μl ampicillin solution, cultivating at 30°C shaker. c. Detect OD600 each hour until it is between 0.3 and 0.4. d. Add L-Arabinose until its concentration is 30mmol/L. Then shake it at 30°C for
---	--	---

		<p>1h.</p> <p>e. Take all of the bacterium solution into 1.5ml centrifuge tubes, ice bath for 10min.</p> <p>f. Centrifuge the bacterium solution in the conditions of 5000rpm and 4°C for 10min, abandoning the supernatant.</p> <p>g. Suspend the sediment with pre-cooling 50% sterile glycerin till 1ml, centrifuging in the conditions of 5000rpm and 4°C for 2min, repeat for 3 times.</p> <p>h. Suspend the sediment with 300µl glycerin.</p>
2	Electroporation	<p>a. Take 100µl bacterium solution and 10µl kan targeting vector into 1.5ml sterile centrifuge tube, mix, then adding into a pre-cooling electroporation cup with no bubbles.</p> <p>b. The size of electroporation cup:1mm or 2mm.</p> <p>c. Electroporation parameters: 2000kV, 25µF, controller 200Ω</p>
3	Follow- up	<p>a. Add 1 ml no resistant liquid LB medium immediately after electroporation, mix.</p> <p>b. Suck all the liquid out of the electroporation cup into a new sterile 1.5ml centrifuge tube, shaking it at 30°C, 180rpm for 1-1.5h.</p> <p>c. Coat the ampicillin resistant solid LB medium with the bacterium solution, cultivated overnight at 37°C.</p>

D. The Verification of Strain

1	Bacteria Shaking Cultivating	Picking a single bacterial colony into the ampicillin resistant liquid LB medium, cultivated overnight at 30°C.												
2	Bacteria Keeping	add 700µl bacterium solution and 700µl 50% sterile glycerin into a 1.5ml sterile centrifuge tube, storing at -20°C.												
3	Plasmid Extraction	TIANprep Mini Plasmid Kit.												
4	PCR	<p>The PCR system(using the detection primers):</p> <table> <tr> <td>PCR MIX</td> <td>12.5µl</td> </tr> <tr> <td>plasmid pKD4</td> <td>1µl</td> </tr> <tr> <td>upstream primer(10umol/L)</td> <td>1µl</td> </tr> <tr> <td>downstream primer(10umol/L)</td> <td>1µl</td> </tr> <tr> <td>sterile water</td> <td>9.5µl</td> </tr> <tr> <td>paraffin oil</td> <td>10µl</td> </tr> </table>	PCR MIX	12.5µl	plasmid pKD4	1µl	upstream primer(10umol/L)	1µl	downstream primer(10umol/L)	1µl	sterile water	9.5µl	paraffin oil	10µl
PCR MIX	12.5µl													
plasmid pKD4	1µl													
upstream primer(10umol/L)	1µl													
downstream primer(10umol/L)	1µl													
sterile water	9.5µl													
paraffin oil	10µl													
5	AGE	To analyze if the strain is correct.												