

Lab tasks Schedule

Task 1: Lab work. September 19th 2017

Tasks:

- 1- Reconstitution of the fragments
- 2- Reconstitution of the primers
- 3- PCR amplification for circ RNA circuit and CRISPR circuit

1- Reconstitution of the fragments :

Following the instructions written in the specification sheet of the IDT of each fragment. We made reconstitution of the lyophilized fragments by deionized RNase free water.

N.B : the volume of added RNase free water was calculated by : http://www.molbiol.edu.ru/eng/scripts/01_07.html

when performing the calculation we found that 105 uL of the deionized water was added to each fragment to reach the desired concentration 10ng / uL

2- Reconstitution of primers :

Following the instructions written in the specification sheet of the IDT of each fragment. We made reconstitution of the lyophilized fragments by deionized RNase free water

3- Preparation of cas9 from distribution kit according to IGEM's instructions (fragment 3 in CRISPR circuit).

4- Preparation for PCR according to protocol and amplification
PCR products

1-5 fragments of regular circuit

1-4 fragments of CRISP circuit

PCR reaction mix : according to kit taq PCT master mix kit / Cat No. 201443 (QIAGEN)

Serial	Reagent	Volume
1	Taq Master Mix	25 uL
2	RNAase free water	19 uL
3	Primer Mix (25 pmol)	5 uL
4	gDNA (1 ug / Reaction)	1 uL
	Total volume	50 uL

PCR protocol:

Was done according to the instruction manual

Annealing temperature: 55 C

No. of cycles: 12 cycles (as recommended by IDT instructions) as most of our gDNA is less than 1 KB

Task 2: Lab work 21st September

Steps:

Real Time {PCR amplification for all fragments from IdT:

Purpose: To ensure that all fragments used in Gibson assembly is pure DNA with satisfactory amplification.

- 1) Run PCR for fragments for fragments to see curves instead of gel electrophoresis (ensure quantification).
- 2) Reconstitution of primes to reach a 100 μ l concentration.
- 3) Dilution before PCR by Rnase free water to reach a concentration of 20 μ l.
- 4) Preparation of the reaction mix for PCR

Serial	Reagent	Volume
1	QuantiTect Master Mix	10 μ L
2	RNAase free water	3 μ L
3	Primer Mix (25 pmol)	2 μ L
4	gDNA (1 ug / Reaction)	5 μ L
	Total volume	20 μ L

Working on 2 parallel circuits

-determine cyclic condition of PCR

-Red rotter in ram (volume 20 μ l)

-Protocol of kit in Q-PCR

Technique optimized from after the non-specific 1st 10-15 cycles initial, log exponential, plateau.

Output: 9 ependorph 0.2

Labled from 1-9

1-5 circ-RNA

6-9 CRISPR

Gibson Assembly Protocol E5510

Gibson for target + positive control

Target: single band on gel electrophoresis

Equation: $\text{weight in ng} * 1000 / (\text{bp} * 650 \text{ dalton})$

By Gibson online calculator (NE, bio calculator) we calculated the total amount of fragments

Fragment no.	Mass in pmol	Volume in μl
Circ-RNA		
Fragment 1	1.378	1.5
Fragment 2	1.46	1.5
Fragment 3	1.63	1.5
Fragment 4	1.45	1.5
Fragment 5	2.2	2.0
CRISPR		
Fragment 6	2.26	3.0
Fragment 7	1.33	2.0
Fragment 8 (cas9)	0.00096	2.0
Fragment 9	1.62	2.0

-insert into thermal cycles to assemble the fragments into gibson master mix to do function

Output 3 ependorph 0.2 μl

1= Gibson assembly for circ-RNA circuit (fragment 1,2,3,4,5)

2= Gibson assembly for CRISPR circuit (fragment 6,7,8,9)

3= positive control for Gibson assembly (provided by the kit)

Task 3: IGEN 2017 Technical Work – 25th September, 2017

IGEM Participants:

Tasks:

- Download and analyze pictures of Gel Electrophoresis of assembled fragments
- Download and analyze pictures of PCR curves
- Transformation of assembled fragment into competent cells

Methods:

1) Preparation of a semisolid media- LB agar plates via application of (Protocol 6)

Main steps: preparation of 0.75 gm agar powder in 50 ml of LB media

Incubate in 37 degrees Celsius

Solvate ampicillin Na salt in 20 ml distilled H₂O

Pour LB agar media into 3 plates to culture a transformed Gibson PC, Crisp, and Circ

http://2011.igem.org/LB_Broth

P.S. to each 1ml of broth, add 1ul of 50 mg/mol ampicillin; therefore, for a 50 ml of broth add 250ul of 10 mg/mol ampicillin

2) Transformation into competent cell (Bioengland e.coli)

Application of Protocol 5

Gibson Assembly[®] Chemical Transformation Protocol (E5510)

<https://www.neb.com/protocols/2012/12/11/gibson-assembly-transformation-protocol-e5510>

Task 4: IGEM 2017 Technical Work – 27th September, 2017

Tasks:

1- We checked the culture media and took pictures and the results were as follows:

Culture	Results
+ve control culture plate	No bacterial Colonies
Circular RNA culture plate	Bacterial colonies Present
CRISPR culture plate	Bacterial colonies Present and of more density

2- We have followed protocol 7 and for DNA extraction “Monarch® Plasmid Miniprep Kit” with inoculation of cells with phosphate buffer saline.

3- We worked on step 3 of protocol 2 to run a new gel of Cas9 for improving characterization

4- We documented output samples and associated them to previous output samples.