The specific binding of AFP to AP-273

- Coated: take 20μL (10mg / ml) Purified protein(AFP), then take 1980μL coated buffer, 100μL per well to add in 96-wells plate, Incubate in a 37 °C incubator for 2 hours.
- 2. Emptied and dry the residual liquid, and washed the 96-wells plate twice with 300µl of washing buffer.
- 3. Added 300µL (5 mg / ml) of BSA block Buffers to each well and allowed to stand at 4°C for 1 hour
- 4. Discard the liquid, wash the plate with 300µl washing-buffer 3 times, during the last time after washing empty the liquid and pat dry
- 5. The aptamer was arranged into diluted solution of 0.2,0.4,0.6,0.8,1.0μmol/L in concentration, 100μl of the appropriate concentration of the aptamer was added to 96-wells .The blank group was added with TE buffer, and the above substances were incubated at 37 °Cfor 2h.
- 6. Discard the liquid, wash the plate with 300µl washing-buffer 3 times, during the last time after washing empty the liquid and pat dry.
- 7. Added 100µl 1: 2000 diluted HRP- streptavidin to the 96-wellss and incubated at 37 ° C for 1.5 h.
- 8. discard the liquid, wash the plate with 300 µl of washing buffer three times, then soak with 300µl of rinse buffer for 5 min, then empty the liquid and pat dry, and then wash the plate with 300µl washing buffer 2 times, during the last time after washing empty the liquid and pat dry.
- Added 10μL of TMB solution to 96-wells. After 20 mins of color development in the dark, added 50μL stop solution to terminate the reaction immediately.
 After the reaction was stopped, read the absorbance value at 450 nm.