


DNA and Genes











<u>Genetic Engineering (GE)</u> is a process where genetic material is taken from one organism and inserted into the cells of another organism

Examples of GE:

Steps of GE: 1. ________ - DNA is extracted from an organism known to have a desired trait 2. _______ - the gene of interest is located and copied 3. _______ - the gene is modified to perform in the desired way (by replacing and altering regions) 4. ______ - the gene is delivered into the tissue of the cells. We hope for it to insert into the nucleus. 5. ______ - the modified cells are allowed to grow.

Genetic Engineering: what are the advantages and disadvantages?

Advantages	Disadvantages

Strawberry DNA Extraction Experiment

Materials:

- 1. Zip bag
- 2. Strawberry
- 3. Water
- 4. Dish-washing liquid
- 5. Salt
- 6. Isopropyl alcohol





Instructions:

- 1. Put a bottle of isopropyl alcohol in a freezer.
- 2. Extraction mixture preparation:
 - a. Measure 90mL of water into a small glass container.
 - b. Measure and add 10mL dish-washing liquid to the water.
 - c. Stirl in 1/4-tsp salt and mix until the salt dissolves.
- 3. Put strawberry into the zip bag and add extraction mixture.
- 4. Remove as much air from the bag as possible and seal it closed. Smash the strawberry so that no big pieces remain.
- 5. Pour the resulting strawberry pulb and extraction mixture through a strainer and into a medium glass bowl or similar container. Use a spoon to press the mashed bits of strawberry against the strainer, forcing even more of the mixture into the container.
- 6. Pour the extracted mixture into a 50-100mL glass container. This will help to isolate the DNA on the surface of the mixture.
- 7. Add 5mL of the chilled isopropyl alcohol to the solution and hold the mixture at eye level. You're looking for a separation of material that shows up as a white layer on top.
- 8. Congrats! You have extracted the strawberry DNA!

Why do you think we add dishwashing soap?

Why do you think we add salt?

Why do you think we add isopropyl alcohol?