

**Laboratory Records: Cloning of Synthetic Catabolic Pathways of
Fluorene and Phenanthrene
Under the Control of Constitutive Promoters**

Table of Contents

1- PURPOSE	3
2- CLONING STRATEGY	3
2.1. BACKGROUND.....	3
2.2. STRATEGY.....	3
3- GENE AND PROMOTER AREAS CHARACTERISTICS	4
3.1. FLUORENE	4
3.2. PHENANTHRENE	5
3.3. DESIGN OF PROMOTER AREAS	6
4- DNA SOURCE	7
4.1. MATERIALS.....	7
4.2. DNA PREPARATION	7
5- CLONING OF INDIVIDUAL POLYCISTRONIC F1, F2 (FLUORENE) AND P1, P2 (PHENANTHRENE) BEHIND CONSTITUTIVE PROMOTERS	8
5.1. MATERIALS.....	8
5.2. METHODS.....	9
5.2.1. Preparation of promoter areas with RBS.....	9
5.2.2. Preparation of insert F1, F2, P1, and P2	9
5.2.3. Preparation of recipient vector.....	10
5.2.3. Gel purification of DNA fragments.....	11
5.2.4. Ligation.....	12
5.2.5. Transformation	14
5.3. RESULTS.....	14

5.3.1. Transformation Results.....	14
5.4. Clone Verification.....	18
5.4. CLONES DESCRIPTION	22
6- CLONING OF FULL LENGTH (SAME PLASMID) POLYCYSTRONIC F1, F2 (FLUORENE) AND P1, P2 (PHENANTHRENE) BEHIND CONSTITUTIVE PROMOTERS	22
6.1. PHENANTHRENE	22
6.1.1. Cloning of Phenanthrene Full Length (P1 and P2 on same plasmid).....	22
6.1.1.1. Preparation of fragments for Phenanthrene Cloning.....	22
6.1.1.2. Ligation Phenanthrene Cloning.....	24
6.1.2. Transformation	26
6.1.3. Transformation Results.....	26
6.1.4. Clone Verification.....	28
6.2. FLUORENE	32
6.2.1. Cloning of Fluorene Full Length (F1 and F2) into a Vector with p15a Ori.....	32
6.2.1.1. Preparation of fragments for Fluorene Cloning.....	32
6.2.1.2. Ligation Fluorene Cloning	34
6.2.1.3. Transformation	35
6.2.1.4. Transformation Results.....	36
6.2.1.5. Clone Verification.....	37
7- SEQUENCING	39
8- TRANSFORMATION OF FULL LENGTH PHENANTHRENE AND FLUORENE CLONES INTO E.COLI BL-21	40
8.1. MATERIALS.....	40
8.2. METHODS.....	40
9- CLONING OF FLUORENE FULL LENGTH F1 AND F2 INTO pSB1C3 FOR DEPOSIT IN REGISTRY	42
9.1. PREPARATION OF DNA FRAGMENTS FOR CLONING OF FLUORENE F1 AND F2 INTO pSB1C3 FOR DEPOSIT IN REGISTRY	42
9.2. LIGATION OF DNA FRAGMENTS F1 AND F2 INTO pSB1C3.....	43
9.3. TRANSFORMATION	44
9.4. TRANSFORMATION RESULTS	45
9.5. CLONE VERIFICATION.....	46
9.6. VERIFICATION BY SEQUENCING OF FULL LENGTH CLONE F1_F2	48
10- GLYCEROL STOCK.....	49
11- SEQUENCING FILES	54

1- PURPOSE

- Cloning of fluorene and phenanthrene catabolic pathways under the control of constitutive promoters.
- Testing of 3 different constitutive promoters driving the catabolic pathways in *E.coli* and in other bacteria selected for gene augmentation.

2- CLONING STRATEGY

2.1. BACKGROUND

The catabolic pathway for fluorene and phenanthrene were synthesized as polycistronic operons with the codon optimized for expression in *E.coli*.

The polycistronic catabolic pathways for fluorene and phenanthrene were split into 2 parts for several reasons, F1 and F2 for fluorene and P1 and P2 for phenanthrene:

- (i) To facilitate the synthesis of the genes (cost-effective and in a timely manner);
- (ii) To ensure a good level of expression of the polycistronic genes;
- (iii) To determine if there are orientations of the polycistronic operons that may be more favorable for expression; and
- (iv) To minimize toxicity issues that may arise when the full pathway is synthesized with all the genes in single bacteria.

We have tested the full pathway P1_P2 and F1_F2 under the control of the T7LacZ inducible promoter. From that experiment, we have learned that the expression of the enzymes involved in the fluorene catabolic pathway may be slightly toxic as illustrated by slow bacterial growth. Therefore, the fluorene pathway will be cloned into a low copy number plasmid (origin of replication: p15a). The expression of the enzymes involved in the phenanthrene catabolic pathway did not appear to alter bacterial growth when the pathway is under the control of the T7LacZ inducible promoter. Therefore, the phenanthrene catabolic pathway will be cloned into pSB1C3 (high copy number plasmid). The 2 plasmids are compatible in *E.coli*, so they can be transformed into the same bacterial strain for the degradation of both fluorene and phenanthrene

P1: Synthetic phnF, phnE, phnC, phnD

P2: Synthetic phnAc, phnAd, phnB

F1: Synthetic flnB, dbfA1, dbfA2

F2: Synthetic flnE, flnD1, ORF16, flnC

2.2. STRATEGY

Phenanthrene Catabolic Pathway:

Step 1: Design 3 promoter area containing the prefix sequence, followed by one of the 3 constitutive promoters of part BBa_J23100 or BBa_J23101 or BBa_J23110 followed by a Ribosome Binding Site of part BBa_B0034, followed by the suffix sequence.

Step 2: Order sequence at IDT to be cloned into pIDT_kanamycin vector to facilitate subsequent cloning into vectors with different antibiotic resistance genes (=pIDT is a pUC plasmid)

Step 3: Prepare, ligate and transform the following Linearized DNA fragments to obtain each pathway P1 and P2 under the control of 3 different constitutive promoters on two separate plasmids:

Linearized promoter area fragments as EcoRI/SpeI
 Linearized each of the pathway fragment (P1 and P2) as XbaI/PstI
 Linearized vector pSB1C3 as EcoRI/PstI

Step 4: Check clones by digestion for insertion of correct fragments

Step 5: Prepare, ligate, and transform the following Linearized to obtain pathway P1 and P2 under the control of 3 different constitutive promoters on the same plasmid:

Linearized promoter area + P1 fragments + vector pSB1C3 as SpeI/PstI
 Linearized promoter area + P2 fragments as XbaI/PstI

Step 6: Check clones by digestion and sequencing for insertion of correct fragments

Step 7: Transform correct clones into E.coli BL-21 for gene expression and growth in LB, minimal medium supplemented with various sources of carbons (glucose and/or phenanthrene).

Fluorene Catabolic Pathway:

Step 1: Use the promoter areas designed above [3 promoter area containing the prefix sequence, followed by one of the 3 constitutive promoters of part BBa_J23100 or BBa_J23101 or BBa_J23110 followed by a Ribosome Binding Site of part BBa_B0034, followed by the suffix sequence].

Step 2: Prepare, ligate and transform the following Linearized DNA fragments to obtain each pathway F1 and F2 under the control of 3 different constitutive promoters on two separate plasmids:

Linearized promoter area fragments as EcoRI/SpeI
 Linearized each of the pathway fragment (F1 and F2) as XbaI/PstI
 Linearized vector pSB1C3 as EcoRI/PstI

Step 3: Check clones by digestion for insertion of correct fragments

Step 4: Prepare, ligate, and transform the following Linearized to obtain pathway F1 and F2 under the control of 3 different constitutive promoters on the same plasmid:

Linearized promoter area + F1 fragments as EcoRI/SpeI
 Linearized promoter area + F2 fragments as XbaI/PstI
 Linearized vector p15a as EcoRI/PstI [pSB3T5 Tetracycline]


Step 5: Check clones by digestion and sequencing for insertion of correct fragments

Step 6: Transform correct clones into E.coli BL-21 for gene expression and growth in LB, minimal medium supplemented with various sources of carbons (glucose and/or phenanthrene).


3- GENE AND PROMOTER AREAS CHARACTERISTICS

3.1. FLUORENE


Map of Native Fluorene (Upper pathway) – 11859 bp



Map of Synthetic Fluorene Catabolic Pathway (Insert 1) – 3219 bp:




Map of Synthetic Fluorene Catabolic Pathway (Insert 2) – 3545 bp:




3.2. PHENANTHRENE


Map of Native Phenanthrene (Upper pathway) – 11451 bp



Map of Synthetic Phenanthrene Catabolic Pathway (Insert 1) – 4227bp:



Map of Synthetic Phenanthrene Catabolic Pathway (Insert 2) – 3174bp:



Plasmid Designation	Description	Size of Insert	Vector
FLUORENE-Insert 1=F1	Synthetic flnB, dbfA1, dbfA2	3219 bp	pUC57 Ampicillin
FLUORENE-Insert 2=F2	Synthetic flnE, flnD1, ORF16, flnC	3545 bp	pUC57 Ampicillin
PHENANTHRENE-Insert 1=P1	Synthetic phnF, phnE, phnC, phnD	4227 bp	pUC57 Ampicillin
PHENANTHRENE-Insert 2=P1	Synthetic phnAc, phnAd, phnB	3174 bp	pUC57 Ampicillin

3.3. DESIGN OF PROMOTER AREAS

3.3.1. Assembly using Constitutive_Promoter_BBa_J23100

>BBa_J23100 Part-only sequence (35 bp)
ttgacggctagctcagtcctaggtacagtgctagc

Sequence:

>BBa_J23100_RBS

GAATTCGCGGCCGCTTCTAGAttgacggctagctcagtcctaggtacagtgctagc**aaagaggagaaa**ACTAGTAGCGGCCGCT**GCAG**

Description:

EcoRI NotI XbaI Constitutive_Promoter_BBa_J23100 RBS_BBa_B0034 SpeI NotI PstI

3.3.2. Assembly using Constitutive_Promoter_BBa_J23101

>BBa_J23101 Part-only sequence (35 bp)
tttacagctagctcagtcctaggtattatgctagc

Sequence:

>BBa_J23101_RBS

GAATTCGCGGCCGCTTCTAGAtttacagctagctcagtcctaggtattatgctagc**aaagaggagaaa**ACTAGTAGCGGCCGCT**GCAG**

Description:

EcoRI NotI XbaI Constitutive_Promoter_BBa_J23101 RBS_BBa_B0034 SpeI NotI PstI

3.3.3. Assembly using Constitutive_Promoter_BBa_J23110

>BBa_J23101 Part-only sequence (35 bp)
tttacagctagctcagtcctaggtattatgctagc

Sequence:

>BBa_J23110_RBS

GAATTCGGCGCCGCTTCTAGAttacggctagctcagtcctaggtacaatgctagc**aaagaggagaaaACTAGTAGCGGCCGCTGCAG**

Description:

EcoRI NotI XbaI Constitutive_Promoter_BBa_J23110 RBS_BBa_B0034 SpeI NotI PstI

3.3.4. Alignment of 3 Promoter Areas

CLUSTAL O(1.2.4) multiple sequence alignment

```

BBa_J23101_RBS      GAATTCGGCGCCGCTTCTAGAttacagctagctcagtcctaggtattatgctagcaaag
BBa_J23100_RBS      GAATTCGGCGCCGCTTCTAGAttgacggctagctcagtcctaggtacagtgctagcaaag
BBa_J23110_RBS      GAATTCGGCGCCGCTTCTAGAttacggctagctcagtcctaggtacaatgctagcaaag
*****
BBa_J23101_RBS      aggagaaaACTAGTAGCGGCCGCTGCAG
BBa_J23100_RBS      aggagaaaACTAGTAGCGGCCGCTGCAG
BBa_J23110_RBS      aggagaaaACTAGTAGCGGCCGCTGCAG
*****

```

3.3.5. Promoter Areas Designation

Plasmid Designation	Description	Size of Insert	Vector
100	BBa_J23100_RBS	88 bp	pIDT_Kanamycin
101	BBa_J23101_RBS	88 bp	pIDT_Kanamycin
110	BBa_J23110_RBS	88 bp	pIDT_Kanamycin

4- DNA SOURCE

4.1. MATERIALS

Promoters (100, 101, 110):

- Reagent grade water
- The promoter areas were designed by CCA-IGEM-Team 2017.
- The promoter areas were synthesized by IDT.
- The synthetic promoter areas were delivered to us lyophilized.

Polycistronic codons (F1, F2, P1, P2):

- Reagent grade water
- The mini-genes were designed by CCA-IGEM-Team 2017.
- The mini-genes were synthesized by Genscript.
- The genes were delivered to us lyophilized.

4.2. DNA PREPARATION

The vials containing the lyophilized DNA (~4 µg) were spun down before opening the vials for the first time 16 µL of 0.2µm filtered water was added to the lyophilized powder using a P20 pipet.

After closing the tubes, they were vortexed for 2-3 minutes and were allowed to sit at 60-65°C for 15 minutes to resuspend the DNA. The tubes were then spun. Aliquots were taken to start cloning. All stocks are stored at -20°C.

Date of preparation for promoter synthetic areas: 1-Aug-2017

Date of preparation for Fluorene synthetic genes: 20-Jul-2017

Date of preparation for Phenanthrene synthetic genes: 24-July-2017

5- CLONING OF INDIVIDUAL POLYCISTRONIC F1, F2 (FLUORENE) AND P1, P2 (PHENANTHRENE) BEHIND CONSTITUTIVE PROMOTERS

5.1. MATERIALS

- a. Synthetic DNA [Fluorene F1, F2] and [Phenanthrene P1 and P2]
 - Synthetic F1 : flnB, dbfA1, dbfA2: 3219 bp
 - Synthetic F2: flnE, flnD1, ORF16, flnC: 3545 bp
 - Synthetic P1: phnF, phnE, phnC, phnD: 4227 bp
 - Synthetic P2: phnAc, phnAd, phnB: 3174 bp
- b. Synthetic DNA [Promoter areas 100, 101, 110]
 - Synthetic Promoter 100: 88 bp
 - Synthetic Promoter 101: 88 bp
 - Synthetic Promoter 110: 88 bp
- c. Agarose, 100 g, Fisher, Cat No. BP-164-100
- d. 50XTAE Electrophoresis Buffer, 1L, (1X: 40 mM Tris, 20mM Acetic Acid, 1 mM EDTA), Thermofisher, Cat No. B49
- e. Sybr Safe DNA Gel Stain, Invitrogen, Cat No. S33102
- f. 1kb plus DNA ladder DNA marker, Thermo Scientific, Cat No. SM1334
- g. 10X Blue Juice DNA loading buffer, Invitrogen, 10816-015
- h. Reagent grade water
- i. Gel DNA Recovery Kit, Zymo Research, Cat No. D4007
- j. 1.5 mL tubes
- k. Vortex
- l. Ice bucket and ice
- m. Cryo box for restriction enzyme
- n. P1000 and P200 with corresponding tips

5.2. METHODS

5.2.1. Preparation of promoter areas with RBS

Date: 1-Aug-2017

- a. The promoters/RBS regions were designed by the CCA_IGEM team and synthesized by IDT.
- b. For the writing on plates, the designations for the promoters + RBS were as follows:
 - 1) BBa_J23100 + RBS =100
 - 2) BBa_J23101 + RBS =101
 - 3) BBa_J23110 + RBS =110
- c. Double Digestion of promoter + RBS region with restriction enzymes EcoRI and SpeI
- d. Restriction Digest Set-up (20µl reaction) in a 1.5 mL tube
- e. Add the reagents as described below.
- f. Incubate the tube at 37°C for 15 min
- g. Note: there is no need to add a loading buffer because the digestion buffer already includes it.
- h. Load the reaction on a 1.5 % Agarose gel, TAE
- i. After running the electrophoresis for 1 hour at 80V, cut the linearized band.
- j. Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit
- k. Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

Component	Volume	Final condition
Reagent grade water	12.0 µl	
10X Buffer	2.0 µl	1x
Promoter +RBS region	4.0µl	~ 1 µg Note: the fragment to recover is of small size (<88bp), so needs more DNA for a reasonable recovery from the gel
EcoRI	1.0 µl	
SpeI	1.0 µl	

5.2.2. Preparation of insert F1, F2, P1, and P2

Date: 31-Jul-2017

- a. The following polycistronic sequences are cloned behind the 3 different promoter regions:
 1. For fluorene catabolic pathway
 - a. FLUO-Insert 1 or F1 [Synthetic flnB, dbfA1, dbfA2]
 - b. FLUO-Insert 2 or F2 [Synthetic flnE, flnD1, ORF16, flnC]
 2. For phenanthrene catabolic pathway
 - a. PHE-Insert 1 or P1 [Synthetic phnF, phnE, phnC, phnD]
 - b. PHE-Insert 2 or P2 [Synthetic phnAc, phnAd, phnB]
- b. The fragments are digested with 2 enzymes: XbaI and Pst I

- c. Digest plasmid containing synthetic gene designated F1, F2, P1, and P2 with the restriction enzyme XbaI and PstI
- d. Set-up restriction digestion (15 µl reaction) in a 1.5 mL tube
- e. Turn on water bath at 37°C
- f. Add the reagents in the order and with volume described in the table below.
- g. Spin the tube briefly for 15 seconds at 10,000 rpm
- h. Incubate the tube at 37°C for 15 min
- i. Note: there is no need to add a loading buffer because the digestion buffer already has it.
- j. Load the digestion reaction on a 1% Agarose gel, TAE
- k. After running the electrophoresis for 2 hours at 80V, cut out with a razor blade the linearized band of the desired size.
- l. Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit.
- m. Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

Component	Volume	Final condition
Reagent grade water	10.0µl	
10X Buffer	1.5 µl	1 X
Preparation of Insert F1, F2, P1 and P2	1.5µl	~0.4 µg
XbaI	1.0 µl	
PstI	1.0 µl	

5.2.3. Preparation of recipient vector


- a. The recipient vector is pSB1C3 (carrying the chloramphenicol resistance gene)
- b. Double Digestion of promoter + RBS region with EcoRI and PstI
- c. Restriction Digest Set-up (15 µl reaction) in 3 X 1.5 mL tube
- d. Add the reagents as described below.
- e. Incubate the tube at 37°C for 15 min
- f. Note: there is no need to add a loading buffer because the digestion buffer already has it.
- g. Load the reaction on a 1.0 % Agarose gel
- h. After running the electrophoresis for 2 hours at 80V, cut the linearized band.
- i. Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit
- j. Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

Component	Volume	Final condition
Reagent grade water	10.0 µl	


10X Buffer	1.5 μ l	1x
pSB1C3	1.5 μ l	\sim 0.2 μ g
EcoRI	1.0 μ l	
SpeI	1.0 μ l	

5.2.3. Gel purification of DNA fragments

Date: 31-Jul-2017



Date: 1-Aug-2017



Gel 1


- 1) pSB1C3 vector digested with EcoR1 and Pst 1
- 2) pSB1C3 vector digested with EcoR1 and Pst 1
- 3) pSB1C3 vector digested with EcoR1 and Pst 1
(We did 3 digestions in parallel to prepare a large stock of this linearized vector as it is needed for multiple cloning)
- 4) 1 kb molecular weight ladder
- 5) Fluorene insert 1 XbaI/PstI
- 6) Fluorene insert 2 XbaI/PstI
- 7) Phenanthrene insert 1 XbaI/PstI
- 8) Phenanthrene insert 2 XbaI/PstI

Gel 2

- 9) Promoter [BBa_J23100](#) digested with EcoRI and SpeI
- 10) Promoter [BBa_J23101](#) digested with EcoRI and SpeI
- 11) Promoter [BBa_J23110](#) digested with EcoRI and SpeI
- 12) 1 Kb molecular weight ladder

Arrows indicate the fragments that were cut out and gel-purified and used for ligation. Fragments were visible on the gels but not necessarily on the picture.

DNA Ladder



5.2.4. Ligation

Date: 1-Aug-2017 and 2-Aug-2017

Materials

- T4 DNA Ligation Buffer. Invitrogen, Cat No. 46-0114
- T4 DNA Ligase, Thermo Scientific, Cat No. K1231
- Reagent grade water, NERL, Cat No. 98555
- 1.5 mL tube
- Vortex
- Pipet and tips
- Ice bucket and ice

Methods

- Set ligation as shown below in 1.5 mL tube.
- The tubes were incubated at room temperature for 1 hour.
- The tubes were then transferred to ice.

Ligation Condition (with insert)

Component	Volume (μ l)
Reagent grade water	4.5 μ l
10X T4 ligation buffer	1.0 μ l
EcoR1-Pst1 Linearized - pSB1C3	0.5 μ l
EcoR1/SpeI - Linearized –promoter BBa_J23100 BBa_J23101 BBa_J23110	2.0 μ l
XbaI/Pst1 - Linearized -Insert [Catabolic pathway] Fluorene 1 Or Fluorene 2 Or Phenanthrene 1 Or Phenanthrene 2	1.5 μ l
T4 DNA Ligase (5Weiss/ μ l)	0.5 μ l

Control Ligation Condition (no insert)

Component	Volume (μ l)
Reagent grade water	8.0 μ l
10X T4 ligation buffer	1.0 μ l
EcorR1-Pst1 pSB1C3	0.5 μ l

T4 DNA Ligase (5Weiss/ μ l)	0.5 μ l
---------------------------------	-------------

5.2.5. Transformation

Date: 1-Aug-2017 and 2-Aug-2017

Materials

- a. LB Chloramphenicol 34 μ g/mL agar plates, Cat No. Teknova, L1017
- b. DH5a competent cells, Invitrogen, Cat 18265-017
- c. SOC (Recovery Medium), Lucigen, Cat No. F98226
- d. 1.5 mL tube
- e. Vortex
- f. Pipet and tips
- g. Ice bucket and ice
- h. Water bath (42°C)
- i. Incubator (37°C)

Method

- a. Transform ligated DNA into E.coli DH5a chemically competent cells
- b. Turn on incubator-shaker at 37°C.
- c. Turn on incubator for plates at 37°C.
- d. Set up water bath at 42°C.
- e. Bring to room temperature S.O.C medium.
- f. Bring LB plates supplemented with appropriate antibiotic (chloramphenicol 34 μ g/mL) at room temperature.
- g. Thaw competent cells on ice.
- h. Aliquots competent cells in as many tubes as needed.
- i. Add 5 μ l ligation mix to 50 μ l competent cells to DNA and swirl gently to mix.
- j. Incubate on ice for 20 minutes
- k. Heat shock at 42°C in heat block for 30 seconds. Quickly return to ice and let it sit for 2 min.
- l. Add 200 μ l of 18-25°C SOC medium and transfer mixture to 15mL Falcon tube
- m. Incubate in shaker at 37°C, 225 rpm for 30 min
- n. Plate 2 volumes (50 μ L and ~100 μ L) of the mixture onto 2 different plates of LB agar plates supplemented with Chloramphenicol 34 μ g/mL
- o. Incubate plates at 37°C overnight
- p. Count colonies
- q. Initiate cultures in LB with the appropriate antibiotic from individual clones to analyze clones

5.3. RESULTS

5.3.1. Transformation Results

Date: 25-Jul-2017

Transformation results of mixture plated on 1-Aug-2017 and 2-Aug-2017; Readout on 2-Aug-2017 and 3-Aug-2017

Description	Number of colonies (50 µL volume plating)	Number of clones analyzed by digestion
Phenanthrene		
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23100/RBS_BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD_Ter_BBa_B0015: (4115 bp, XbaI/PstI)		2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23101/RBS_BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD_Ter_BBa_B0015: (4115 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23110/RBS_BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD_Ter_BBa_B0015: (4115 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23100/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Phenanthrene Synthetic P2: phnAc, phnAd, phnB_Ter_BBa_B0015: (3062 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23101/RBS_BBa_B0034 (68 bp, EcoRI/SpeI) +Phenanthrene Synthetic P2: phnAc, phnAd, phnB_Ter_BBa_B0015: (3062 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23110/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Phenanthrene Synthetic P2: phnAc, phnAd, phnB_Ter_BBa_B0015: (3062 bp, XbaI/PstI)	>300	2
Fluorene		
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23100/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F1: flnB, dbfA1, dbfA2_Ter_BBa_B0015: (3107 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23101 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F1: flnB, dbfA1, dbfA2_Ter_BBa_B0015: (3107 bp, XbaI/PstI)	>300	2

Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23110/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F1: flnB, dbfA1, dbfA2_Ter_BBa_B0015: (3107 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23100/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F2: flnE, flnD1, ORF16, flnC_Ter_BBa_B0015: (3433 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23101/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F2: flnE, flnD1, ORF16, flnC_Ter_BBa_B0015: (3433 bp, XbaI/PstI)	>300	2
Vector (2080 bp, EcorI/PstI) + Promoter BBa_J23110/RBS_BBa_B0034 (68 bp, EcoRI/SpeI)+ Fluorene Synthetic F2: flnE, flnD1, ORF16, flnC_Ter_BBa_B0015: (3433 bp, XbaI/PstI)	>300	2
Vector alone		
Vector_pSB1C3_ (2080 bp, EcorI/PstI)	3	




Figure. LB-Chloramphenicol plates with transformation of E.coli DH5a cells with ligation of vector pSB1C3 alone, vector pSB1C3 + promoter 100 or 101 or 110 + Phenanthrene-1 fragment.

Promoter area 100: part BBa_J23100 +RBS_BBa_B0034
 Promoter area 101: part BBa_J23101 +RBS_BBa_B0034
 Promoter area 110: part BBa_J23110 +RBS_BBa_B0034

Phenanthrene-1= P1: Synthetic phnF, phnE, phnC, phnD




Figure. LB-Chloramphenicol plates with transformation of E.coli DH5a cells with ligation of vector pSB1C3 alone, vector >300+ promoter 100 or 101 or 110 + Phenanthrene-2 fragment.

Promoter area 100: part BBa_J23100 +RBS_BBa_B0034

Promoter area 101:part BBa_J23101 +RBS_BBa_B0034

Promoter area 110:part BBa_J23110 +RBS_BBa_B0034

Phenanthrene-2= P2: Synthetic phnAc, phnAd, phnB



Figure. LB-Chloramphenicol plates with transformation of E.coli DH5a cells with ligation of vector pSB1C3 alone, vector + promoter pSB1C3 100 or 101 or 110 + Fluorene-1 fragment.

Promoter area 100: part BBa_J23100 +RBS_BBa_B0034

Promoter area 101:part BBa_J23101 +RBS_BBa_B0034

Promoter area 110:part BBa_J23110 +RBS_BBa_B0034

Fluorene-1= F1: Synthetic flnB, dbfA1, dbfA2




Figure. LB-Chloramphenicol plates with transformation of E.coli DH5a cells with ligation of vector pSB1C3 alone, vector pSB1C3 + promoter 100 or 101 or 110 + Fluorene-2 fragment.

Promoter area 100: part BBa_J23100 +RBS_BBa_B0034

Promoter area 101:part BBa_J23101 +RBS_BBa_B0034

Promoter area 110:part BBa_J23110 +RBS_BBa_B0034

Fluorene-2= F2: Synthetic flnE, flnD1, ORF16, flnC

5.4. Clone Verification

Date: 4-Aug-2017

Materials

- DNA extraction kit, Zymo Research
- LB liquid medium
- Antibiotic stock solution (chloramphenicol 34 mg/mL)
- Vortex
- Pipet and tips
- 15 mL culture tube
- Incubator-Shaker
- 10 mL pipette
- Pipet aid
- Vortex
- Rack
- Toothpick

Culture Set-up

- Grow selected number of colonies in 3 mL LB medium supplemented with Chloramphenicol (35 μ g/mL) overnight at 37°C in the incubator/shaker at 220 rpm.
- Spin down 1.5 mL of culture from each culture and isolate plasmid DNA using the Zymo Research DNA Isolation kit.
- DNA preparation is resuspended in a final volume of **35 μ L** reagent grade water.
- Store remaining 1mL of culture for glycerol stock preparation.
- Check clones by digestion.

Set-up Digestion for Clone Verification

Date: 5-Aug-2017

- a. Set up digestion as shown below.
- b. This is a 1% gel agarose gel, TAE, prepared with Sybr green dye.
- c. Double digestion of miniprep DNA with EcoRI/SpeI for F1+promoter area and P1+promoter area clones
- d. Double digestion of miniprep DNA with XbaI/PstI for F2+promoter area and P2+promoter area clones
- e. Set-up restriction digestion (20µl reaction) in a 1.5 mL tube
- f. Turn on water bath at 37°C
- g. Add the reagents in the order and with volume described in the table below.
- h. Spin the tube briefly for 15 seconds at 10,000 rpm
- i. Incubate the tube at 37°C for 15 min
- j. The reaction already have a loading buffer
- k. Load the digestion reaction on a 1% Agarose gel, TAE
- l. After running the electrophoresis for 2 hours at 80V, cut out with a razor blade the linearized band of the desired size.
- m. Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit.
- n. Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

Component	Volume	Final condition
Reagent grade water	10.5 µl	
10X BufferFastDigest Green Buffer	1.5 µl	1 X
Plasmid miniprep	1.5µl	~0.2-0.4 µg
Restriction Enzyme 1	1.0 µl	
Restriction Enzyme 2	1.0 µl	

Expected Fragments Size for Verification

Phenanthrene:


EcoR/SpeI digestion for P1 based clones

XbaI/PstI digestion for P2 based clones

Description	Expected size
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23100 _BBa_B0034 (68 bp, EcoRI/SpeI) +	Vector + 4183 bp

Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD: (4115 bp, XbaI/PstI)	
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23101 _BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD: (4115 bp, XbaI/PstI)	Vector + 4183 bp
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23110 _BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P1: phnF, phnE, phnC, phnD: (4115 bp, XbaI/PstI)	Vector + 4183 bp
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23100 _BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P2: phnAc, phnAd, phnB: (3062 bp, XbaI/PstI)	Vector + 4130 bp
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23101 _BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P2: phnAc, phnAd, phnB: (3062 bp, XbaI/PstI)	Vector + 4130 bp
Vector (2080 bp, EcoRI/PstI) + Promoter BBa_J23110 _BBa_B0034 (68 bp, EcoRI/SpeI) + Phenanthrene Synthetic P2: phnAc, phnAd, phnB: (3062 bp, XbaI/PstI)	Vector + 4130 bp

Gel Electrophoresis Picture for P1 derived Clones



- 1) Clone CCA-23 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23100](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 2) Clone CCA-24 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23100](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 3) Clone CCA-25 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23100](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 4) Clone CCA-26 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23101](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 5) Clone CCA-27 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23101](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 6) Clone CCA-28 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23101](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 7) Clone CCA-29 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23110](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 8) Clone CCA-30 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23110](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 9) Clone CCA-31 [Phenanthrene insert 1 XbaI/PstI + Promoter [BBa_J23110](#) EcoRI/SpeI+ pSB1C3 EcoRI/PstI]
- 10) BBa_J04450, pSB3T5 Tetracycline, reporter gene, EcoRI- PstI
- 11) DNA Ladder

Result for P1 cloning:

CCA-23 , CCA-26, and CCA29 were first checked for presence of insert using a EcoRI/PstI digestion. As an insert of the correct size was observed on the gel, they were selected for the second step of cloning and digested with SpeI/PstI.


Gel Electrophoresis Picture for F1 and F2 and P2 derived ClonesElectrophoresis Gel 1

Upper Gel 1

- 1) DNA Miniprep Clone CCA-30 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23100 + pSB1C3]
- 2) DNA Miniprep Clone CCA-31 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23100 + pSB1C3]
- 3) DNA Miniprep Clone CCA-32 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23101 + pSB1C3]
- 4) DNA Miniprep Clone CCA-33 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23101 + pSB1C3]
- 5) DNA Miniprep Clone CCA-34 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23110 + pSB1C3]
- 6) DNA Miniprep Clone CCA-35 Digestion EcoRI/SpeI [Fluorene insert 1 + Promoter BBa_J23110 + pSB1C3]
- 7) DNA Miniprep Clone CCA-36 Digestion XbaI/PstI [Fluorene insert 2 + Promoter BBa_J23100 + pSB1C3]
- 8) DNA Miniprep Clone CCA-37 Digestion XbaI/PstI [Fluorene insert 2 + Promoter BBa_J23100 + pSB1C3]
- 9) DNA Miniprep Clone CCA-38 Digestion XbaI/PstI [Fluorene insert 2 + Promoter BBa_J23101 + pSB1C3]
- 10) DNA Miniprep Clone CCA-39 Digestion XbaI/PstI [Fluorene insert 2 + Promoter BBa_J23101 + pSB1C3]
- 11) DNA Ladder (Molecular Weight Marker)

Lower Gel 1

- 1) DNA Miniprep Clone CCA-40 Digestion EcoRI/SpeI [Fluorene insert 2 + Promoter BBa_J23110 + pSB1C3]
- 2) DNA Miniprep Clone CCA-41 Digestion XbaI/PstI [Fluorene insert 2 + Promoter BBa_J23110 + pSB1C3]
- 3) DNA Miniprep Clone CCA-42 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23100 + pSB1C3]
- 4) DNA Miniprep Clone CCA-43 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23100 + pSB1C3]
- 5) DNA Miniprep Clone CCA-44 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23101 + pSB1C3]
- 6) DNA Miniprep Clone CCA-45 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23101 + pSB1C3]
- 7) DNA Miniprep Clone CCA-46 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23110 + pSB1C3]
- 8) DNA Miniprep Clone CCA-47 Digestion XbaI/PstI [Phenanthrene insert 2 + Promoter BBa_J23110 + pSB1C3]
- 9) DNA Ladder (Molecular Weight Marker)



40 41 42 43 44 45 46 47 48 49 MW



5.4. CLONES DESCRIPTION

The clones that were kept are the following:

Clone CCA-23 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-26 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-29 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-42 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-44 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-46 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-30 [Fluorene insert 1 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-32 [Fluorene insert 1 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-34 [Fluorene insert 1 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-36 [Fluorene insert 2 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-38 [Fluorene insert 2 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-40 [Fluorene insert 2 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

6- CLONING OF FULL LENGTH (SAME PLASMID) POLYCISTRONIC F1, F2 (FLUORENE) AND P1, P2 (PHENANTHRENE) BEHIND CONSTITUTIVE PROMOTERS

6.1. PHENANTHRENE

6.1.1. Cloning of Phenanthrene Full Length (P1 and P2 on same plasmid)

6.1.1.1. Preparation of fragments for Phenanthrene Cloning

Source of clones:

The 2 clones CCA-23 and CCA-42 obtained as shown above were used to prepare a plasmid DNA containing the full length phenanthrene pathway under the control of promoter BBa_J23100

Miniprep DNA of CCC-23 was digested with SpeI/PstI and miniprep DNA of CCC-42 was digested with XbaI/PstI to release BBa-J23100_RBS_P2.

Clone CCA-23 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-42 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23100 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

The 2 clones CCA-26 and CCA-44 obtained as shown above were used to prepare a plasmid DNA containing the full length phenanthrene pathway under the control of promoter BBa_J23101

Miniprep DNA of CCC-26 was digested with SpeI/PstI and miniprep DNA of CCC-44 was digested with XbaI/PstI to release BBa-J23101_RBS_P2.

Clone CCA-26 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-44 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23101 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

The 2 clones CCA-29 and CCA-46 obtained as shown above were used to prepare a plasmid DNA containing the full length phenanthrene pathway under the control of promoter BBa_J23110.

Miniprep DNA of CCC-29 was digested with SpeI/PstI and miniprep DNA of CCC-46 was digested with XbaI/PstI to release BBa-J23110_RBS_P2.

Clone CCA-29 [Phenanthrene insert 1 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-46 [Phenanthrene insert 2 XbaI/PstI + Promoter BBa_J23110 EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

- a. The clones listed above were the sources of DNA for these steps of cloning.
- b. Double Digestion of recipient vector already containing P1 with SpeI and PstI
- c. Double Digestion of clones containing P2 with XbaI/PstI
- d. Restriction Digest Set-up (20µl reaction) in a 1.5 mL tube
- e. Add the reagents as described below.
- f. Incubate the tube at 37°C for 15 min
- g. Note: there is no need to add a loading buffer because the digestion buffer already has it.
- h. Load the reaction on a 1% Agarose gel, TAE
- i. After running the electrophoresis for 2 hours at 80V, cut the linearized band.
- j. Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit
- k. Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

P1 digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 µl	
10X Buffer	1.5 µl	1x
P1 Clones (CCA-23, CCA-26 and CCA-29)	2.0 µl	~ 0.3 µg
SpeI	1.0 µl	
PstI	1.0 µl	


P2 digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
P2 Clones (CCA-42, CCA-44 and CCA-46)	2.0 μ l	\sim 0.3 μ g
XbaI	1.0 μ l	
PstI	1.0 μ l	

Electrophoresis Gel 1

Preparation of recipient plasmid:

- 1) DNA Miniprep Clone CCA-23 Digestion SpeI/PstI [Phenanthrene insert 1 + Promoter are [BBa_J23100_RBS](#) + pSB1C3]
- 2) DNA Miniprep Clone CCA-26 Digestion SpeI/PstI [Phenanthrene insert 1 + Promoter [BBa_J23101_RBS](#) + pSB1C3]
- 3) DNA Miniprep Clone CCA-29 Digestion SpeI/PstI [Phenanthrene insert 1 + Promoter [BBa_J23110_RBS](#) + pSB1C3]
- 4) DNA Ladder (Molecular Weight Marker)
- 5) No gel loading



6.1.1.2. Ligation Phenanthrene Cloning

Date: 1-Aug-2017 and 2-Aug-2017

Materials

- a. T4 DNA Ligation Buffer. Invitrogen, Cat No. 46-0114
- b. T4 DNA Ligase, Thermo Scientific, Cat No. K1231
- c. Reagent grade water, NERL, Cat No. 98555
- d. 1.5 mL tube
- e. Vortex
- f. Pipet and tips
- g. Ice bucket and ice

Methods

- a. Set ligation as shown below in 1.5 mL tube.
- b. The tubes were incubated at room temperature for 1 hour.
- c. The tubes were then transferred to ice.

Ligation Condition (with insert)

Component	Volume (μ l)
Reagent grade water	6.0 μ l
10X T4 ligation buffer	1.0 μ l
SpeI/PstI Linearized recipient vector- pSB1C3-P1 under the control of promoters and RBS BBa_J23100 BBa_J23101 BBa_J23110	0.5 μ l
XbaI/PstI - fragment pSB1C3-P2 under the control of promoters and RBS: BBa_J23100 BBa_J23101 BBa_J23110	2.0 μ l
T4 DNA Ligase (5Weiss/ μ l)	0.5 μ l

Control Ligation Condition (no insert)

Component	Volume (μ l)
Reagent grade water	8.0 μ l
10X T4 ligation buffer	1.0 μ l
SpeI/PstI Linearized recipient vector- pSB1C3-P1 under the control of promoters and RBS BBa_J23100	0.5 μ l

BBa_J23101 BBa_J23110	
T4 DNA Ligase (5Weiss/ μ l)	0.5 μ l

6.1.2. Transformation

Date: 2-Aug-2017 and 3-Aug-2017

Materials

- a. LB Chloramphenicol 34 μ g/mL agar plates, Cat No. Teknova, L1017
- b. DH5a competent cells, Invitrogen, Cat 18265-017
- c. SOC (Recovery Medium), Lucigen, Cat No. F98226
- d. 1.5 mL tube
- e. Vortex
- f. Pipet and tips
- g. Ice bucket and ice
- h. Water bath (42°C)
- i. Incubator (37°C)

Method

- a. Transform ligated DNA into E.coli DH5a chemically competent cells
- b. Turn on incubator-shaker at 37°C.
- c. Turn on incubator for plates at 37°C.
- d. Set up water bath at 42°C.
- e. Bring to room temperature S.O.C medium.
- f. Bring LB plates supplemented with appropriate antibiotic at room temperature.
- g. Thaw competent cells on ice.
- h. Aliquots competent cells in as many tubes as needed.
- i. Add 5 μ l ligation mix to 50 μ l competent cells to DNA and swirl gently to mix.
- j. Incubate on ice for 20 minutes
- k. Heat shock at 42°C in heat block for 30 seconds. Quickly return to ice and let it sit for 2 min.
- l. Add 200 μ l of 18-25°C SOC medium and transfer mixture to 15mL Falcon tube
- m. Incubate in shaker at 37°C, 225 rpm for 30 min
- n. Plate 2 volumes (50 μ L and ~100 μ L) of the mixture onto 2 different plates of LB agar plates supplemented with Chloramphenicol 34 μ g/mL
- o. Incubate plates at 37°C overnight
- p. Count colonies and estimate transformation efficiency
- q. Initiate cultures in LB with the appropriate antibiotic from individual clones to analyze clones.

6.1.3. Transformation Results

Date: 2-Aug-2017 and 3-Aug-2017

Description [ligation]	Number of colonies (50 μ L volume plating)	Number of clones analyzed by digestion	Clones Designation
<u>Control:</u> Recipient vector containing P1_BBa_J23100_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI]	0		
<u>Ligation:</u> Recipient vector containing P1_BBa_J23100_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI] + Synthetic phnAc, phnAd, phnB [XbaI/PstI]	>100	3	CCA-57 CCA-58 CCA-59
<u>Control:</u> Recipient vector containing P1_BBa_J23101_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI]	0		
<u>Ligation:</u> Recipient vector containing P1_BBa_J23101_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI] + Synthetic phnAc, phnAd, phnB [XbaI/PstI]	>100	3	CCA-60 CCA-61 CCA-62
<u>Control:</u> Recipient vector containing P1_BBa_J23110_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI]	0		
<u>Ligation:</u> Recipient vector containing P1_BBa_J23110_RBS_Synthetic phnF, phnE, phnC, phnD_ pSB1C3 [SpeI/PstI] + Synthetic phnAc, phnAd, phnB [XbaI/PstI]	>100	3	CCA-63 CCA-64 CCA-65




Figure. LB-Chloramphenicol plates with E.coli DH5a cells transformed with product of ligation of pSBB1C3 containing Phenanthrene-1 (SpeI/PstI) ligated to Phenanthrene-2 (XbaI/PstI) under the control of the constitutive promoter [BBa_J23100](#) [100] or [BBa_J23101](#) [101] or [BBa_J23110](#)[110] and RBS.

Phenanthrene-1= P1: Synthetic phnF, phnE, phnC, phnD

Phenanthrene-2= P2: Synthetic phnAc, phnAd, phnB

6.1.4. Clone Verification

Date: 5-Aug-2017/6-Aug-2017

Clone verification is done by digestion. Cultures are set up to extract DNA.

Materials

- DNA extraction kit, Zymo Research
- LB liquid medium
- Antibiotic stock solution

- d. Vortex
- e. Pipet and tips
- f. 15 mL culture tube
- g. Incubator-Shaker
- h. 10 mL pipette
- i. Pipet aid
- j. Vortex
- k. Rack
- l. Toothpick

Set-up Culture

- a. Grow selected number of colonies in 3 mL LB medium supplemented with chloramphenicol (34 $\mu\text{g}/\text{mL}$) overnight at 37°C in the incubator/shaker at 220 rpm.
- b. Spin down 1.5 mL of culture from each culture and isolate plasmid DNA using the Zymo Research DNA Isolation kit.
- c. DNA preparation is resuspended in a final volume of **35 μL** reagent grade water.
- d. Store remaining 1mL of culture for glycerol stock preparation of sequence confirmed clones.
- e. Check clones by digestion.

Clones

Analysis of Individual clones: CCA_57; CCA_58; CCA_59

BBa_J23100_RBS_P1_Terminator + BBa_J23100_RBS_P2_Terminator in vector pSB1C3

Analysis of Individual clones: CCA_60; CCA_61; CCA_62

BBa_J23101_RBS_P1_Terminator + BBa_J23101_RBS_P2_Terminator in vector pSB1C3

Analysis of Individual clones: CCA_63; CCA_64; CCA_65

BBa_J23110_RBS_P1_Terminator + BBa_J23110_RBS_P2_Terminator in vector pSB1C3

P1: Synthetic phnF, phnE, phnC, phnD_Terminator

P2: Synthetic phnAc, phnAd, phnB-Terminator

Set-up Digestion for Clone Verification

Date: 5-Aug-2017

- a. Set up digestion as shown below.
- b. This is a 1% gel agarose gel, TAE, prepared with Sybr green dye.
- c. Digest clones with EcoRI/PstI
- d. Set-up restriction digestion (15 μl reaction) in a 1.5 mL tube
- e. Turn on water bath at 37°C
- f. Add the reagents in the order and with volume described in the table below.
- g. Spin the tube briefly for 15 seconds at 10,000 rpm
- h. Incubate the tube at 37°C for 15 min
- i. The reaction already have a loading buffer

- j. Load the digestion reaction on a 1% Agarose gel, TAE
- k. Set up the power for electrophoresis to 120 V and let it run for ~ 1 hour.

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X BufferFastDigest Green Buffer	1.5 μ l	1 X
Plasmid miniprep	1.5 μ l	~0.2-0.4 μ g
Restriction Enzyme 1	1.0 μ l	
Restriction Enzyme 2	1.0 μ l	

Expected Fragments Size for Verification

Phenanthrene clones:

EcoRI/PstI digestion: 2 fragments should be observed: the vector pSB1C3 and P1+P2 (4165 bp+3057bp)

P1_Terminator: XbaI/Pst1: 4115 bp

Promoter_RBS_P1_Terminator: 4183 bp

Promoter_RBS_P1_Terminator: 4165 bp EcoRI/SpeI

P2_Terminator: XbaI/Pst1: 3022 bp

Promoter_RBS_P2_Terminator: 3090 bp

Promoter_RBS_P2_Terminator: 3057 bp EcoRI/SpeI

6.2. FLUORENE

6.2.1. Cloning of Fluorene Full Length (F1 and F2) into a Vector with p15a Ori

6.2.1.1. Preparation of fragments for Fluorene Cloning

Source of clones:

The 2 clones **CCA-30** and **CCA-36** obtained as shown above were used to prepare a plasmid DNA containing the full length fluorene pathway under the control of promoter **BBa_J23100**.

Miniprep DNA of CCC-30 was digested with EcoRI/SpeI to release BBa-J23100_RBS_F1 and miniprep DNA of CCC-36 and was digested with XbaI/PstI to release BBa-J23100_RBS_F2.

Clone CCA-30 [Fluorene insert 1 XbaI/PstI + Promoter **BBa_J23100** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-36 [Fluorene insert 2 XbaI/PstI + Promoter **BBa_J23100** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

The 2 clones **CCA-32** and **CCA-38** obtained as shown above were used to prepare a plasmid DNA containing the full length fluorene pathway under the control of promoter **BBa_J23101**.

Miniprep DNA of CCC-32 was digested with EcoRI/SpeI to release BBa-J23101_RBS_F1 and miniprep DNA of CCC-38 and was digested with XbaI/PstI to release BBa-J23101_RBS_F2.

Clone CCA-32 [Fluorene insert 1 XbaI/PstI + Promoter **BBa_J23101** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-38 [Fluorene insert 2 XbaI/PstI + Promoter **BBa_J23101** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

The 2 clones **CCA-34** and **CCA-40** obtained as shown above were used to prepare a plasmid DNA containing the full length fluorene pathway under the control of promoter **BBa_J23110**.

Miniprep DNA of CCC-34 was digested with EcoRI/SpeI to release BBa-J23110_RBS_F1 and miniprep DNA of CCC-38 and was digested with XbaI/PstI to release BBa-J23110_RBS_F2.

Clone CCA-34 [Fluorene insert 1 XbaI/PstI + Promoter **BBa_J23110** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Clone CCA-40 [Fluorene insert 2 XbaI/PstI + Promoter **BBa_J23110** EcoRI/SpeI+ pSB1C3 EcoRI/PstI]

Plasmid pSB3T5 (available with part BBa_J04450 in 2017 distribution plate) which contain p15a as origin of replication was digested with EcoRI and PstI.

Methods:

- Double Digestion of F1 under the control of one of the 3 promoters with EcoRI/SpeI
- Double Digestion of F2 under the control of one of the 3 promoters with XbaI/PstI
- Double Digestion of recipient vector with EcoRI/PstI
- Restriction Digest Set-up (20µl reaction) in a 1.5 mL tube
- Add the reagents as described below.
- Incubate the tube at 37°C for 15 min
- Note: there is no need to add a loading buffer because the digestion buffer already has it.
- Load the reaction on a 1% Agarose gel
- After running the electrophoresis for 2 hours at 80V, cut the linearized band.
- Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit
- Elute the DNA from the Zymoclean column with 12 µl of sterile filtered H₂O

F1 digestion condition:

Component	Volume	Final condition

Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
F1 Clones (CCA-30, CCA-32 and CCA-34)	2.0 μ l	~ 0.3 μ g
EcoRI	1.0 μ l	
SpeI	1.0 μ l	

F2 digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
F1 Clones (CCA-36, CCA-38 and CCA-40)	2.0 μ l	~ 0.3 μ g
XbaI	1.0 μ l	
PstI	1.0 μ l	

Vector digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
Vector	2.0 μ l	~ 0.3 μ g
EcoRI	1.0 μ l	
PstI	1.0 μ l	

Electrophoresis Gel 1


- 4) DNA Ladder (Molecular Weight Marker)
- 5) No gel loading

Preparation of vector in duplicate:

- 6) pSB3T5 (p15 a) vector EcoRI/PstI

7) pSB3T5 (p15 a) vector EcoRI/PstI

MW567



6.2.1.2. Ligation Fluorene Cloning

Date: 1-Aug-2017 and 2-Aug-2017

Materials

- T4 DNA Ligation Buffer. Invitrogen, Cat No. 46-0114
- T4 DNA Ligase, Thermo Scientific, Cat No. K1231
- Reagent grade water, NERL, Cat No. 98555
- 1.5 mL tube
- Vortex
- Pipet and tips
- Ice bucket and ice

Methods

- Set ligation as shown below in 1.5 mL tube.
- The tubes were incubated at room temperature for 1 hour.
- The tubes were then transferred to ice.

Ligation Condition (with insert)

Component	Volume (µl)
Reagent grade water	5.0 µl
10X T4 ligation buffer	1.0 µl
Linearized Recipient Vector EcoRI/PstI (Tetracycline)	0.5 µl
EcoRI/SpeI - F1 insert under the control of promoters and RBS	1.5 µl

BBa_J23100 BBa_J23101 BBa_J23110	
XbaI/PstI – F2 insert under the control of promoters and RBS: BBa_J23100 BBa_J23101 BBa_J23110	1.5µl
T4 DNA Ligase (5Weiss/µl)	0.5 µl

Control Ligation Condition (no insert)

Component	Volume (µl)
Reagent grade water	8.0 µl
10X T4 ligation buffer	1.0 µl
Linearized Recipient Vector EcoRI/PstI (Tetracycline)	0.5 µl
T4 DNA Ligase (5Weiss/µl)	0.5µl

6.2.1.3. Transformation

Date: 4-Aug-2017

Materials

- a. LB Tetracycline 12.5 µg/mL agar plates, Cat No. Teknova, L5072
- b. DH5a competent cells, Invitrogen, Cat 18265-017
- c. SOC (Recovery Medium), Lucigen, Cat No. F98226
- d. 1.5 mL tube
- e. Vortex
- f. Pipet and tips
- g. Ice bucket and ice
- h. Water bath (42°C)
- i. Incubator (37°C)

Method

- a. Transform ligated DNA into E.coli DH5a chemically competent cells
- b. Turn on incubator-shaker at 37°C.
- c. Turn on incubator for plates at 37°C.
- d. Set up water bath at 42°C.
- e. Bring to room temperature S.O.C medium.
- f. Bring LB plates supplemented with appropriate antibiotic at room temperature.
- g. Thaw competent cells on ice.
- h. Aliquots competent cells in as many tubes as needed.

- i. Add 5 μ l ligation mix to 50 μ l competent cells to DNA and swirl gently to mix.
- j. Incubate on ice for 20 minutes
- k. Heat shock at 42°C in heat block for 30 seconds. Quickly return to ice and let it sit for 2 min.
- l. Add 200 μ l of 18-25°C SOC medium and transfer mixture to 15mL Falcon tube
- m. Incubate in shaker at 37°C, 225 rpm for 30 min
- n. Plate 2 volumes (50 μ L and \sim 100 μ L) of the mixture onto 2 different plates of LB agar plates supplemented with tetracycline 12.5 μ g/mL
- o. Incubate plates at 37°C overnight
- p. Count colonies and estimate transformation efficiency
- q. Initiate cultures in LB with the appropriate antibiotic from individual clones to analyze clones

6.2.1.4. Transformation Results

Date: 5-Aug-2017

Description [ligation]	Number of colonies (50 μ L volume plating)	Number of clones analyzed by digestion	Clones Designation
<u>Control</u> : Recipient vector [EcoRI/PstI]	0		
<u>Ligation</u> : Recipient vector p15a + F1 +F2 each under control of promoter BBa_J23100	>100	3	CCA-48 CCA-49 CCA-50
<u>Control</u> : Recipient vector [EcoRI/PstI]	0		
<u>Ligation</u> : Recipient vector p15a + F1 +F2 each under control of promoter BBa_J23101	>100	3	CCA-51 CCA-52 CCA-53
<u>Control</u> : Recipient vector [EcoRI/PstI]	0		
<u>Ligation</u> : Recipient vector p15a + F1 +F2 each under control of promoter BBa_J23110	>100	3	CCA-54 CCA-55 CCA-56




Figure. LB-Tetracycline plates with E.coli DH5a cells transformed with product of ligation of p15aTet (vector pSB3T5) (EcoRI/PstI) ligated to Fluorene-1 (EcoRI/SpeI) and ligated to Fluorene-2 (XbaI/PstI) under the control of the constitutive promoter [BBa_J23100](#) [100] or [BBa_J23101](#) [101] or [BBa_J23110](#)[110].

Full length Synthetic fluorene catabolic pathway = Synthetic flnB, dbfA1, dbfA2+ Synthetic flnE, flnD1, ORF16, flnC

Fluorene-1= Promoter_RBS_Synthetic flnB, dbfA1, dbfA2
 Fluorene-2= Promoter_RBS_Synthetic flnE, flnD1, ORF16, flnC

6.2.1.5. Clone Verification

Date: 6-Aug-2017/7-Aug-2017

Clone verification is done by digestion. Cultures are set up to extract DNA.

Materials

- DNA extraction kit, Zymo Research
- LB liquid medium
- Antibiotic stock solution
- Vortex
- Pipet and tips
- 15 mL culture tube
- Incubator-Shaker
- 10 mL pipette
- Pipet aid
- Vortex
- Rack
- Toothpick

Set-up Culture

- Grow selected number of colonies in 3 mL LB medium supplemented with tetracycline (12.5µg/mL) overnight at 37°C in the incubator/shaker at 220 rpm.

- b. Spin down 1.5 mL of culture from each culture and isolate plasmid DNA using the Zymo Research DNA Isolation kit.
- c. DNA preparation is resuspended in a final volume of **35 µL** reagent grade water.
- d. Store remaining 1mL of culture for glycerol stock preparation of sequence confirmed clones.
- e. Check clones by digestion.

Clones

Analysis of Individual clones: CCA_48; CCA_49; CCA_50

BBa_J23100_RBS_F1_Terminator + BBa_J23100_RBS_F2_Terminator in vector pSB3T5

Analysis of Individual clones: CCA_51; CCA_52; CCA_53

BBa_J23101_RBS_F1_Terminator + BBa_J23101_RBS_F2_Terminator in vector pSB3T5

Analysis of Individual clones: CCA_54; CCA_55; CCA_56

BBa_J23110_RBS_F1_Terminator + BBa_J23110_RBS_F2_Terminator in vector pSB3T5

F1: Synthetic flnB, dbfA1, dbfA2

F2: Synthetic flnE, flnD1, ORF16, flnC

Set-up Digestion for Clone Verification

Date: 7-Aug-2017

- a. Set up digestion as shown below.
- b. This is a 1% gel agarose gel, TAE, prepared with Sybr green dye.
- c. Digest clones with EcorRI/PstI
- d. Set-up restriction digestion (20µl reaction) in a 1.5 mL tube
- e. Turn on water bath at 37°C
- f. Add the reagents in the order and with volume described in the table below.
- g. Spin the tube briefly for 15 seconds at 10,000 rpm
- h. Incubate the tube at 37°C for 15 min
- i. The reaction already have a loading buffer
- j. Load the digestion reaction on a 1% Agarose gel, TAE
- k. Set up the power for electrophoresis to 120 V and let it run for ~ 1 hour.

Component	Volume	Final condition
Reagent grade water	15.0 µl	
10X BufferFastDigest Green Buffer	1.5 µl	1 X
Plasmid miniprep	1.5 µl	~0.2-0.4 µg
Restriction Enzyme 1	1.0 µl	
Restriction Enzyme 2	1.0 µl	

Expected Fragments Size for Verification

Fluorene: EcoRI/PstI digestion:

2 fragments should be observed: the vector pSB3T5 and F1+F2 (3157bp+3468 bp)

F1_Terminator: XbaI/PstI: 3107 bp

Promoter_RBS_F1_Terminator: 3175 bp

Promoter_RBS_F1_Terminator: 3157 bp EcoRI/Spel

F2_Terminator: XbaI/PstI: 3433 bp

Promoter_RBS_F_Terminator2: 3501 bp


Promoter_RBS_F2_Terminator: 3468 bp EcoRI/Spel

Results

The DNA digestion of clones CCA-48 to CCA-58 indicated that the fragments of the correct size were obtained.

Clone CCA-48, CCA51, and CCA54 will be sent out for sequencing for verification.

CCA-48 49 50 5152 5354 55 56 MW

**7- SEQUENCING**

Date: 8-Aug-2017

An aliquot of the DNA plasmid mini-prep was sent out for sequencing.

All files are attached below.

The sequencing data indicated that the sequence of the clones was correct.

8- TRANSFORMATION OF FULL LENGTH PHENANTHRENE AND FLUORENE CLONES INTO E.COLI BL-21

Date: 10-Aug-2017

8.1. MATERIALS

- a. DNA miniprep
- b. LB agar plates, Cat No. Teknova, appropriate antibiotic
- c. E coli **BL21 DE3**, Life Technology, Cat No. 60106-1
- d. SOC (Recovery Medium), Lucigen, Cat No. F98226
- e. 15 mL culture tube
- f. Vortex
- g. Tooth pick
- h. Incubator shaker
- i. 42°C Water bath
- j. Ice and ice bucket
- k. Pipet

8.2. METHODS

- a. Transform 1 μL of DNA mini-preparation into E coli BL21 DE3 chemically competent cells
- b. For double transformation, transform with of 1 μL of each DNA mini-preparation into E coli BL21 DE3 chemically competent cells
- c. Turn on incubator-shaker at 37°C.
- d. Turn on incubator for plates at 37°C.
- e. Set up water bath at 42°C.
- f. Bring to room temperature S.O.C medium.
- g. Bring LB plates supplemented with appropriate antibiotic at room temperature.
- h. Thaw competent cells on ice.
- i. Aliquots competent cells in as many tubes as needed.
- j. Add 1.0 μl DNA preparation to 40 μl competent cells to DNA and swirl gently to mix
- k. Incubate on ice for 20 minutes
- l. Heat shock at 42°C in heat block for 30 seconds. Quickly return to ice and let it sit for 2 min.
- m. Add 200 μl of 18-25°C SOC medium and transfer mixture to 15mL Falcon tube
- n. Incubate in shaker at 37°C, 225 rpm for 30 min
- o. Plate 1 volumes ($\sim 100\mu\text{L}$) of the mixture onto one LB agar plates supplemented with Chloramphenicol 35 $\mu\text{g}/\text{mL}$ or Tetracycline 12.5 $\mu\text{g}/\text{mL}$ or both Chloramphenicol 35 $\mu\text{g}/\text{mL}$ and Tetracycline 12.5 $\mu\text{g}/\text{mL}$ depending on the clone
- p. Incubate plates at 37°C overnight
- q. Initiate cultures in LB with the appropriate antibiotic from individual clones to analyze clones by digestion




Figure. LB-Chloramphenicol plates with E.coli BL-21 cells transformed with product mini-prep DNA of clone CCA-57, CCA 60, and CCA 64.

Clone_57: BBa_J23100_RBS_P1_Terminator +
BBa_J23100_RBS_P2_Terminator in vector pSB1C3

CCA_60: BBa_J23101_RBS_P1_Terminator +
BBa_J23101_RBS_P2_Terminator in vector pSB1C3

CCA_64: BBa_J23110_RBS_P1_Terminator +
BBa_J23110_RBS_P2_Terminator in vector pSB1C3

Phenanthrene-1= Synthetic phnF, phnE, phnC, phnD
Phenanthrene-2= Synthetic phnAc, phnAd, phnB




Figure. LB-Tetracycline plates with E.coli BL-21 cells transformed with product mini-prep DNA of clone CCA-48, CCA 51, and CCA 54.

CCA_48=BBa_J23100_RBS_F1_Terminator +
BBa_J23100_RBS_F2_Terminator in vector pSB3T5

CCA_51=BBa_J23101_RBS_F1_Terminator +
BBa_J23101_RBS_F2_Terminator in vector pSB3T5

CCA_54=BBa_J23110_RBS_F1_Terminator +
BBa_J23110_RBS_F2_Terminator in vector pSB3T5

Fluorene-1= Promoter_RBS_Synthetic flnB, dbfA1, dbfA2
Fluorene-2= Promoter_RBS_Synthetic flnE, flnD1,
ORF16, flnC

9- CLONING OF FLUORENE FULL LENGTH F1 AND F2 INTO pSB1C3 FOR DEPOSIT IN REGISTRY

9.1. PREPARATION OF DNA FRAGMENTS FOR CLONING OF FLUORENE F1 AND F2 INTO pSB1C3 FOR DEPOSIT IN REGISTRY

Source of clones:

DNA of CCA-48, CCA-51, and CCA 54 containing the full fluorene pathway in the pSB3T5 vector are the source of materials.

Methods:

- Double Digestion of CCA-48, CCA-51, and CCA 54 under the control of one of the 3 promoters with EcoRI/PstI
- Double Digestion of recipient vector pSB1C3 with EcoRI/PstI
- Restriction Digest Set-up (20 μ l reaction) in a 1.5 mL tube
- Add the reagents as described below.
- Incubate the tube at 37°C for 15 min
- Note: there is no need to add a loading buffer because the digestion buffer already has it.
- Load the reaction on a 1% Agarose gel
- After running the electrophoresis for 2 hours at 80V, cut the linearized band.
- Purify the DNA from the gel using the Zymoclean Gel DNA Recovery Kit
- Elute the DNA from the Zymoclean column with 12 μ l of sterile filtered H₂O

digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
CCA-48, CCA-51, and CCA 54	2.0 μ l	~ 0.3 μ g
EcoRI	1.0 μ l	
PstI	1.0 μ l	

Vector digestion condition:

Component	Volume	Final condition
Reagent grade water	9.5 μ l	
10X Buffer	1.5 μ l	1x
Vector	2.0 μ l	~ 0.3 μ g

EcoRI	1.0 μ l	
PstI	1.0 μ l	


Electrophoresis Gel 1

- 4) DNA Ladder (Molecular Weight Marker)
- 5) No gel loading

Preparation of vector in duplicate:

- 6) pSB3T5 (p15 a) vector EcoRI/PstI
- 7) pSB3T5 (p15 a) vector EcoRI/PstI

MW567



9.2. LIGATION OF DNA FRAGMENTS F1 AND F2 INTO pSB1C3

Date: 1-Sep-2017 and 2-Sep-2017

Materials

- h. T4 DNA Ligation Buffer. Invitrogen, Cat No. 46-0114
- i. T4 DNA Ligase, Thermo Scientific, Cat No. K1231
- j. Reagent grade water, NERL, Cat No. 98555
- k. 1.5 mL tube
- l. Vortex
- m. Pipet and tips
- n. Ice bucket and ice

Methods

- d. Set ligation as shown below in 1.5 mL tube.
- e. The tubes were incubated at room temperature for 1 hour.

- f. The tubes were then transferred to ice.

Ligation Condition (with insert)

Component	Volume (μ l)
Reagent grade water	5.5 μ l
10X T4 ligation buffer	1.0 μ l
Linearized Recipient Vector EcoRI/PstI (Chloramphenicol)	0.5 μ l
Linearized insert	2.0 μ l
T4 DNA Ligase (5Weiss/ μ l)	1 μ l

Control Ligation Condition (no insert)

Component	Volume (μ l)
Reagent grade water	8.0 μ l
10X T4 ligation buffer	1.0 μ l
Linearized Recipient Vector EcoRI/PstI (Chloramphenicol)	0.5 μ l
T4 DNA Ligase (5Weiss/ μ l)	0.5 μ l

9.3. TRANSFORMATION

Materials

- j. LB Chloramphenicol 34 μ g/mL agar plates, Teknova
- k. DH5a competent cells, Invitrogen, Cat 18265-017
- l. SOC (Recovery Medium), Lucigen, Cat No. F98226
- m. 1.5 mL tube
- n. Vortex
- o. Pipet and tips
- p. Ice bucket and ice
- q. Water bath (42°C)
- r. Incubator (37°C)

Method

- r. Transform ligated DNA into E.coli DH5a chemically competent cells
- s. Turn on incubator-shaker at 37°C.
- t. Turn on incubator for plates at 37°C.
- u. Set up water bath at 42°C.
- v. Bring to room temperature S.O.C medium.
- w. Bring LB plates supplemented with appropriate antibiotic at room temperature.

- x. Thaw competent cells on ice.
- y. Aliquots competent cells in as many tubes as needed.
- z. Add 5 μ l ligation mix to 50 μ l competent cells to DNA and swirl gently to mix.
- aa. Incubate on ice for 20 minutes
- bb. Heat shock at 42°C in heat block for 30 seconds. Quickly return to ice and let it sit for 2 min.
- cc. Add 200 μ l of 18-25°C SOC medium and transfer mixture to 15mL Falcon tube
- dd. Incubate in shaker at 37°C, 225 rpm for 30 min
- ee. Plate 2 volumes (50 μ L and ~100 μ L) of the mixture onto 2 different plates of LB agar plates supplemented with tetracycline 12.5 μ g/mL
- ff. Incubate plates at 37°C overnight
- gg. Count colonies and estimate transformation efficiency
- hh. Initiate cultures in LB with the appropriate antibiotic from individual clones to analyze clones

9.4. TRANSFORMATION RESULTS

Description [ligation]	Number of colonies (50 μ L volume plating)	Number of clones analyzed by digestion	Clones Designation
<u>Control</u> : Recipient vector [EcoRI/PstI]	0		
<u>Ligation</u> : Recipient vector pSB1C3+ F1_F2 each under control of promoter BBa_J23100	>100	2	CCA-66 CCA-67
<u>Ligation</u> : Recipient vector pSB1C3+ F1_F2 each under control of promoter BBa_J23101	>100	2	CCA-68 CCA-69
<u>Ligation</u> : Recipient vector pSB1C3+ F1_F2 each under control of promoter BBa_J23110	>100	2	CCA-70 CCA-71




Figure. LB-Chloramphenicol plates with E.coli DH5a cells transformed with product of ligation of pSB1C3 ligated to Fluorene-1 + Fluorene-2 (XbaI/PstI) under the control of the constitutive promoter BBa_J23100 [100] or BBa_J23101 [101] or BBa_J23110[110].

Full length Synthetic fluorene catabolic pathway = Synthetic flnB, dbfA1, dbfA2+ Synthetic flnE, flnD1, ORF16, flnC

Fluorene-1= Promoter_RBS_Synthetic flnB, dbfA1, dbfA2
 Fluorene-2= Promoter_RBS_Synthetic flnE, flnD1, ORF16, flnC

9.5. CLONE VERIFICATION

Clone verification is done by digestion. Cultures are set up to extract DNA.

Materials

- DNA extraction kit, Zymo Research
- LB liquid medium
- Antibiotic stock solution
- Vortex
- Pipet and tips
- 15 mL culture tube
- Incubator-Shaker
- 10 mL pipette
- Pipet aid
- Vortex
- Rack
- Toothpick

Set-up Culture

- f. Grow selected number of colonies in 3 mL LB medium supplemented with tetracycline (12.5µg/mL) overnight at 37°C in the incubator/shaker at 220 rpm.
- g. Spin down 1.5 mL of culture from each culture and isolate plasmid DNA using the Zymo Research DNA Isolation kit.
- h. DNA preparation is resuspended in a final volume of **35 µL** reagent grade water.
- i. Store remaining 1mL of culture for glycerol stock preparation of sequence confirmed clones.
- j. Check clones by digestion.

Clones

Analysis of Individual clones: CCA_66 and CCA_67

BBa_J23100_RBS_F1_Terminator + BBa_J23100_RBS_F2_Terminator in vector pSB1C3

Analysis of Individual clones: CCA_68 and CCA_69

BBa_J23101_RBS_F1_Terminator + BBa_J23101_RBS_F2_Terminator in vector pSB1C3

Analysis of Individual clones: CCA_70 and CCA_71

BBa_J23110_RBS_F1_Terminator + BBa_J23110_RBS_F2_Terminator in vector pSB1C3

F1: Synthetic flnB, dbfA1, dbfA2

F2: Synthetic flnE, flnD1, ORF16, flnC

Set-up Digestion for Clone Verification

Date: 7-Aug-2017

- l. Set up digestion as shown below.
- m. This is a 1% gel agarose gel, TAE, prepared with Sybr green dye.
- n. Digest clones with EcorRI/PstI
- o. Set-up restriction digestion (20µl reaction) in a 1.5 mL tube
- p. Turn on water bath at 37°C
- q. Add the reagents in the order and with volume described in the table below.
- r. Spin the tube briefly for 15 seconds at 10,000 rpm
- s. Incubate the tube at 37°C for 15 min
- t. The reaction already have a loading buffer
- u. Load the digestion reaction on a 1% Agarose gel, TAE
- v. Set up the power for electrophoresis to 120 V and let it run for ~ 1 hour.

Component	Volume	Final condition
Reagent grade water	15.0 µl	
10X BufferFastDigest Green Buffer	1.5 µl	1 X
Plasmid miniprep	1.5 µl	~0.2-0.4 µg
Restriction Enzyme 1	1.0 µl	

Restriction Enzyme 2	1.0 μ l	
----------------------	-------------	--

Expected Fragments Size for Verification

Fluorene: EcoRI/PstI digestion:

2 fragments should be observed: the vector pSB1C3 and F1+F2 (3157bp+3468 bp)

F1_Terminator: XbaI/PstI: 3107 bp

Promoter_RBS_F1_Terminator: 3175 bp

Promoter_RBS_F1_Terminator: 3157 bp EcoRI/SpeI

F2_Terminator: XbaI/PstI: 3433 bp

Promoter_RBS_F_Terminator2: 3501 bp


Promoter_RBS_F2_Terminator: 3468 bp EcoRI/SpeI

Results

Clones are ok except for CCA-71.

Clone CCA-66, CCA68, and CCA70 will be sent out for sequencing for verification.

CCA-66 67 68 69 70 71



9.6. VERIFICATION BY SEQUENCING OF FULL LENGTH CLONE F1_F2

An aliquot of the DNA plasmid mini-prep was sent out for sequencing.

All files are attached below.

The data showed that the sequence of the clones was correct.

10- GLYCEROL STOCK

Date: 14-AUG-2017

Materials

- a. Glycerol
- b. LB liquid medium
- c. Antibiotic stock solution
- d. 15 mL culture tube
- e. Incubator-Shaker
- f. 10 mL pipette
- g. Pipet aid
- h. Vortex
- i. Rack
- j. Toothpick

Set-up Culture and Prepare Glycerol Stocks

- a. Grow selected clones that have been checked (so they are correct) in 3 mL LB medium supplemented with appropriate antibiotics overnight at 37°C in the incubator/shaker at 220 rpm.
- b. Prepare LB medium with 40% glycerol and add 0.5 mL to a cryogenic vial
- c. Add 0.5 mL of culture sample to be stored
- d. Gently vortex the cryogenic vial to ensure the culture and glycerol is well mixed
- e. Label tube with date and identifier
- f. Organize in a freezer box and label box
- g. Prepare excel spreadsheet with all information
- h. Store freezer box at -80°C

Glycerol Stock	Date of Glycerol Stock	Clone	Description	Cell Description	Selection Marker	Vector
CCA-1000	29-Jul-17	CCA-1	Full Synthetic catabolic pathway Fluorene Inducible promoter	E.coli DH5a	Ampicillin	pUC57
CCA-1001	29-Jul-17	CCA-12	Full Synthetic catabolic pathway Phenanthrene Inducible promoter	E.coli DH5a	Ampicillin	pUC57
CCA-1002	29-Jul-17	CCA-1	Full Synthetic catabolic pathway Fluorene Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1003	29-Jul-17	CCA-12	Full Synthetic catabolic pathway Phenanthrene Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1004	29-Jul-17	Synthetic Fluorene-1	Synthetic flnB, dbfA1, dbfA2 Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1005	29-Jul-17	Synthetic Fluorene-2	Synthetic flnE, flnD1, ORF16, flnC Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1006	29-Jul-17	Synthetic Phenanthrene 1	Synthetic phnF, phnE, phnC, phnD Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1007	29-Jul-17	Synthetic Phenanthrene 2	Synthetic phnAc, phnAd, phnB Inducible promoter	E.coli BL-21	Ampicillin	pUC57
CCA-1008	29-Jul-17	pUC19		E.coli BL-21	Ampicillin	pUC19
CCA-1009	14-Aug-17	CCA-23	Promoter BBa_J23100 /RBS BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI + pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1010	14-Aug-17	CCA-26	Promoter BBa_J23101 /RBS BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI + pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1011	14-Aug-17	CCA-29	Promoter BBa_J23110 /RBS BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI + pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1012	14-Aug-17	CCA-30	Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100 BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1013	14-Aug-17	CCA-32	Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter 14-Aug-17 BBa_J23101 BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1014	14-Aug-17	CCA-34	Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110 BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1015	14-Aug-17	CCA-36	Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100 BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1035	14-Aug-17	CCA-38	Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter	E.coli DH5a	Chloramphenicol	pSB1C3

			BBa_J23101_BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI			
CCA-1016	14-Aug-17	CCA-40	Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1017	14-Aug-17	CCA-42	Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1018	14-Aug-17	CCA-44	Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1019	14-Aug-17	CCA-46	Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI+ pSB1C3 EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1020	14-Aug-17	CCA-48	CCA-30 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-36 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Tetracycline	pSB3T5
CCA-1021	14-Aug-17	CCA-51	CCA-30 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-38 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Tetracycline	pSB3T5
CCA-1022	14-Aug-17	CCA-54	CCA-30 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-40 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Tetracycline	pSB3T5
CCA-1023	14-Aug-17	CCA-57	CCA-23 [Promoter BBa_J23100 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-42 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1024	14-Aug-17	CCA-60	CCA-26 [Promoter BBa_J23101 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-44 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1025	14-Aug-17	CCA-64	CCA-29 [Promoter BBa_J23110 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-46 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1026	14-Aug-17	CCA-48	CCA-30 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter	E.coli BL-21	Tetracycline	pSB3T5

			BBa_J23100_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-36 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI			
CCA-1027	14-Aug-17	CCA-51	CCA-32 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-38 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli BL-21	Tetracycline	pSB3T5
CCA-1028	14-Aug-17	CCA-54	CCA-34 [Fluorene insert 1_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-40 [Fluorene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli BL-21	Tetracycline	pSB3T5
CCA-1029	14-Aug-17	CCA-57	CCA-23 [Promoter BBa_J23100 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-42 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23100_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli BL-21	Chloramphenicol	pSB1C3
CCA-1030	14-Aug-17	CCA-60	CCA-26 [Promoter BBa_J23101 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-44 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23101_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli BL-21	Chloramphenicol	pSB1C3
CCA-1031	14-Aug-17	CCA-64	CCA-29 [Promoter BBa_J23110 /RBS_BBa_B0034 EcoRI/SpeI + Phenanthrene insert 1_Ter_BBa_B0015 XbaI/PstI] as SpeI/PstI CCA-46 [Phenanthrene insert 2_Ter_BBa_B0015 XbaI/PstI + Promoter BBa_J23110_BBa_B0034 EcoRI/SpeI] as XbaI/PstI	E.coli BL-21	Chloramphenicol	pSB1C3
CCA-1032	14-Aug-17	CCA-48 + CCA-57		E.coli BL-21	Chloramphenicol, Tetracycline	pSB1C3, pSB3T5
CCA-1033	14-Aug-17	CCA-51 + CCA-60		E.coli BL-21	Chloramphenicol, Tetracycline	pSB1C3, pSB3T5
CCA-1034	14-Aug-17	CCA-54 + CCA-64		E.coli BL-21	Chloramphenicol, Tetracycline	pSB1C3, pSB3T5
CCA-1036	14-Aug-17			E.coli BL-21	Tetracycline	pSB3T5

CCA-1037				E.coli BL-21	Chloramphenicol	pSB1C3
CCA-1038	14-Aug-17			E.coli BL-21	Chloramphenicol, Tetracycline	pSB1C3, pSB3T5
CCA-1039	1-Sep-17	CCA-66	CCA-30 [Fluorene insert 1_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23100_BB_a_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-36 [Fluorene insert 2_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23100_BB_a_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1040	1-Sep-17	CCA-68	CCA-32 [Fluorene insert 1_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23101_BB_a_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-38 [Fluorene insert 2_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23101_BB_a_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3
CCA-1041	1-Sep-17	CCA-70	CCA-34 [Fluorene insert 1_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23110_BB_a_B0034 EcoRI/SpeI] as EcoRI/SpeI CCA-40 [Fluorene insert 2_Ter_BB_a_B0015 XbaI/PstI + Promoter BBa_J23110_BB_a_B0034 EcoRI/SpeI] as XbaI/PstI Vector pSB3T5 as EcoRI/PstI	E.coli DH5a	Chloramphenicol	pSB1C3

11- SEQUENCING FILES

Primers:

Name	Description	Sequence	%GCTm	Length
BBa_G00100	Forward primer for sequencing/amplifying BioBrick parts (VF2)	tgccacctgacgtctaagaa50	60C	20
BBa_G00101	Reverse primer for sequencing/amplifying BioBrick parts (VR)	attaccgcctttgagtgagc50	60C	20

Sequence verification

Fluorene clones CCA-48, CCA-51, and CCA-54.

Sequence Alignment of clones CCA-58, CCA-51, and CCA-54.

CLUSTAL O(1.2.4) multiple sequence alignment

Forward sequence

```

54      TGCGAGAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGT
48      --GAAGAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGT
51      -----CCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGT
          *****

54      TAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCCTAGG
48      TAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCCTAGG
51      TAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCCTAGG
          ***** ** *****

54      TACAATGCTAGCAAAGAGGAGAAAAC TAGATGTCCGAATCAGGGGGTGGGACTGTTGCTA
48      TACAGTGTAGCAAAGAGGAGAAAAC TAGATGTCCGAATCAGGGGGTGGGACTGTTGCTA
51      TATTATGCTAGCAAAGAGGAGAAAAC TAGATGTCCGAATCAGGGGGTGGGACTGTTGCTA
          ** *****

54      CCGCACGTCACGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGG
48      CCGCACGTCACGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGG
51      CCGCACGTCACGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGG
          *****

54      GTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTCGCAGTTTATGTGAAGGTC
48      GTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTCGCAGTTTATGTGAAGGTC
51      GTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTCGCAGTTTATGTGAAGGTC
          *****
    
```



```
GHRNTLDTTDRDWEVMFGVNVFGLTKAIRRFIEPMRAQQRGSIVNVVSSGVLAVAAGGGY
HGLRPWTVEMPYQATKAAVMALTFYLAEEVRGDGVAVNAIMPGHTRASWFDATARAFNEQ
GIAYFMRPAIPEHLLPISLFLAAQDSAGASGRLYYVPEWNYDHGYGDYAXWQDHELPPDM
EEIYSRLEAATXXYEXXGVXXLPFDXXGXLXAXGMANLGXXXXXIXMXXIXXXLXXX
```

Blast of clone CCA-48 (forward)
 fluoren-9-ol dehydrogenase [Janibacter terrae]

Alignment statistics for match #1

Score	Expect	Method	Identities	Positives	Gaps
623 bits(1606)	0.0	Compositional matrix adjust.	305/306(99%)	305/306(99%)	0/306
Query 1	MSESGGGTVATARQRQLVERALGEWQGEVAGRIVVVTGGARGIGRSLCEGLLRAGAKVVA		60		
Sbjct 1	MSESGGGTVATARQRQLVERALGEWQGEVAGRIVVVTGGARGIGRSLCEGLLRAGAKVVA		60		
Query 61	ADLTWDDADDFRKQLESDGSGMAVDMITDDDALDAARDAVIDRFGTVDVLVNNASLVSE		120		
Sbjct 61	ADLTWDDADDFRKQLESDGSGMAVDMITDDDALDAARDAVIDRFGTVDVLVNNASLVSE		120		
Query 121	TLFPPTGHRNTLDTTDRDWEVMFGVNVFGLTKAIRRFIEPMRAQQRGSIVNVVSSGVLAV		180		
Sbjct 121	TLFPPTGHRNTLDTTDRDWEVMFGVNVFGLTKAIRRFIEPMRAQQRGSIVNVVSSGVLAV		180		
Query 181	AAGGGYHGLRPWTVEMPYQATKAAVMALTFYLAEEVRGDGVAVNAIMPGHTRASWFDATA		240		
Sbjct 181	AAGGGYHGLRPWTVEMPYQATKAAVMALTFYLAEEVRGDGVAVNAIMPGHTRASWFDATA		240		
Query 241	RAFNEQGIAYFMRPAIPEHLLPISLFLAAQDSAGASGRLYYVPEWNYDHGYGDYAXWQDH		300		
Sbjct 241	RAFNEQGIAYFMRPAIPEHLLPISLFLAAQDSAGASGRLYYVPEWNYDHGYGDYAAWQDH		300		
Query 301	ELPPDM 306				
Sbjct 301	ELPPDM 306				

Translation

>48-BBaG00100-R 1274 289 277 0.05

```
XXXXXXXXLXXXXAAXXXLLGXXFXLWDEDAIPLLEDYLGADPLWYLXLAFEAPLGGLEVIG
PTMKFRIKGNWKLAAENFAGXDYHGLYTHGATFQIGFLPDYDTLXDYIAHFDHGHGMGDI
GKPGRAYQNDVGMAGFLGPEAIEYVNAVHERLKRVSPLQAEMHGLGQGNIFPNLSWIKF
GVFHVFLFQXHPRGPEIEVWXTALXDRDAPXSXKDXAXXMSQXXAXAGIFXXDDGEN
FEXITEXXXGVXSHTXDFHYAMXXGHEGXXXXXXTPVHLXXXYSXNXRISXGIGXELMT
XSXNXXXXPSXXXXDHLLNEFTXXXXXXXXVLXXQXXXXXXXXXXXX-SXXXXTSLXXAHX
```

Blast

aromatic ring-hydroxylating dioxygenase subunit alpha [Janibacter terrae]

Sequence ID: [WP_032491530.1](#) Length: 443 Number of Matches: 1

Alignment statistics for match #1

Score	Expect	Method	Identities	Positives	Gaps
328 bits(842)	1e-109	Compositional matrix adjust.	160/180(89%)	161/180(89%)	0/180
Query 1	PTMKFRIKGNWKLAAENFAGXDYHGLYTHGATFQIGFLPDYDTLXDYIAHFDHGHGMGDI		60		
Sbjct	PTMKFRIK NWKLAAENFAG DYH LYTHG+ FQIGFLPDYDTL DYIAHFDHGHGMGDI				

Sbjct	198	PTMKFRIKANWKLAAENFAGDDYHVLYTHGSAFQIGFLPDYDTLGDYIAHFDHGHGMGDI	257
Query	61	GKPGRAYQNDVGMQAQFLGPEAIEYVNAVHERLKARVSPLQAEMHGLGQGNIFPNLSWIKF GKPGRAYQNDVGMQAQFLGPEAIEYVNAVHERLKARVSPLQAEMHGLGQGNIFPNLSWIKF	120
Sbjct	258	GKPGRAYQNDVGMQAQFLGPEAIEYVNAVHERLKARVSPLQAEMHGLGQGNIFPNLSWIKF	317
Query	121	GVFHVFLFQXHPRGPEIEVWXTALXDRDAPXSXKDXAXXXMSQXXAXAGIFXXDDGEN GVFHVFLFQ HPRGPGEIEVW TAL DRDAP S KD A MSQ A AGIF DDGEN	180
Sbjct	318	GVFHVFLFQWHPRGPGEIEVWQTALFDRDAPQSVKDVARTQMSQENAAAGIFGQDDGEN	377

DNA sequences of clones CCA-48, CCA-51, and CCA-54 obtained using primers BBaG00100-F and BBaG00100-R

>48-BBaG00100-F

GAAGAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGTTAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTGACGGCT
AGCTCAGTCTTAGGTACAGTGTAGCAAAGAGGAGAAAAGTAGATGTCCGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGGAAACGTG
CATTGGGCGAGTGGCAAGGAGAAGTGGCGGGTTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTTACGCGCAGGT
GCCAAGGTCGTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGA
CGATGATGCCTTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCCAC
CAACGGGGCACCGTAATACCTGGACACGACAGATCGCGACTGGGAGGTAATGTTGGTGTGAATGTCTTGGAACTTAAGGCGATTTCGTCGCTTCATCGAG
CCAATGCGCGCTCAACAGCGCGGTTGATTGTCAACGTGGTAAGCAGTGGCGTCTTCGAGTCGCGAGTGGCGGGGATAACCATGGCTTCGCTCCATGGACCGT
TGAGATGCCCTATCAGGCTACTAAAGCAGCTGTCTATGGCTCTTACATTCTACTTGGCCGAAGAGGTGCGCGCGATGGGGTGGCGGTCAATGCTATCATGCCTG
GTCACACCCGCGCTTCTTGGTTTGTATGCGACCGCTCGTGCCTTTAATGAGCAGGGGATCGCATACTTCATGCGCCCTGCTATTCCCGAGCACTTGCTTCTC
TCTTGTCTCTTCGAGCGCAGGACTCCGCTGGCGCCTCTGGGCGTCTTTACTATGTGCCGAGTGGAACTACGACCACGGGTACGGCGATTATGCGGNGTGGCA
GGATCATGAATTACCACCGGACATGGAGGAAATCTATTGCGCTCTGGAAGCAGCAACGCCNANTTATGAACGNACAGGTGTANNNCNNTTGGCGTTTGTATGCNN
AAGGAACNTTGNNTGCCNAGGAATGGCTAATCTTGGAAANANNGAANNNTNGGATTNCAATGNNTNNNNNATAANTNNAANTTAANCCGNNA

>51-BBaG00100-F

CCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGTTAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACAGCTAGTCTAGTCTCTA
GGTATTATGCTAGCAAAGAGGAGAAAAGTAGATGTCCGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGGAAACGTGCATTGGGCGAGT
GGCAAGGAGAAGTGGCGGGTTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTTACGCGCAGGTGCCAAGGTCGTG
GCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCCTT
GGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCAC
GTAATACCTGGACACGACAGATCGCGACTGGGAGGTAATGTTGGTGTGAATGTCTTGGAACTTAAGGCGATTTCGTCGCTTCATCGAGCCAATGCGCGCT
CAACAGCGCGGTTGATTGTCAACGTGGTAAGCAGTGGCGTCTTCGAGTCGCGAGTGGCGGGGATAACCATGGCTTCGCTCCATGGACCGTTGAGATGCCCTA
TCAGGCTACTAAAGCAGCTGTCTATGGCTCTTACATTCTACTTGGCCGAAGAGGTGCGCGCGATGGGGTGGCGGTCAATGCTATCATGCCTGGTCACACCCGCG
CTTCTTGGTTTGTATGCGACCGCTCGTGCCTTTAATGAGCAGGGGATCGCATACTTCATGCGCCCTGCTATTCCCGAGCACTTGCTTCTCCTATCTCCTTGTCTC
GCAGCGCAGGACTCCGCTGGCGCCTCTGGGCGTCTTTACTATGTGCCGANTGGAACACGACCACGGGTACGGCGATTATGCGGNGTGGCAGGATCATG

>54-BBaG00100-F 1279 19 1006 0.05

TGCGAGAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGTTAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACGG
CTAGCTCAGTCTTAGGTACAATGCTAGCAAAGAGGAGAAAAGTAGATGTCCGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGGAAACG
TGCAATGGGCGAGTGGCAAGGAGAAGTGGCGGGTTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTTACGCGCAG
GTGCCAAGGTCGTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACA
GACGATGATGCCTTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCC
ACCAACGGGGCACCGTAATACCTGGACACGACAGATCGCGACTGGGAGGTAATGTTGGTGTGAATGTCTTGGAACTTAAGGCGATTTCGTCGCTTCATCG
AGCCAATGCGCGCTCAACAGCGCGGTTGATTGTCAACGTGGTAAGCAGTGGCGTCTTCGAGTCGCGAGTGGCGGGGATAACCATGGCTTCGCTCCATGGACC
GTTGAGATGCCCTATCAGGCTACTAAAGCAGCTGTCTATGGCTCTTACATTCTACTTGGCCGAAGAGGTGCGCGCGATGGGGTGGCGGTCAATGCTATCATGCC
TGGTCACACCCGCGCTTCTTGGTTTGTATGCGACCGCTCGTGCCTTTAATGAGCAGGGGATCGCATACTTCATGCGCCCTGCTATTCCCGAGCACTTGCTTCTC
TCTCCTTGTCTCTTCGAGCGCAGGACTCCGCTGGGNCCTCTGGGCGTCTTTACTATGNTGCCGAGTGGAACTACCAACNCGGGTTACGGCGAAT

>54-BBaG00100-F 1279 19 1006 0.05

TGCGAGAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGCTACGTAATTAGTTAGTTAGCCCTTAGTGACTCGAATTCGCGGCCGCTTCTAGATTTACGG
 CTAGCTCAGTCCTAGGTACAATGCTAGCAAAGAGGAGAAAAGTACGATGTCGGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGGAAACG
 TGCAATGGGGCAGTGGCAAGGAGAAGTGGCGGGTCGCGTAATTGTGGTAAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTACGCGCAG
 GTGCCAAGGTGCTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACA
 GACGATGATGCCTTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCC
 ACCAACGGGGCACCGTAATACCTTGGACACGACAGATCGGACTGGGAGGTAATGTTTGGTGTGAATGCTTTTGGAACTTAAGGCGATTGCTGCTTCATCG
 AGCCAATGCGCGCTCAACAGCGCGGTTTCGATTGTCAACGTGGTAAAGCAGTGGCGTCCCTGACAGTTCGCGAGTGGCGGGGATACCATGGCTTGGCTCCATGGACC
 GTTGAGATGCCCTATCAGGCTACTAAAGCAGCTGTCAATGGCTTACATTCTACTTGGCCGAAGAGGTGCGCGGGCATGGGGTGGCGGTCAATGCTATCATGCC
 TGGTCACACCCGCGCTTCTTGGTTTGTATGCGACCGCTCGTGCCTTTAATGAGCAGGGGATCGCATACTTCATGCGCCCTGCTATTCCCGAGCACTTGCTTCTCA
 TCTCTTGTTCCTTGCAGCGCAGGACTCCGCTGGGNCCTCTGGGGCGTCTTTACTATGNTGCCCGAGTGGAACTACCAACNCGGGTTACGGCGAATNNN

>48R

ATTNCNCGNGGCCGCTGNGNTNNACTACTAGGCNTNNATTTTNCNCTTTGGGACGAGGACGCTATTTCCCTGGAGGATTATCTGGGGGCAGATCCCTTATGGT
 ACTTANACTTAGCCTTCGAAGCACCCTGGGGGGTTAGAGGTGATTGGTCCGACGATGAAATTTTCGATTAAGGGCAATTGGAACTTGGCGCGGAGAACTTC
 GCCGANACGACTACCACGGTTTGTACACCCATGGGGCAACTTCCAGATTGGGTTCTGCCCGATTACGACACACTTGGNGATTACATCGCTCATTTTGACCA
 CGGACACGGAATGGGCGACATTTGTAAGCCTGGACGCGCTATCAAAACGACGTGGGGATGGCGCAATTTCTGGGCCCTGAGGCTATCGAATATGTTAACGCGG
 TGCATGAACGTTTGAAGCCCGTGTCTCACCGTTGCAAGCGGAGATGCACGGCTTAGGGCAAGGAAATATTTTCCAAATTTATCATGGATTAAATTCGGCGTG
 TTTACGTTGTTGGATTGTTTTCAGNGGCATCCAGTGGGCCGGGGAAATCGAAGTCTGNNAGACTGCCTTGNTTGTATCGTGTATGCTCCGANAGTGNCAAGGA
 CGTNGCTCANNANTCANNATGAGTCAGGANAANGCGNNTGCCGAATCTTTGNCANGACGACGAGAGAAAATTTGAACANATCACCGAATCANNNTCGNNGCGTTG
 NCAGTCACACTCGNGACTTCCACTACGCTATGGNNTGGGGCATGAGGGGNAANCNTGANGAAGGTNACCCCTGTCCATTTANGNCCNNACTATTCCGAGCAN
 AACNCTCGCATTTCTANCGGTATTGGCTGNGAGCTGATGACANCGTCGNANAACAANNANTTGNCCCTCGGANGNTANNNNCGACCACCTCCTTAACGAATT
 TACGNNNNNACGNANANGANTNTN

>51R

TCTACNGGATTTTTTTTTCCNATTTGGGACGAGGACGCTATTTCCCTGTAGGATTATCTGGGGGCAGATCCCTTATGGNCTGAGACTTAGCCTTCGAAGCACC
 ACTGGGGGGNTAGAGGTGATTGGACCGACGATGAAATTTTCGATTAAGGGCAATTTGTAACCTTGGCGCGGAAAACCTTCGCCGGAGACGACTACCACGATTTGT
 ACACCCATGGGTCAACTTCCAGATTGGCTTCTGCCCGATTACGACACACTTGGNGATTACATCGCTCATTTTGACCACGGACACGGAATGGGCGACATTGGT
 AAGCCTGGACGCGCTATCAAAACGACGAGGGGATGGCGCAATTTCTGGGCCCTGAGGCTATCGAATATGTTAACGCGGTGCATGAACGTTTGAAGCCCGTGT
 CTCACCGTTGCAAGCGGAGATGCACGGCTTAGGGCAAGGAAATATTTTCCAAATTTATCATGGATTAAATTCGGCGTGTTCACGTTTGGATTGTTTCAGT
 GGCATCCAGTGGGCCGGGGAAATCGAAGTCTGGCAGACTGCCTTGTGATCGTGTATGCTCCGAAAAGTGTCAAGGACGTNGCTCGTACTCAGATGAGTCAG
 GAGAATGCGGCTGCCGAATCTTTGGCCAGNACGACGAGANAATAATTTGAACAAATCACCGAATCAGCTCGNNGCGTGNAGNACAGACTCGCGACTTCNCTAC
 GCTATGGGTTNGGGCATGAGGGGAAATCCATGAGGAAGGTACCCTGNCAATTTANGTCCGCACTATTTCCGAGCANAACCATCGCAATTTCTACCGNTATTGGN
 NGGAGTTGANGACAACGCCAGGAGACCAANTANGCCCCATCGGAGGATANTAT

>54R

GCCCGNGNGGAAAACATAAGGGATTTATTTTCGCAANNNTGGGACGAGGACGCTATTTCCCTGGTGGATTATCTGGGGGCATATCCTTATGGNACTTAGACTT
 AGCCTTCGAAGCACCCTGGGGGGTTAGAGGTGATTGGACCGACGATGAAATTTTCGATTAAGGGCAATTTGGAACTTGGCGCGGAGAACTTCGCCGGAGACG
 ACTACCACGTTTGTACACCCATGGGTGAGCTTCCAGATTGGCTTCTGCCCGATTACGACACACTTGGCGATTACATCGCTCATTTTGACCACGGACACGGA
 ATGGGCGACATTTGTAAGCCTGGACGCGCTATCAAAACGACGTGGGGATGGCGCAATTTCTGGGCCCTGAGGCTATCGAATATGTTAACGCGGTGCATGAACG
 TTTGAAGCCCGTGTCTCACCGTTGCAAGCGGAGATGCACGGCTTAGGGCAAGGAAATATTTTCCAAATTTATCATGGATTAAATTCGGCGTGTTCACGTTG
 TTGGATTGTTTTCAGTGGCATCCAGTGGNCGNNGGAAATCGAAGTCTGGCAGACTGCCTTGTGATCGTGTATGCTCCGANAGNNGCAAGGACGTTGCTCGT
 ACTCAGATGANTCNGGANAATGCGGNTGCCGAATCTTTGNCAGNACGACNGAGAANTTNTGAACAAATCACCGAATCAGCTCANNCGTGNAGTCANACTC
 GCGACTCCNCTACNNTATGGGTTTGGGGCATGANGGGGAAATCCATGAGGAAAGNTACCCTGGTCAATTTAGGNCCGCACTATCCGAGNANGATNCATCGCAAT
 TNC

Sequence verification

Phenanthrene clones CCA-57, CCA-60, and CCA-64.

Sequence Alignment of clones CCA-57, CCA-60, and CCA-64.

Clustal

Forward

```
64F          -----TCACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGAT
57-BBaG00100-F  -GGCGTATCACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGAT
```

```

60F          AGGCGTATNACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGAT
              * *****
64F          TTCTGGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCCTAGGTACAATGCTAG
57-BBaG00100-F TTCTGGAATTCGCGGCCGCTTCTAGATTGACGGCTAGCTCAGTCCTAGGTACAGTGCTAG
60F          TTCTGGAATTCGCGGCCGCTTCTAGATTTACAGCTAGCTCAGTCCTAGGTATTATGCTAG
              ***** ** *****
64F          CAAAGAGGAGAAAAC TAGATGGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCG
57-BBaG00100-F CAAAGAGGAGAAAAC TAGATGGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCG
60F          CAAAGAGGAGAAAAC TAGATGGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCG
              *****
64F          ACTGCTGCCGCGACCTTTGAACGTCGTAGTCCTACAACCGGCGAATTAGTGACTCGCGCC
57-BBaG00100-F ACTGCTGCCGCGACCTTTGAACGTCGTAGTCCTACAACCGGCGAATTAGTGACTCGCGCC
60F          ACTGCTGCCGCGACCTTTGAACGTCGTAGTCCTACAACCGGCGAATTAGTGACTCGCGCC
              *****
64F          GCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGACTCTGCTGCCGCAGCTTATCGT
57-BBaG00100-F GCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGACTCTGCTGCCGCAGCTTATCGT
60F          GCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGACTCTGCTGCCGCAGCTTATCGT
              *****
64F          TCCTGGAGCACTACTGGGCCACCGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTA
57-BBaG00100-F TCCTGGAGCACTACTGGGCCACCGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTA
60F          TCCTGGAGCACTACTGGGCCACCGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTA
              *****

```

Clustal Reverse sequences

```

57R          --ACGACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTT
60R          CGACGACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTT
64R          ---GACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTT
              *****
57R          TCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACACCTTGCCCTTTTTTGCCGACTGCAGC
60R          TCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACACCTTGCCCTTTTTTGCCGACTGCAGC
64R          TCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACACCTTGCCCTTTTTTGCCGACTGCAGC
              *****
57R          GGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTA
60R          GGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTA
64R          GGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTA
              *****
57R          GTAGAGAGCGTTCACCGACAAACAACAGATAAAACGAAAGGCCAGTCTTTCGACTGAGC
60R          GTAGAGAGCGTTCACCGACAAACAACAGATAAAACGAAAGGCCAGTCTTTCGACTGAGC
64R          GTAGAGAGCGTTCACCGACAAACAACAGATAAAACGAAAGGCCAGTCTTTCGACTGAGC
              *****
57R          CTTTCGTTTTATTTGATGCCTGGTTATTAATGCAAAATCATCACCACCACACACCTTCATA
60R          CTTTCGTTTTATTTGATGCCTGGTTATTAATGCAAAATCATCACCACCACACACCTTCATA
64R          CTTTCGTTTTATTTGATGCCTGGTTATTAATGCAAAATCATCACCACCACACACCTTCATA
              *****

```

```

57R      ATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAA
60R      ATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAA
64R      ATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAA
          *****

57R      TTTTCGCGGGAGGCCAACAGGACATAAAGTCCACAGTAATCGGCCGGAACGGGGGCCATA
60R      TTTTCGCGGGAGGCCAACAGGACATAAAGTCCACAGTAATCGGCCGGAACGGGGGCCATA
64R      TTTTCGCGGGAGGCCAACAGGACATAAAGTCCACAGTAATCGGCCGGAACGGGGGCCATA
          *****

```

DNA sequences of clones CCA-57, CCA-60, and CCA-64 obtained using primers BBaG00100-F and BBaG00100-R

>57-BBaG00100-F

```

GGCGTATCACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGATTTCTGGAATTCGCGGCCGCTTCTAGATTGACGGC
TAGCTCAGTCCTAGGTACAGTGCTAGCAAAGAGGAGAAAAGTATGATGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCGACTGCTG
CCGCGACCTTTGAACGTCGTAGTCTACAACCGGCGAATTAGTGACTCGCGCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGACT
CTGCTGCCGACGCTTATCGTTCCGGAGCACTACTGGGCCACCAGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTATTAGAAGCCCGCA
CCCCTGAGTTTTACGCGTCATGGCCTTGAAGTGGGTGCATCAGACCTTTGGGCTGGCGTTAACGTAATGCTTGCAGCTAACTTATTTTCGCG
AGGCTGCTGCATTAACAACGCAAATTCAGGGGGAGACCATTCCAACAGACAAGGCTGGTGTCTTATCGATGACTGTTTCGTCAGCCGTTGGCG
TGATTTTATCTATTGCCCGTGAATGGACCTGTTGTATTAGCTGCTCGTGCCATCGCGTACCCATTGGTCTGCGGAAACACAGTAGTGTTC
GTGCGTCGGAGTTATCTCCAAGACGCACATGCTTATCGTGGATGTCTTACGTGACGCTGGCCTGCCGCCAGGAGTCTTGAATGCAGTTACCA
ACGCACCTCAAGATGCACCAGAAGTAGTGGACGCGTTAATCGCACACCCTGCGGTGCGTCGCATTAACTTTACGGGGAGTACCCGTGTGGGGC
GTGTAATTGCAGAAAAAGCAGCTCGTCACTTAAAGCGTTGCCCTTCTTGGAGCTTGGAGGCAAGGCCCCCTGGTANTTTTAAATGATGCGGATA
TCGACGAAGCTGTTAAAGCCGCCGCTTTCGGTGCCTTCTTGTACCAAGGGCAGATTTGNNNTGTCGACAGAACGTATTGTTGTGCGANNAAAAT
TGCTGATACTTTTCGTTGCCCGTTTCGCCGCC

```

>60F

```

AGGCGTATNACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGATTTCTGGAATTCGCGGCCGCTTCTAGATTTACAG
CTAGCTCAGTCCTAGGTATTATGCTAGCAAAGAGGAGAAAAGTATGATGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCGACTGCT
GCCGCGACCTTTGAACGTCGTAGTCTACAACCGGCGAATTAGTGACTCGCGCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGAC
TCTGCTGCCGACGCTTATCGTTCCGGAGCACTACTGGGCCACCAGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTATTAGAAGCCCGC
ACCCCTGAGTTTTACGCGTCATGGCCTTGAAGTGGGTGCATCAGACCTTTGGGCTGGCGTTAACGTAATGCTTGCAGCTAACTTATTTTCG
GAGGCTGCTGCATTAACAACGCAAATTCAGGGGGAGACCATTCCAACAGACAAGGCTGGTGTCTTATCGATGACTGTTTCGTCAGCCGTTGGC
GTGATTTTATCTATTGCCCGTGAATGGACCTGTTGTATTAGCTGCTCGTGCCATCGCGTACCCATTGGTCTGCGGAAACACAGTAGTGTTC
CGTGCCTCGGAGTTATCTCCAAGACGCACATGCTTATCGTGGATGTCTTACGTGACGCTGGCCTGCCGCCAGGAGTCTTGAATGCAGTTACC
AACGCACCTCAAGATGCACCAGAAGTAGTGGACGCGTTAATCGCACACCCTGCGGTGCGTCGCATTAACTTTACGGGGAGTACCCGTGTGGGG
CGTGTAAATTGCAGAAAAAGCAGCTCGTCACTTAAAGCGTTGCCCTTCTTGGAGCTTGGAGGCAAGGCCCCCTGGTAGTTTTAGATGATGCGGAT
ATCGACGAAGCTGTTAAAGCCGCCGCTTTCGGTGCCTTCTTGTACCAAGGGCAGATTTGTATGTCGACAGAACGTATTGTTGTGCGACGAAAA
ATTGCTGATACTTTTCGTTGCCCGTTTCNCCGCCCGCCNNAATTA

```

>64F

```

TCACGAGGCAGAATTTTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGATGATTTCTGGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTC
AGTCCTAGGTACAATGCTAGCAAAGAGGAGAAAAGTATGATGATACGCAACTTATCATTGATAACGCAGACGTTCCGGCGACTGCTGCCGCGA
CCTTTGAACGTCGTAGTCTACAACCGGCGAATTAGTGACTCGCGCCGCCGCCAGCGTCGCTGACGCAATTGCAGCCGCTGACTCTGCTG
CCGACGCTTATCGTTCCGGAGCACTACTGGGCCACCAGAGCGCCGCCGCATCTTGTGAAAGCCGCCGATTTATTAGAAGCCCGCACCCCTG
AGTTTTACGCGTCATGGCCTTGAAGTGGGTGCATCAGACCTTTGGGCTGGCGTTAACGTAATGCTTGCAGCTAACTTATTTTCGCGAGGCTG
CTGCATTAACAACGCAAATTCAGGGGGAGACCATTCCAACAGACAAGGCTGGTGTCTTATCGATGACTGTTTCGTCAGCCGTTGGCGTGATTT
TATCTATTGCCCGTGAATGGACCTGTTGTATTAGCTGCTCGTGCCATCGCGTACCCATTGGTCTGCGGAAACACAGTAGTGTTCGTCGCT
CGGAGTTATCTCCAAGACGCACATGCTTATCGTGGATGTCTTACGTGACGCTGGCCTGCCGCCAGGAGTCTTGAATGCAGTTACCAACGCAC
CTCAAGATGCACCAGAAGTAGTGGACGCGTTAATCGCACACCCTGCGGTGCGTCGCATTAACTTTACGGGGAGTACCCGTGTGGGGCGTGTA
TTGCAGAAAAAGCAGCTCGTCACTTAAAGCGTTGCCCTTCTTGGAGCTTGGAGGCAAGGCCCCCTGGTAGTTTTAGATGATGCGGATATCGACG
AAGCTGTTAAAGCCGCCGCTTTCGGTGCCTTCTTGTACCAAGGGCAGATTTGTATGTCGACAGAACGTATTGTTGTGCGACGAAAAAATTGCTG
ATACTTTCGTTGCCCGTTTCNCCNCCNGCGCCCGCAATTACCTGNGGGCGATCTGNACCTGNNGNNGN

```

>57R

ACGACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTTTCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACACC
 TTGCCCTTTTTTGCCGGACTGCAGCGGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTAGTAGAGAG
 CGTTCACCGACAAACAACAGATAAAAACGAAAGGCCAGTCTTTTCGACTGAGCCTTTCGTTTTATTGATGCCTGGTTATTAATGCAAATCATC
 ACCACCACACACCTTCATAATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAATTTTCGCGGGAGGC
 CAACAGGACATAAGGTCCACAGTAATCGGCCGGAACGGGGGCCATACGCAGGGGAGTATTTTCCACCAGCATAGCCTCCAAATTTGTCAACTTG
 ATCAAGGGTTTTGAGCTGACGTACCAGTTGCCGACAAACCTCCCAAATTTGTGCGCATCCCCCTGGAGCGACGCCATTTACACGTACTTTTGG
 GGCCAACTCGTGAGCCAATTCACGAATCAGTCCCACCAAAGCATGTTTACTTGCAGTGTACAAGGGGCCACCCCCCGGGGTAGAAACCTGC
 ATTAGAAACGGTAAAGATTAACGAACCCTGGGATTTTACCAGTTCAGCTAACGCTGCTTTGGCCCCAAGCAGTGCAGCTTTTACGTTCACTGC
 GAATAATTCGTCGAATGCGCGTGAGATGCTCCCGCGTCCATCTGGGGCAATGTCTGGAAAAAATCGAAGACTCCGGCGTTTCCGACAAAGTT
 GTCAGGGCGACCGAAACGTGCAACTGTTTTCTGGACCACAAGAAGCTTATCTTCATACAACGTCACGTCTCCTACCCTACTTCTACTGCGTC
 CCCAAATTCACGAGCAAGCTCACGGGCGGTTCCGCCAGAGCGTTCAAGGACCCCCACACGTCCACCCTCATTGATGAAACGTTCCACTAATGC
 TTTGCCAGTCCA

>60R

CGACGACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTTTCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACA
 CCTTGCCCTTTTTTGCCGGACTGCAGCGGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTAGTAGAG
 AGCGTTCACCGACAAACAACAGATAAAAACGAAAGGCCAGTCTTTTCGACTGAGCCTTTCGTTTTATTGATGCCTGGTTATTAATGCAAATCA
 TCACCACCACACACCTTCATAATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAATTTTCGCGGGAG
 GCCAACAGGACATAAGGTCCACAGTAATCGGCCGGAACGGGGGCCATACGCAGGGGAGTATTTTCCACCAGCATAGCCTCCAAATTTGTCAACT
 TGATCAAGGGTTTTGAGCTGACGTACCAGTTGCCGACAAACCTCCCAAATTTGTGCGCATCCCCCTGGAGCGACGCCATTTACACGTACTTTT
 GGGGCCAACTCGTGAGCCAATTCACGAATCAGTCCCACCAAAGCATGTTTACTTGCAGTGTACAAGGGGCCACCCCCCGGGGTAGAAACCT
 GCATTAGAAACGGTAAAGATTAACGAACCCTGGGATTTTACCAGTTCAGCTAACGCTGCTTTGGCCCCAAGCAGTGCAGCTTTTACGTTCACT
 GCGAATAATTCGTCGAATGCGCGTGAGATGCTCCCGCGTCCATCTGGGGCAATGTCTGGAAAAAATCGAAGACTCCGGCGTTTCCGACAAAG
 TTGTCAAGGGCGACCGAAACGTGCAACTGTTTTCTGGACCACAAGAAGCTTATCTTCATACAACGTCACGTCTCCTACCCTACTTCTACTGCG
 TCCCCAAATTCACGAGCAAGCTCACGGGCGGTTCCGCCAGAGCGTTCAAGGACCCCCACACGTCCACCCTCATTGATGAAACGTTCCACTAAT
 GCTTTGCCAGTCCANAGCCNCCACCCGGNATGATTGGCACTGGATTATTCAACCAGC

>64R

GACCGAGCGCAGCGAGTCAGTGAGCGAGGAAGCCTGCATAACGCGAAGTAATCTTTTCGGTTTTAAAGAAAAAGGGCAGGGTGGTGACACCTT
 GCCCTTTTTTGCCGGACTGCAGCGGCCGCTACTAGTATATAAACGCAGAAAGGCCACCCGAAGGTGAGCCAGTGTGACTCTAGTAGAGAGCG
 TTCACCGACAAACAACAGATAAAAACGAAAGGCCAGTCTTTTCGACTGAGCCTTTCGTTTTATTGATGCCTGGTTATTAATGCAAATCATCAC
 CACCACACACCTTCATAATGCCACGAACACCGAATCCGCCATCTGTATTGATGACCACCCCTGTCATCGGACCAGAATTTTCGCGGGAGGCCA
 ACAGGACATAAGGTCCACAGTAATCGGCCGGAACGGGGGCCATACGCAGGGGAGTATTTTCCACCAGCATAGCCTCCAAATTTGTCAACTTGAT
 CAAGGGTTTTGAGCTGACGTACCAGTTGCCGACAAACCTCCCAAATTTGTGCGCATCCCCCTGGAGCGACGCCATTTACACGTACTTTTGGGG
 CCAACTCGTGAGCCAATTCACGAATCAGTCCCACCAAAGCATGTTTACTTGCAGTGTACAAGGGGCCACCCCCCGGGGTAGAAACCTGCAT
 TAGAAACGGTAAAGATTAACGAACCCTGGGATTTTACCAGTTCAGCTAACGCTGCTTTGGCCCCAAGCAGTGCAGCTTTTACGTTCACTGCGA
 ATAATTCGTCGAATGCGCGTGAGATGCTCCCGCGTCCATCTGGGGCAATGTCTGGAAAAAATCGAAGACTCCGGCGTTTCCGACAAAGTTGT
 CAAGGCGACCGAAACGTGCAACTGTTTTCTGGACCACAAGAAGCTTATCTTCATACAACGTCACGTCTCCTACCCTACTTCTACTGCGTCCC
 CAAATTCACGAGCAAGCTCACGGGCGGTTCCGCCAGAGCGTTCAAGGACCCCCACACGTCCACCCTCATTGATGAAACGTTCCACTAATGGT

DNA sequence of clone CCA-57-Forward blast

57Fdehydrogenase [Burkholderia sp. Bk]

Sequence ID: [WP_088176951.1](#) Length: 483 Number of Matches: 1

Alignment statistics for match #1

Score	Expect	Method	Identities	Positives	Gaps
572 bits(1474)	0.0	Compositional matrix adjust.	289/296(98%)	291/296(98%)	0/296
Query 1		MDTQLIIDNADVPATAAATFERRSPTTGELVTRAAAASVADAIAAADSAAAAYRSWSTTG	60		
		MDTQLIIDNADVPATAAATFERRSPTTGELVTRAAAASVADAIAAADSAAAAYRSWSTTG			
Sbjct 1		MDTQLIIDNADVPATAAATFERRSPTTGELVTRAAAASVADAIAAADSAAAAYRSWSTTG	60		

Query	61	PTERRRILLKAADLLEARTPEFSRVMALEVGASDLWAGVNVMLAANLFREAAALTTQIQG	120
		PTERRRILLKAADLLEARTPEFSRVMALEVGASDLWAGVNVMLAANLFREAAALTTQIQG	
Sbjct	61	PTERRRILLKAADLLEARTPEFSRVMALEVGASDLWAGVNVMLAANLFREAAALTTQIQG	120
Query	121	ETIPTDKAGVLSMTVRQPVGVILSIAPWNGPVVLAARAIAIYPLVCGNTVVFRASELSPKT	180
		ETIPTDKAGVLSMTVRQPVGVILSIAPWNGPVVLAARAIAIYPLVCGNTVVFRASELSPKT	
Sbjct	121	ETIPTDKAGVLSMTVRQPVGVILSIAPWNGPVVLAARAIAIYPLVCGNTVVFRASELSPKT	180
Query	181	HMLIVDVLRDAGLPPGVLNAVNTNAPQDAPEVVDALIAHPAVRRINFTGSTRVGRVIAEKA	240
		HMLIVDVLRDAGLPPGVLNAVNTNAPQDAPEVVDALIAHPAVRRINFTGSTRVGRVIAEKA	
Sbjct	181	HMLIVDVLRDAGLPPGVLNAVNTNAPQDAPEVVDALIAHPAVRRINFTGSTRVGRVIAEKA	240
Query	241	ARHLKRCLLELGGKAPLVXLNDADIDEAVKAAVFGAFLYQGQIXXSTERIVVXXKL	296
		ARHLKRCLLELGGKAPLV L+DADIDEAVKAAVFGAFLYQGQI STERIVV K+	
Sbjct	241	ARHLKRCLLELGGKAPLVVLDADIDEAVKAAVFGAFLYQGQICMSTERIVVDEKI	296

DNA sequence of clone CCA-57-Forward Reverse

3-(cis-5,6-dihydroxycyclohexa-1,3-dien-1-yl)propanoate dehydrogenase [Burkholderia sp. Bk]

Sequence ID: [WP_088176957.1](#) Length: 272 Number of Matches: 1

Alignment statistics for match #1

Score	Expect	Method	Identities	Positives	Gaps
511 bits(1316)	0.0	Compositional matrix adjust.	256/256(100%)	256/256(100%)	0/256
Query 1		GLGKALVERFINEGGRVGVLSRGERARELAREFGDAVEVVGDVTLYEDNVLVVQKTVA	60		
		GLGKALVERFINEGGRVGVLSRGERARELAREFGDAVEVVGDVTLYEDNVLVVQKTVA			
Sbjct 17		GLGKALVERFINEGGRVGVLSRGERARELAREFGDAVEVVGDVTLYEDNVLVVQKTVA	76		
Query 61		RFGRLDNFVGNAGVFDFQTLPQMDAGSISRAFDELFAVNVKAALLGAKAALAEVLKSQG	120		
		RFGRLDNFVGNAGVFDFQTLPQMDAGSISRAFDELFAVNVKAALLGAKAALAEVLKSQG			
Sbjct 77		RFGRLDNFVGNAGVFDFQTLPQMDAGSISRAFDELFAVNVKAALLGAKAALAEVLKSQG	136		
Query 121		SLIFTVSNAGFYPPGGGPLYTASKHALVGLIRELAHELAPKVRVNGVAPGGMRTNLGGLS	180		
		SLIFTVSNAGFYPPGGGPLYTASKHALVGLIRELAHELAPKVRVNGVAPGGMRTNLGGLS			
Sbjct 137		SLIFTVSNAGFYPPGGGPLYTASKHALVGLIRELAHELAPKVRVNGVAPGGMRTNLGGLS	196		
Query 181		ATG TSAQTLDQVDNLEAMLVENTPLRMAPVPADYCGPYVLLASRENSGPMTGVVINTDGG	240		
		ATG TSAQTLDQVDNLEAMLVENTPLRMAPVPADYCGPYVLLASRENSGPMTGVVINTDGG			
Sbjct 197		ATG TSAQTLDQVDNLEAMLVENTPLRMAPVPADYCGPYVLLASRENSGPMTGVVINTDGG	256		
Query 241		FGVRGIMKVCGGDDLH 256			
		FGVRGIMKVCGGDDLH			
Sbjct 257		FGVRGIMKVCGGDDLH 272			

Translation

>57F

MDTQLIIDNADVPAT
 AAATFERRSPPTGELVTRAAAASVADAIAAADSAAAAYSWSSTTGPTEERRRILLKAADLL
 EARTPEFSRVMALEVGASDLWAGVNVMLAANLFREAAALTTQIQGETIPTDKAGVLSMTV
 RQPVGVILSIAPWNGPVVLAARAIAIYPLVCGNTVVFRASELSPKTHMLIVDVLRDAGLPP
 GVLNAVNTNAPQDAPEVVDALIAHPAVRRINFTGSTRVGRVIAEKAARHLKRCLLELGGKA
 PLVXLNDADIDEAVKAAVFGAFLYQGQIXXSTERIVVXXKLLILSLPVSP

>57R

GLGKALVERFINEGGRVGVLEERSGERARELAREFGDAVEVVVGDVTLYEDNVLVVQKTVA
 RFGRLDNFVGNAGVFDFQTLPOMDAGSISRFAFDELFAVNVKAALLGAKAALAEVLKSOQ
 SLIFTVSNAGFYPPGGGPLYTASKHALVGLIRELAHELAPKVRVNGVAPGGMRTNLGGLS
 ATG TSAQTLDQVDNLEAMLVENTPLRMAPVPADYCGPYVLLASRENSGPMTGVVINTDGG
 FGVRGIMKVCGGDDLH

Sequence request

Primers to be synthesized by Retrogen

[BBa_G00100_F](#) 5'-tgccacctgacgtctaagaa-3"

[BBa_G00101_R](#) 5'- attaccgcctttgagtgagc-3"


Sample Number	Sample Description	Primer
48	DNA miniprep	BBa_G00100_F
49	DNA miniprep	BBa_G00100_F
51	DNA miniprep	BBa_G00100_F
53	DNA miniprep	BBa_G00100_F
54	DNA miniprep	BBa_G00100_F
55	DNA miniprep	BBa_G00100_F
57	DNA miniprep	BBa_G00100_F
59	DNA miniprep	BBa_G00100_F
60	DNA miniprep	BBa_G00100_F
61	DNA miniprep	BBa_G00100_F
64	DNA miniprep	BBa_G00100_F
65	DNA miniprep	BBa_G00100_F
48	DNA miniprep	BBa_G00101_R
49	DNA miniprep	BBa_G00101_R
51	DNA miniprep	BBa_G00101_R
53	DNA miniprep	BBa_G00101_R
54	DNA miniprep	BBa_G00101_R
55	DNA miniprep	BBa_G00101_R
57	DNA miniprep	BBa_G00101_R
59	DNA miniprep	BBa_G00101_R
60	DNA miniprep	BBa_G00101_R
61	DNA miniprep	BBa_G00101_R
64	DNA miniprep	BBa_G00101_R
65	DNA miniprep	BBa_G00101_R

SeqId	Sample	Primer	Date	Phred Q20	Comments
-------	--------	--------	------	-----------	----------

40776-59	48	BBaG00100-F	Aug 7 2017	phd_qual 1133 fasta	Great Sequence
40776-60	48	BBaG00100-R	Aug 7 2017	phd_qual 629 fasta	black epi, use 18723-BBa-R 40776-60 D8 REDO GC
40776-61	49	BBaG00100-F	Aug 7 2017	phd_qual 1122 fasta	Great Sequence
40776-62	49	BBaG00100-R	Aug 7 2017	phd_qual 675 fasta	black epi, use 18723-BBa-R 40776-62 F8 REDO GC
40776-63	51	BBaG00100-F	Aug 7 2017	phd_qual 1102 fasta	Great Sequence
40776-64	51	BBaG00100-R	Aug 7 2017	phd_qual 653 fasta	black epi, use 18723-BBa-R 40776-64 H8 REDO GC
40776-65	53	BBaG00100-F	Aug 7 2017	phd_qual 1121 fasta	Great Sequence
40776-66	53	BBaG00100-R	Aug 7 2017	phd_qual 480 fasta	black epi, use 18723-BBa-R 40776-66 B9 REDO GC
40776-67	54	BBaG00100-F	Aug 7 2017	phd_qual 1060 fasta	Great Sequence
40776-68	54	BBaG00100-R	Aug 7 2017	phd_qual 663 fasta	black epi, use 18723-BBa-R 40776-68 D9 REDO GC
40776-69	55	BBaG00100-F	Aug 7 2017	phd_qual 1055 fasta	Great Sequence
40776-70	55	BBaG00100-R	Aug 7 2017	phd_qual 482 fasta	black epi, use 18723-BBa-R 40776-70 F9 REDO GC
40776-71	57	BBaG00100-F	Aug 7 2017	phd_qual 1086 fasta	Great Sequence
40776-72	57	BBaG00100-R	Aug 7 2017	phd_qual 1112 fasta	black epi, use 18723-BBa-R 40776-72 H9 REDO GC
40776-73	59	BBaG00100-F	Aug 7 2017	phd_qual 1111 fasta	Great Sequence
40776-74	59	BBaG00100-R	Aug 7 2017	phd_qual 1090 fasta	black epi, use 18723-BBa-R 40776-74 B10 REDO GC
40776-75	60	BBaG00100-F	Aug 7 2017	phd_qual 1115 fasta	Great Sequence
40776-76	60	BBaG00100-R	Aug 7 2017	phd_qual 1103 fasta	black epi, use 18723-BBa-R 40776-76 D10 REDO GC
40776-77	61	BBaG00100-F	Aug 7 2017	phd_qual 1084 fasta	Great Sequence
40776-78	61	BBaG00100-R	Aug 7 2017	phd_qual 1099 fasta	black epi, use 18723-BBa-R 40776-78 F10 REDO GC
40776-79	64	BBaG00100-F	Aug 7 2017	phd_qual 1122 fasta	Great Sequence
40776-80	64	BBaG00100-R	Aug 7 2017	phd_qual 1072 fasta	black epi, use 18723-BBa-R 40776-80 H10 REDO GC
40776-81	65	BBaG00100-F	Aug 7 2017	phd_qual 1110 fasta	Great Sequence
40776-82	65	BBaG00100-R	Aug 7 2017	phd_qual 1095 fasta	black epi, use 18723-BBa-R 40776-82 B11 REDO GC

There are 24 samples.

Plasmid map



Sequence verification

Fluorene clones for deposit in registry

CCA-66, CCA-68, and CCA-70.

Sequence Alignment of clones CCA-66, CCA-68, and CCA-70.

Results:

- All 3 clones are matching the template sequence of the first gene of the fluorene pathway as shown in the alignment below performed using clustal omega. The difference between the 3 clones is the promoter (strong, weak, and moderate) but the ORF is the same (flnB).
- The forward primer BBaG00100-F worked well. The reverse primer BBaG00100-R did not work (as reported on IGEM for some sequences).
- The translation matches the flnB amino acid sequence.

Clustal Alignment

CLUSTAL O(1.2.4) multiple sequence alignment

```

68      ---AGGCGTATCACGAGGCAGAATTTTCAGATAAAAAAAAAATCCTTAGCTTTCGCTAAGGAT
66      --TAGGCGTATCACGAGGCAGAATTTTCAGATAAAAAAAAAATCCTTAGCTTTCGCTAAGGAN
70      AATAGGCGTATCACGAGGCAGAATTTTCAGATAAAAAAAAAATCCTTAGCTTTCGCTAAGGAT
          *****
68      GATTTCTGGAATTCGCGCCGCTTCTAGATTTACAGCTAGCTCAGTCCTAGGTATTATGC
    
```

```

66      GATTTCTGGAATTCGCGGCCGCTTCTAGATTGACGGCTAGCTCAGTCCTAGGTACAGTGC
70      GATTTCTGGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCCTAGGTACAATGC
          ***** ** .***** :.***

68      TAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGTGGGACTGTTGCTACCGCACGT
66      TAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGTGGGACTGTTGCTACCGCACGT
70      TAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGTGGGACTGTTGCTACCGCACGT
          *****

68      CAGCGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGGGTCGCGTA
66      CAGCGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGGGTCGCGTA
70      CAGCGCCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAGAAGTGGCGGGTCGCGTA
          *****

68      ATTGTGGTAAACAGGTGGGGCTCGCGGGATCGGTTCGCAGTTTATGTGAAGGTCTTTTACGC
66      ATTGTGGTAAACAGGTGGGGCTCGCGGGATCGGTTCGCAGTTTATGTGAAGGTCTTTTACGC
70      ATTGTGGTAAACAGGTGGGGCTCGCGGGATCGGTTCGCAGTTTATGTGAAGGTCTTTTACGC
          *****

68      GCAGGTGCCAAGGTCGTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAA
66      GCAGGTGCCAAGGTCGTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAA
70      GCAGGTGCCAAGGTCGTGGCCGCTGATTTAACCTGGGACGACGCCGATGACTTCCGCAA
          *****

68      CAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCC
66      CAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCC
70      CAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCC
          *****

68      TTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAAT
66      TTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAAT
70      TTGGACGCTGCCCGTGACGCAGTAATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAAT
          *****

68      AACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGAC
66      AACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGAC
70      AACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGAC
          *****

68      ACGACAGATCGCGACTGGGAGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCC
66      ACGACAGATCGCGACTGGGAGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCC
70      ACGACAGATCGCGACTGGGAGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCC
          *****
    
```

>Translation of 70

```

MSESGGTVATAR
QRQLVERALGEWQGEVAGRVIIVVTGGARGIGRSLCEGLLRAGAKVVAADLTWDDADDFRK
QLES DSGSMAVDMIDITDDDALDAARDAVIDRFGTVDVLVNNASLVSETLFPPTGHRNTLD
TTDRDWEVMFVNVFVGTGLKAI RRFIEPMRAQQRGSIVNVVSSGVLAVAAGGGYHGLRPWT
VEMPYQATKAAVMALT FYLAEEVVRGDGVAVNAIMPGHTRASWFDATARAFNEQGIAYFMR
PAIPEHLLPISLFLAAQESAGASGRLYYVPPXXNYDH
    
```

>66_F

```

TAGGCGTATCACGAGGCAGAAATTCAGATAAAAAAATCCTTAGCTTTCGCTAAGGANGATTTCTGGAATTCGCGGCCGCTTCTAGATTGACGGCTAGCTCAGTCCTAGGTA
CAGTGCTAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGTGGGACTGTTGCTACCGCACGTGACGCGCAGTTGGTGGAACGTGCATTGGGCGAGTGGCAAGGAG
AAGTGGCGGGTCGCGTAATTGTGTAACAGGTGGGGCTCGCGGGATCGGTTCGCAGTTTATGTGAAGGTCTTTACGCGCAGGTGCCAAGGTCGTGGCCGCTGATTTAACCT
GGGACGACGCCGATGACTTCCGCAAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCCTTGGACGCTGCCCGTGACGCAGTAA
    
```

TCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGACACGACAGATCGCGACTGG
GAGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCGATTGCTCGCTTCATCGAGCCAATGCGCGCTCAACAGCGCGGGNTCGATTGTCGACNN

>68_F

AGGCGTATCACGAGGCAGAATTTAGATAAAAAAATCCTTAGCTTCGCTAAGGATGATTTCTGGAATTCGCGGCCGCTTCTAGATTTACAGCTAGCTCAGTCTAGGTATT
ATGCTAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGAACGTGCATTGGGCGAGTGGCAAGGAGAA
GTGGCGGGTTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTTACGCGCAGGTGCCAAGGTCGTGGCCGCTGATTTAACCTGG
GACGACGCCGATGACTTCCGCAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCCTTGGACGCTGCCCGTACGCGAGTAATC
GACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGACACGACAGATCGCGACTGGG
AGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCGATTGCTCGCTTCATCGAGCCAATGCGCGCTCAACAGCGCGTTTCGATTGTCAACGTGGTAAGCAGTGGCGT
CCTTGCACTGCGAGCTGGCGGGGATACCATGGCTTGCCTGACCGTTGAGATGCCCTATCAGGCTACTAAAGCAGCTGTCATGGCTCTTACATTCTACTTGGCCGAA
AAGGTGCGCGCGATGGGGTGGCGTCAATGCTATCATGCCTGGTACACCCGCGCTTNTNGGTTTG

>70_F

AATAGGCGTATCACGAGGCAGAATTTAGATAAAAAAATCCTTAGCTTCGCTAAGGATGATTTCTGGAATTCGCGGCCGCTTCTAGATTTACGGCTAGCTCAGTCTAGGT
ACAATGCTAGCAAAGAGGAGAAAACTAGATGTCCGAATCAGGGGGTGGGACTGTTGCTACCGCACGTCAGCGCCAGTTGGTGAACGTGCATTGGGCGAGTGGCAAGGA
GAAGTGGCGGGTTCGCGTAATTGTGGTAACAGGTGGGGCTCGCGGGATCGGTGCGAGTTTATGTGAAGGTCTTTTACGCGCAGGTGCCAAGGTCGTGGCCGCTGATTTAAC
TGGGACGACGCCGATGACTTCCGCAACAATTAGAGTCCGACGGCTCTGGTATGGCCGTAGATATGGATATTACAGACGATGATGCCTTGGACGCTGCCCGTACGCGAGTA
ATCGACCGCTTCGGAACCGTTGATGTCTTGGTGAATAACGCTTCGCTGGTCTCTGAGACTTTGTTTCCACCAACGGGGCACCCTAATACCCTGGACACGACAGATCGCGACTG
GGAGGTAATGTTTGGTGTGAATGTCTTTGGAACACTTAAGGCGATTGCTCGCTTCATCGAGCCAATGCGCGCTCAACAGCGCGTTTCGATTGTCAACGTGGTAAGCAGTGGC
GTCCTTGCAGTTCGAGCTGGCGGGGATACCATGGCTTGCCTGACCGTTGAGATGCCCTATCAGGCTACTAAAGCAGCTGTCATGGCTCTTACATTCTACTTGGCCGA
AGAGGTGCGCGCGATGGGGTGGCGTCAATGCTATCATGCCTGGTACACCCGCGCTTCTTGGTTTGTGCGACCGCTCGTGCCTTAAATGAGCAGGGGATCGCATACTTC
ATGCGCCTGCTATCCGAGCACTTGCTTCTATCTCTTGTTCCTTGCAGCGCAGGAATCCGCTGGCGCTCTGGGCGTCTTACTATGTGCCCGANNGGAACACGACCAC