

# Restriction digestions

Restriction digestions was done using the following basic setup adopted from <http://nebcloner.neb.com/#!/redigest>.

Incubation time, temperature, and buffer used varied depending on the restriction enzymes used as suggested by NEB.

## Protocol

1. For each 10 uL reaction the following was mixed:

Component	10 uL reaction
10x NEB buffer	1.0 $\mu$ L
DNA	4.0 $\mu$ L
Restriction enzyme 1	0.25 $\mu$ L
Restriction enzyme 2 or H <sub>2</sub> O	0.25 $\mu$ L
Nuclease-free water	4.5 $\mu$ L

2. Subsequently samples were incubated as suggested (generally 1 hour at 37 degrees).

## Verification

Verification was done using gel electrophoresis.