

## Cloning Strategy with Restriction Enzymes

- ◆ Digestion: Vector/Backbone
  - ◇ X  $\mu\text{L}$  600 ng Vektor
  - ◇ 3  $\mu\text{L}$  10 x buffer (NEB2.1 or Cutsmart)
  - ◇ 1  $\mu\text{L}$  1. Restriction enzyme
  - ◇ 1  $\mu\text{L}$  2. Restriction enzyme
  - ◇ Add 30  $\mu\text{L}$  ddH<sub>2</sub>O
- Incubation 1 h, 37 °C
- Inactivation 10 min, 37 °C
  
- ◆ Dephosphorylation of Vector/Backbone
  - ◇ Add 2  $\mu\text{L}$  of Shrimp Alkaline Phosphatase (Jena Bioscience)
- Incubation 30 min, 37 °C
- Inactivation 5 min, 75 °C
  
- ◆ Phosphorylation of Insert
  - ◇ Add 3  $\mu\text{L}$  10x T4 Ligase Buffer (Thermo Fisher) -> contains right amount of ATP
  - ◇ Add 1  $\mu\text{L}$  PNK (NEB)
- Incubation 30 min, 37 °C
- Inactivation 20 min, 65 °C
  
- ◆ Purification of Vector and Insert
- Any given Gel and PCR Purification Kit
- We recommend purification after gel electrophoresis

- ◆ Ligation

Molar ratio: 1:5; Vector: XXX bp; Insert: XXX bp

Calculation of needed Insert in ng:

$$XX \text{ ng Insert} = 50 \text{ ng Vector} \cdot 1/5 \cdot (\text{insert (bp)} / \text{vector (bp)})$$

- ◇ X  $\mu\text{L}$  50 ng Vector
- ◇ X  $\mu\text{L}$  XX ng Insert
- ◇ 2  $\mu\text{l}$  T4 DNA Ligase Puffer (Thermo Fisher)
- ◇ 1  $\mu\text{L}$  T4 DNA Ligase (Thermo Fisher)
- ◇ Ad 30  $\mu\text{L}$  ddH<sub>2</sub>O

- ➔ Incubation 1 h, RT
- ➔ Inactivation 20 min, 80 °C

- ◆ Transformation

We recommend Heat Shock Transformation of chemically competent cells.

If you prefer transformation with electro competent cells we recommend purifying the ligation product.