

Protocol: Ligation Reactions

Author: Akshitha Yarrabothula

Date: April 8, 2011

Trinh Lab

- A. Obtain the following materials to set up the ligation reaction:
1. Vector
 2. Insert
 3. T4 DNA Ligase
 4. T4 DNA Ligase Buffer
 5. Deionized Water
- B. Measure the concentrations of the vector and insert using the nanodrop measuring device.
- C. Set up the following four reactions, making sure to use 5 ng of vector in each of the reactions for optimization:
1. Vector to insert ratio of 0:1
 2. Vector to insert ratio of 5:1
 3. Vector to insert ratio of 10:1
 4. Vector to insert ratio of 25:1

Ex: Ligation of pCT17 (pETite + lucC) with lucD

Number of base pairs of pCT17 = 3.6 kbp

Number of base pairs of lucD = 0.9 kbp

Normalization of length: $[I] / [P] = [(0.9) / (3.6)]k$

$$n_{\text{lucD}} / n_{\text{pCT17}} = 0.25k$$

$$n_{\text{lucD}} = n_{\text{pCT17}}(0.25k)$$

Using $n_{\text{pCT17}} = 5 \text{ ng}$, n_{pCT17} becomes 1.25k

k	n_{lucD}
5	6.25 ng
10	12.5 ng
25	31.25 ng

Table 1

We want to use the following amounts (in μl) of our reactants:

	V:I Ratio	V	I	Buffer	Enzyme	Water	Total
1:0	0:1	1	0	1	1	7	10
1:5	1:5	1	1	1	1	6	10
1:10	1:10	1	2	1	1	5	10