

iGEM TU/e 2015
Biomedical Engineering

Eindhoven University of Technology
Room: Ceres 0.04
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
Tel. no. +31 50 247 55 59
2015.igem.org/Team:TU_Eindhoven

Measuring Fluorescence

Table of contents

Measuring (bio)luminescence and
fluorescence

1	Preparation samples	3
1.1	Materials	3
1.2	Setup & Protocol	3
2	Measuring fluorescence	4
2.1	Materials	4
2.2	Setup & Protocol	4

1 Preparation samples

Estimated bench time: 20 minutes

Estimated total time: 20 minutes

Purpose: Prepare samples for fluorescence assay

1.1 Materials

- 384 Black wells plate
- dH₂O
- Eppendorf tubes
- Protein samples
- PCR tube

1.2 Setup & Protocol

- Make the desired protein concentrations and mixtures in a PCR tube.
- Load the samples in the 384 wells plate with a maximum of 90 µl of each sample. (The maximum loading volume may depend on the brand that is used, so search for the well-plate specifications beforehand.)

2 Measuring fluorescence

Estimated bench time: 30 minutes

Estimated total time: 30 minutes

Purpose: Analyses of the bioluminescence of the expressed proteins.

2.1 Materials

- Prepared 384 Black wells plate with loaded samples
- Tecan Safire2 plate-reader

2.2 Setup & Protocol

- For the measurement of fluorescence, open the program (Excel-file) of the Tecan Safire2 plate-reader and edit the measurement parameters.
 - Select the fluorescence measurement.
 - Plate: select the plate that is used and the wells that contain the samples.
 - Excitation wavelength: select the excitation wavelength to one suitable for the fluorophore that needs to be measured.
 - Emission wavelength: it's possible to scan over a broad wavelength range and determine the emission of the fluorophore. If you want to only measure the excitation at a certain wavelength, you need to select emission wavelength instead of emission scan.
 - Nr. of reads: select 10 (default)
 - Select endpoints measurement if you only want to perform one measurement.
- Insert the plate into the Tecan Safire2 plate-reader and start the measurement.