

# Gel Electrophoresis

Guide for visualizing PCR products, and for preparations to gel extractions and DNA purification.

## Gel Electrophoresis

---

PCR products are first run on agarose gels using gel electrophoresis where after they the DNA is extracted using DNA gel extraction.

### Creating the gel

- Use 1.0 % agarose from common stock to pour gels. If unavailable or another percentage is desired use 1.0 g Agarose per 100 mL 1x TAE for 1.0 %. Mix and heat in microwave till all agarose powder is dissolved.
- Pour liquid gel into mold and wait ~15 min for it to solidify, longer if using a lower agarose concentration.

### Running the gel

- Use a DNA ladder in first lane (max. 5  $\mu$ L).
- Mix 6x loading dye in 1:5 ratio with each sample (e.g. 5  $\mu$ L sample with 1  $\mu$ L 6x loading dye) and insert in well.
- For ~1.5 kb PCR product 20 min at 100 V is sufficient. Lower voltage = slower run time. Larger PCR product: higher run time for better separation.

### Visualizing the gel

- The gel can be visualized in the UV machine next to the door.
- Unless you are making a gel extraction, make sure to remove it from the tray, as the tray lowers the quality of the picture.
- Save the picture, print and glue it into the lab book.