

Lab Notebook Week 16

Project: NU iGEM 2017 Shared Project

Authors: Katerina Cheronis

Dates: 2017-10-03 to 2017-10-05

TUESDAY, 10/3/17

Overnight cultures started:

PeIB (2 colonies picked)

iPCR for construct correction:

pC92 and pC94

Primers used:

P139 and P140 for pC92

P140 and 141 for pC94

PCR protocol (2 tubes used)

98 deg C for 30s

Repeat 10x:

98 deg C for 15s

54.6 deg C for 30s

72 deg C for 2:20 mins

Repeat 25x:

98 deg C for 15s

65.6 deg C for 30s

72 deg C for 2:20 mins

72 deg C for 4 minutes

4 deg C for inf

Gel Results:

pC92: successful

pC94: unsuccessful

PCR for cjCas9 insertion

Plasmids	
	Plasmid
1	pK95
2	pK96

Primers		A	B	C	D
1	P144	AATTAAGAGGAGAAAAGGTCATGGGCAGCAGCCATCATC A	FW	59.8, 66.2	
2	P145	TAGTGGTGATGGTGATGATGAGGTCCGACTTTTTGAAGT C	REV	55.0, 65.6	

PCR protocol (2 tubes used)

98 deg C for 30s

Repeat 10x:

98 deg C for 15s

58.0 deg C for 30s

72 deg C for 1 mins

Repeat 25x:

98 deg C for 15s

68.6 deg C for 30s

72 deg C for 1 mins

72 deg C for 4 minutes

4 deg C for inf

Insert successfully amplified.

WEDNESDAY, 10/4/17

PCR wash (Do while gel runs)

- Add 250uL of PB and mix well with each PCR reaction
 - **Used the one in the Qiagen kit**
- Pipette everything into a column from the miniprep kit and spin for 1 minute
- Wash with CWC (500 uL) and spin for 1 minute
- Elute into a steril tube using 30uL of ddH2O

Digestion with DpnI (Prepare while gel runs)

- Add DpnI (1uL/(50uL)) to the purified PCR product and incubate for 1-4 hours (37 degC)
 - **Incubated for 60 minutes**

Nanodrop results (nanodrop after PCR wash and DpnI treatment):

Nanodrop results - His addit...		
	A	B
1	Plasmid	Concentration (ng/uL)
2	pC91 A	
3	pC91 B	
4	pC92 A	
5	pC92 B	

Ligation

Procedure:

1. Calculate the volume of the DpnI-treated PCR product containing 50-150 ng and fill in the table below.
2. Calculate the volume of water to make up to 10 μ L.
3. Mix the components below:
4. Incubate for 120 min at room temperature (~22°C).

Reaction Volumes for Ligat...					
	A	B	C	D	E
1	Reagent	pC91 A	pC91 B	pC92 A	pC92 B
2	ddH2O	3	2	3	3
3	Purified PCR product (50 ng)	5	6	5	5
4	10x reaction buffer	1	1	1	1
5	T4 DNA ligase, 5U/uL	1	1	1	1
6	Total:	10	10	10	10

Perform Transformations:

Volumes for iPCR Transformations			
	A	B	C
1	30 mL C-cells	5 uL iPCR ligated product	35 uL rescue media

PCR Wash

Introduction

This protocol outlines the procedure to perform a PCR cleanup.

Materials

- › PCR Product
 - › Take straight from PCR after removing ~3 uL to run a gel
- › 250 uL Buffer PB
 - › This is typically found in the Qiagen miniprep kit, this is a binding buffer
- › 400 uL CWC
 - › This is an ethanol wash, found in the Promega miniprep kit
- › 30 uL NF Water
- › Spin minicolumn and collection tube
- › 1.5 mL microcentrifuge tubes
- ›

Procedure

Procedure

- ✓ 1. Transfer PCR product to a microcentrifuge tube and add 250 uL buffer PB. Vortex to mix.
- ✓ 2. Pipette everything into a minicolumn and collection tube setup and spin at 14000 rpm for 1 min
- ✓ 3. Empty the collection tube and place minicolumn back into the same collection tube.
- ✓ 4. Add CWC and spin for 1 min
- ✓ 5. Transfer minicolumn to a fresh 1.5 mL microcentrifuge tube
- ✓ 6. Add 30 uL NF water. Let stand for 1 minute
- ✓ 7. Spin for 30 seconds to elute DNA

Transformation Protocol (Chelsea's)

Introduction

Transformation protocol as recommended by Chelsea Hu.

Materials

>

- > Competent cells
- > DNA
- > 2 mL microcentrifuge tube
- > Agar plates (with right antibiotic resistance if needed)
- > Ice (in bucket)
- >

Procedure

- ✓ 1. Thaw comp cells (more than 50 uL) on ice
- ✓ 2. 10 uL of competent cells in each tube (20 uL for JC8031)
- ✓ 3. Remove agar plates (containing the appropriate antibiotic) from storage at 4°C and let warm up to room temperature
- ✓ 4. .5 uL of 1 pg-10 ng of each plasmid DNA in tubes, do not disturb in anyway (not even flick)
 - a. Dilute stocks of 50 ng to 10 ng
- ✓ 5. Sit on ice for 20 mins
- ✓ 6. Heat shock at 42 C for 60 secs
- ✓ 7. Chill on ice for 5 mins
- ✓ 8. Add 10 uL (20 uL if using JC8031) of SOC (no antibiotics) and rescue for an hour at 37 C in shaker
- ✓ 9. Plate on appropriate antibiotic -agar plates