

***Chlamydomonas* transformation**

1. Preparing *Chlamydomonas* batch culture for transformation

Inoculation to produce a 4L culture:

- i. Inoculate 50 mL of TAP with a single, fresh colony from a TAP streak plate.

Use a 250mL conical flask for efficient aeration and mixing while shaking.

- i. Grow, while continuously shaking, for 5 days at 50 – 100 μ E until culture density reaches $\sim 1 \times 10^7$ cells/mL.
- ii. Transfer 50 mL into 400 mL of fresh TAP.
- iii. Grow while shaking for another 2 days at 50 – 100 μ E until culture density reaches $\sim 8\text{-}9 \times 10^6$ cells/mL.
- iv. Inoculate 4L of fresh TAP with all 450 mL *Chlamydomonas*, grow while bubble air and continuously stirring for ~ 30 hrs at 100 μ E until culture density reaches $2\text{-}4 \times 10^6$ cells/mL.

Ensure the final inoculation uses an exponential-phase culture.

For smaller volumes (up to 1L), this step can be done in a conical flask with continuous shaking.

The high volume is in case the culture does not reach the desired concentration, a larger volume of cells can be spun down for transformation.

2. Transformation

Work with open cultures should be performed in an appropriate flow hood to decrease contamination risk.

*Use 50 mL *Chlamydomonas*, at $2\text{-}4 \times 10^6$ cells/mL, per transformation.*

Each transformation requires 5×10^7 cells.

Materials to prepare

- i. TAP Sucrose, 40mM sucrose -- 15mL per transformation
- ii. Selection plates

Plate preferences depend on type of screening that will be performed

We use rectangular plates (250 mm x 250 mm)

1.5% TAP agar with appropriate antibiotic

Label on the side to avoid interference with the colony picker camera during selection

- iii. DNA insert cassettes digested with appropriate enzyme
- iv. Electroporation cuvettes
- v. Centrifuge tubes
- vi. Glass L-Shaped spreaders
- vii. Pipette tips and pipettes



Harvesting

1. Aliquot *Chlamydomonas* culture into centrifuge tubes
2. Pellet *Chlamydomonas* culture using centrifuge (4 min at 1000x g)
3. Discard supernatant and avoid disturbing cell pellet

Electroporation

1. Use TAP 40 mM sucrose to dilute cell pellet to 2×10^8 cells/mL.
2. Add 250 μ L of *Chlamydomonas* culture per electroporation cuvette.
4. Incubate at 16°C for 5 mins.
5. Prepare recovery tubes (15 mL falcon tubes) with 8mL TAP sucrose.
Steps 6-8 should be performed as quickly as is possible, especially step 8.
Before electroporation, dry the cuvette.
6. Add transformation cassette DNA (14.5ng/kb cassette per 250 μ L) and mix briefly.
Do not mix by repeatedly pipetting, use pipette tip to mix by stirring.
7. Electroporator should be set to 800V, 25 μ F.
8. Administer pulse and immediately add the mixture from the cuvette to the recovery tube.

Plating and growth

1. After all electroporations are complete, transfer recovery tubes to a shaker, and grow for 6 hours in the dark.
2. Centrifuge the recovery tubes for 4 min at 1000x g.
3. Discard most of the supernatant, leaving about 500 μ L.
4. Resuspend the pellet in the residual supernatant and spread on a transformation plate using a glass spreader.
5. After all plates are adequately dry, place them upside-down in stacks in plastic bags, and return to growth room.
6. Incubate transformation plates at 2-10 μ E for 1-2 days.
9. Increase the light intensity to 25-50 μ E for 8-12 days
This should yield >2mm colonies that can be used for selection