



dynamic

COMPOSITION

100 SHEETS / 200 PAGES

IGEM MANHATTAN

WIDE RULED

~~Procedure 1 for making~~

5/30/17

1. Streak out frozen glycerol stock of bacterial cells (TOP10, DH5 α , etc.) onto an LB plate (no antibiotics since these cells don't have a plasmid in them). Work sterile. Grow plate overnight at 37°C.

5/30/17

1. Autoclave

- 1 L LB

- 1 L of 100 mM CaCl₂

- 1 L of 100 mM MgCl₂

- 100 mL of 85 mM CaCl₂, 15% glycerol v/v

- 4 centrifuge bottles + caps

- Lots of microfuge tubes.

2. Chill overnight at 4°C

- 100 mM CaCl₂

- 100 mM MgCl₂

- 85 mM CaCl₂, 15% glycerol v/v

- Centrifuge rotor

3. Prepare starter culture of cells

- Select a single colony of E. coli from fresh LB plate and inoculate a 10 mL starter culture of LB.

Grow culture @ 37° C in shaker overnight.

Procedure I for making competent cells 5/21/17

1. Inoculate 1 L of LB media with 10 mL starter culture and grow in 37° C shaker. Measure the OD₆₀₀ every hour, then every 15-20 mins when the OD gets above 0.2.
2. When the OD₆₀₀ reaches 0.35-0.4, put the cells on ice immediately. Chill the culture for 20-30 mins, swirling occasionally to ensure even cooling. Place centrifuge bottles on ice at this time.
3. Spin #1: Pour 500 mL into ice cold centrifuge bottles. Harvest the cells by centrifugation at 3000 g (~4000 rpm) in the Beckman JA-10 rotor for 20 mins. at 4° C.
4. Decant the supernatant and gently resuspend each pellet in about 200 mL of ice cold MgCl₂. Make sure to prepare a blank bottle as a balance.

2
8

a.

5. Spin #2: Harvest the cells by centrifugation at 2000 g (~3000 rpm in the Beckman JA-10 rotor) for ~~15~~²⁰ mins. at 4°C

6. Decant the supernatant and resuspend the pellet in about 100 mL of ice cold CaCl₂. Keep this suspension on ice for at least 20 mins. Start putting 1.5 mL microfuge tubes on ice, if not already chilled.

7. Spin #3: Harvest the cells by centrifugation at 2000 g (~3000 rpm in the Beckman JA-10 rotor) for ~~15~~²⁰ mins at 4°C. Chill a 50 mL conical tube on ice.

8. Decant the supernatant and resuspend the pellet in ~25 mL of ice cold 85 mM CaCl₂, 15% glycerol. Transfer the suspension to the 50 mL conical tube

9. Spin #4: Harvest the cells by centrifugation at 1000 g (~2100 rpm

in the Beckman GH-3.8 rotor
for 18²⁰ mins at 4°C.

10. Decant the supernatant and
resuspend the pellet in 1 mL of
ice cold 85 mM CaCl₂, 15% glycerol.
The final OD₆₀₀ of the suspended
cells should be ~ 200-250.

11. Aliquot 50 μ L into sterile 1.5 mL
microfuge tubes and snap freeze
with liquid nitrogen. Store frozen
cells in the -80°C freezer.

Testing for efficiency:

06/1/17

Results:

No - no cells

low - no cells

medium - no cells

high - 10 cells (red)

$$\text{efficiency (cfu/ug)} = \frac{\text{colonies on plate (cfu)}}{\text{amount of DNA plated (ng)} \times 1000 \text{ (ng/ug)}}$$

↓
colony forming units

high - 100 pg/mL of DNA \Rightarrow total of 100 pg in transformation

10 pg of DNA was plated

$$10 \text{ pg} \left(\frac{1 \text{ ng}}{1000 \text{ pg}} \right) = 0.01 \text{ ng}$$

$$\text{efficiency} = \left(\frac{10 \text{ cfu}}{0.01 \text{ ng} \times 1000 \text{ (ng/ug)}} \right) = 1.0 \times 10^6$$

ideal efficiency: 1.5×10^8 to 6×10^8

Buffers.

10 mM KOAc ~~pt~~ (10 ml of a 1M stock)

Making of SOB (0.5L)

- 2.5 g yeast extract
- 10 g tryptone
- 0.292 g NaCl
- ~~0.093~~ 0.093 g KCl
- 1.2 g MgSO₄

0.292g
0.0928
1.2g

SOB (1L)

(PH of 7.5)

- 5 g yeast extract
- 20 g tryptone
- 0.584 g NaCl
- 0.186 g KCl
- 2.4 g MgSO₄

↳ Autoclaved the SOB

154g

Process
1)
2)
3)
4)
5)
6)
7)
8)
9)
10)
11)

06/07/17

Sample A

absorbance = 0.045
0.052
0.105
0.156

Sample B

absorbance = ~~0.044~~
0.107
0.160
0.317

Procedure:

- 1) Pick single colony of cells from SOB streaked plate.
- 2) Inoculate 2 mL SOB → grow O/N @ RT
- 3) Inoculate 250 mL SOB w/ 1 mL of O/N culture.
- 4) Grow until $OD_{600} = 0.3$ (RT)
- 5) Fill an ice bucket half way and pre-chill tubes.
- 6) Transfer cells to centrifuge bottles (flat bottom)
- 7) Centrifuge @ $3000 \times g$ @ $4^{\circ}C$ for 10 min
* ↳ (2060 rpm on our centrifuge)
- 8) Remove liquid and then resuspend cells gently in 80 mL of ice cold CCM80 buffer.
- 9) Incubate on ice for 20 min.
- 10) Centrifuge again @ $3000 \times g$ (2060 RPM) @ $4^{\circ}C$ for 10 min.
- 11) Remove liquid and resuspend pellet in 10 mL of ice cold CCM80.

pH of LB = 7.55


Orange cells = 0.15
Black cells = 0.317

X60
9353


ECO,
4561

-e
f

T



06/13/17
1B6



5E8

- 1 mL DNA
 - 1 mL of each enzyme
 - 2.5 mL buffer
 - enough H₂O to make the entire mixture
- 25 mL


1B6

- 3 mL DNA
 - 1 mL each enzyme
 - 2.5 mL buffer
 - enough water (H₂O) to make entire mixture
- 25 mL

Four Reactions

- ① EcoRI / PstI → Buffer 0
- ② EcoRI → EcoRI buffer
- ③ PstI → Buffer 0
- ④ xbaI / SpeI → tango Buffer.

Time the tubes were put in incubator = 12:50 @ 37°C



Plates 1, 2, 3 → pSB1C3
 plasmid backbones: pSB1C3, pSB1A3, pSB1K3, pSB1T3

RFP gene: BBa_j04450 in pSB1C3
 ↳ Plate 1 → well 230

Plaque #7

RFP BBa_j04450

Spec. Res.

Spec. Res.	Plasmid Backbone	Plate	Well
chloramphenicol	pSB1C3	1	230
ampicillin	pSB1A3	4	2H → single error
kanamycin	pSB1K3	4	6B
tetracycline	pSB1T3	4	5B → single error

Date: 06/20/17

- Used sputter machine to sputter gold onto a glass slide until it was 20.1 nm thick.

- Made two 0.1 M solutions of $\text{Cl}_2\text{Ni} \cdot 6\text{H}_2\text{O}$ (200 mL) and NaBH_4 (200 mL).

- Put membrane in DI water and sonicated water.

- Put polycarbonate membrane (7-10 microns) thick in the u-tube and then filled each of the sides w/ each of the solutions simultaneously.

- Waited for reaction for 15 min.

Then repeated procedure but w/ anodized aluminum oxide (60 microns thick) membrane.

And waited 15 min.

06/22/17

- Made .01 M solution (200 mL) of $\text{CaCl}_2 \cdot 6\text{H}_2\text{O}$
and .01 M solution (200 mL) of NaBr .

- Grew nanowires with .01 M solutions and another
membrane w/ .01 M solutions.

- Membrane w/ ~~0.01 M~~ ^{0.1 M} solutions ran for 30 min
and conductivity was good .24 (ohm)

- Dissolved the polycarbonate membrane (^{0.01 M} 30 min)
with *methylene chloride, and centrifuged
at 4000 rpm for 7 min and discarded methylene chloride
two times. And repeated

- *Methylene chloride was discarded into waste
beaker after centrifuge

06/28/17

Digest with EcoRI / PstI

- ① #2 DNA (100 ng)
 - ② #7 DNA (100 ng)
- } in separate tubes

Digest w/ EcoRI / PstI

- ③ PCDF duet plasmid (1100 ng total)
- ④ linearized pSBIC3 (100 ng total)

1/2

Digest w/ EcoRI / PstI (#2 DNA / #7 DNA)

lost

- ↳ 10 μ L DNA
 - ↳ 2.5 μ L buffer 0
 - ↳ 0.5 μ L enzyme 1
 - ↳ 0.5 μ L enzyme 2
 - ↳ 11.5 μ L H₂O
- } total of 25 μ L in rxn

Digest w/ EcoRI / PstI

- ③ PCDF duet plasmid (1100 ng)
 - ↳ 2 μ L of plasmid
 - ↳ 0.5 μ L of enzyme 1
 - ↳ 0.5 μ L of enzyme 2
 - ↳ 12.5 μ L of buffer 0
 - ↳ 9.5 μ L water
- } total of 25 μ L in rxn

Digest w/ EcoRI / PstI

④ linearized pSBIC3 (100 ng total)

↳ 4 μL of pSBIC3

↳ 0.5 μL of enzyme 1

↳ 0.5 μL of enzyme 2

↳ 6.25 μL of buffer 0

↳ 13.75 μL H₂O

total of 25 μL in rxn

Last ←

↳ PstI → not sensitive to heat inactivation
therefore: cleanup kit.

General DNA cleanup from enzymatic reactions protocol

1. Adjust the volume of the reaction mixture to 200 μL w/ water, nuclease-free or TE buffer (not included).
2. Add 100 μL of Binding Buffer. Mix thoroughly by pipetting.
3. Add 300 μL of ethanol (96-100%) and mix by pipetting.
4. Transfer the mixture to the DNA Purification Micro Column preassembled w/ a collection tube. Centrifuge the column for 30-60 seconds at 14,000 \times g. Discard the flow through. Place the DNA Purification Micro Column back into the collection tube.

5. Add 700 μ L of Wash Buffer (supplemented w/ ethanol) to the DNA Purification Micro Column and centrifuge for 30-60 seconds at $14,000 \times g$. Discard the flow-through and place the purification column back into the collection tube.

6. Repeat step 5

7. Centrifuge the empty DNA Purification Micro Column for an additional 1 minute at $14,000 \times g$ to completely remove residual Wash Buffer.

Note: This step is essential to avoid residual ethanol in the purified DNA so the presence of ethanol in the DNA sample may ~~not~~ inhibit downstream enzymatic rxns.

8. Transfer the DNA Purification Micro Column into a clean 7.5 mL microcentrifuge tube.

9. Add 10 μ L of water to the center of the DNA Purification Micro Column membrane. Centrifuge for 1 minute at $14,000 \times g$ to elute DNA.

10. Discard the purification column and store the purified DNA @ $-20^{\circ}C$.

Electrophoresis

PDF
EORL / stl
whole sample. ladder.
5 μ

7/7/17

o Digest of insert fragments

1) iGBM-2017-2 } Both @ 10 ng/mL
iGBM-2017-7 }

- we want to digest 200 ng total in a final volume of 25 μ L

2) $200 \text{ ng} \times \frac{1 \mu\text{L}}{10 \text{ ng}} = 20 \mu\text{L}$ of each fragment needed
per rxn

3) Buffer (Buffer 2) is 10x concentration:

25 μ L final volume $\rightarrow \frac{25}{10} = 2.5 \mu\text{L}$
needed in rxn

4) 0.5 μ L of each enzyme needed \rightarrow 0.5 μ L PstI

0.5 μ L EcoRI-HF

1.0 μ L Total

5) H₂O added to make the final total volume equal to
25 μ L

EXN # 1

20 μ L iGEM-2017-2

2.5 μ L Buffer 2

1.5 μ L H₂O

0.5 μ L PstI

0.5 μ L EcoRI-HF

EXN # 2

20 μ L iGEM-2017-7

2.5 μ L Buffer 2

1.5 μ L H₂O

0.5 μ L PstI

0.5 μ L EcoRI-HF

- 6) Incubate at 37°C for 60 min.
- 7) Incubate at 80°C for 20 min.

PCDF-duet plasmid digest

- digest a total of 1000 ng in a final volume of 20 μ L

- PCDF-duet @ 150 ng/ μ L

1) $1000 \text{ ng} \times \frac{1 \text{ mL}}{150 \text{ ng}} = 6.7 \text{ mL per row}$

2) Buffer (Buffer 2) is 10x concentration:

- 20 μ L final volume $\rightarrow \frac{20}{10} = 2.0 \mu\text{L}$ needed

3) 0.5 μ L of each enzyme: 0.5 μ L PstI

0.5 μ L EcoRI-HF

1.0 μ L Total

4) H₂O added to make the final volume equal to 20 μ L

PND

6.7 μ L pCDF-duet

2.0 μ L Buffer 2

10.3 μ L H₂O

0.5 μ L PstI

0.5 μ L EcoRI - HF

5) Incubate for 60 min @ 37°C

6) Incubate for 20 min @ 80°C

7) Run DNA on gel for electrophoresis separation (45-60 min)

07/11/17.

- Followed Procedure 1 for making competent cells.
with BL21 and DH5 α

Steps done for:


BL21
1 3 5 7 9
2 4 6 8

DH5 α
1 3 5 7 9
2 4 6 8

Microfuge Tubes

BL21 = X

DH5 α = O



competent cell
Transformation
with Expression
vectors

50
30 μ L DNA in each

07/13/17

2 μ L from wells

Ice for 30 min.

In incubator for 45 sec @ 42°C

Ice for 2 min

put in 500 μ L LB.

then in incubator (37°C) for 1 hr.

then put contents from microfuge tubes into chloramphenicol

- 1 = BBa-K1316012 Plate 5-7E
- 2 = ~~Ba~~ BBa-K608002 Plate 1-30
- 3 = BBa-I20270 Plate 7-21B.
- 4 = BBa-R0040 Plate 7-21D
- 5 = BBa-J364000 Plate 7-21F
- 6 = BBa-J364001 Plate 7-21H
- 7 = BBa-J364002 Plate 7-21J
- 8 = BBa-J364003 Plate 7-~~21I~~ 21L
- 9 = BBa-J364004 Plate 7-21N
- 10 = BBa-J364005 Plate 7-21P

RED Kit

7/18/17
took 6uL from DNA and put it into separate tubes.

Rxn Master \rightarrow Digest
Super mix

25 mL Total

2.5 mL buffer $\times 6 = 15$ mL \rightarrow 15 uL buffer

6 mL DNA $\times 6 = 36$ mL DNA + 1 uL PstI

16.5 mL H₂O $\times 6 = 99$ mL H₂O \rightarrow 1 uL EcoRI
97 uL H₂O

25

114 mL

6

19 uL
of master
mix in
each

Incubate at 37°C for 1 hour

Made agarose gel. (1%)

\hookrightarrow [50 mL TBE
with 0.5g agarose DNA]

"silver stuff"

- add red dye to gel (5 μ L)

~~Greening dye~~

Loading Gel (~~10 μ L~~)

3 μ L glycerol

25 ~~75~~ mg bromophenol

dH_2O to 10 mL

- 5 μ L of dye in DNA.

adder	D2	D2	D2	D2	D2	D2	D7	D7	D7	D7
	1	2	3	4	5	6	1	2	3	4
									D7	D7
									5	6

Chosen tube: D7-3

\rightarrow Used 1 μ L, 2 μ L, 3 μ L for the

transformation procedure from 07/13/17.

- Transferred cells from Spec (~~Dr. Bryan~~
~~and David~~) D2 and B2 into plates
(into culture tubes 5 mL LB and 5 μ L of
spec)

7/19/17
↳ DNA extraction using GeneJET Plasmid
Miniprep kit.

from D2 and B2

↳ Cut with $PstI$ / $ECORI$

→ β = 1

RAN gel for 1:20

OR 1000

Dilution: 2.14

7/20/17

↳ Absorbance of LB/spec w/ BL21 cells /

PCDF- iGem-7#3 spec I. was .575

~~After 20 min absorbance:~~

Started w/ OD_{600} of .75

① ↳ took out 1 ml out of LB/spec w/ BL21
cells PCDF- iGem-7#3 and labeled it
microfuge tube (X) negative, put it in
iGem box:

② ↳ ~~Mass~~ Put 40 μ l of \pm 1M IPTG in
flask w/ rest of LB/spec w/ BL21 cells

③ ↳ we let it incubate @ 37°C for 1 hr

④ Took 1 ml from step 3. Spun down cells and put tube in Igem box.

⑤ Repeated steps 3-4 for after 2 hrs.

- Used Thermo Scientific Gene Jet Gel extraction kit. to extract DNA

1 - OD₂₆₀ → 0.037 A 37 ng/μl

2 - OD₂₆₀ → 0.045 A 45 ng/μl

ratio:

Ligation

1: 40 ng vector

3 58.4 ng insert

1 μL vector DNA

6 μL insert (iGEM-X)

1 μL T4 DNA Ligase Buffer.

1 μL H₂O

1 μL T4 DNA ligase

10 μL total.
RT 30 min.

→ transform into competent cells

07/25/17

- 2 500 mL of LB w/ spec in them and ^{iGem-7} PCDF
and 1 1000 mL of LB w/ spec ^{iGem-2} PCDF were
incubated until an absorbance of

Gel electrophoresis: w/ PCDF dnet E/P
Comparison of cut vs uncut DNA

ladder

PCDF
dnet
E/P

control

- incubated another PCDF dnet E/P at 80 °C.

EcoRI / PstI

Weight of DNA gel slice: .05g

- Follow protocol A. DNA extraction from the gel
using centrifuge from Thermo Scientific GeneJET
GEL Extraction Kit. until step 9.

Procedure:

1. Excise gel slice containing DNA fragments. ~~use~~

Record weight of the gel slice. .05g.

2. Add 1:1 volume of Binding Buffer. 100ul

3. Incubate gel mixture @ 50-60 °C for 10 min.

Vortex gel mixture.

nanopip
time
of plasmid

4. Transfer up to 800 μL of gel solution to the GeneJET purification column. Centrifuge for 1 min. Discard flow-through.
5. Add 700 μL of Wash Buffer to the column. Spin for 1 min. Discard flow through.
6. Centrifuge again for 1 min and Discard flow through.
7. Add 15 μL of H_2O ^{\rightarrow wait 1 min}. Centrifuge for 3 min at 8000 rpm and collect in microfuge tube. \rightarrow wait 1 min.
8. Repeat w/ 10 μL of H_2O .
9. Discard the GeneJET purification column and store DNA at -20°C .

4 sets of α .

- ① Negative
- ② Control PCDF 40 ng / μL
- ③ Purified DNA from Gel
- ④ Heated @ 80°C .

1 mL of Control
1 mL of treated
4 mL of Purified Cell.

Transformation to Competent Cells

Ice for 20 min.

Incubate for 45 sec @ 42°C

Ice for 2 min.

put in 500 μ L LB.

Incubate (37°C) for 1 hr.

Then put contents from microfuge tubes onto spec. plates.

wait
(min)

7/26/17
Ligation w/ insert 2, 8, 9, and a
negative control.

1 40 ng/ml vector

58.4 ng/ml insert

1 μ L vector DNA

6 μ L insert (igem-7)

1 μ L T4 DNA Ligase Buffer

1 μ L H₂O

1 μ L T4 DNA ligase

Total: 10 μ L

Transformation to competent cells:

Ice for 20 min

Incubate for 45 sec @ 42°C

Ice for 2 min

put in 500 μ L of LB

Incubate (37°C) for 1 hour

Then put contents from microfuge tubes

onto spec plates.

50ml

7/27/17

↳ SDS-PAGE on iGem-7 PCDF w/ IPTG
and iGem-7 Neg.

↳ Put 100 μ L of SDS-PAGE buffer (blue buffer) in each tube
and boil @ 95°C for 30 min.

Making Gold Anode

50 mM

and 50 mM

1 M ~~50~~ Sodium borohydride

NaBH_4

MW: 37.83 g/mol

$$\left(\frac{.1 \text{ mol}}{\text{L}} \right) \left(.05 \text{ L} \right) = .005 \text{ mol}$$

$$\left(\frac{37.83 \text{ g}}{\text{mol}} \right) \left(.005 \text{ mol} \right) = .189 \text{ g}$$

.189 g NaBH_4

.05 M tetrachloroauric
(III) acid
trihydrate

$\text{AuCl}_4 \cdot 3\text{H}_2\text{O}$

MW: 393.83 g/mol

$$\left(\frac{.05 \text{ mol}}{\text{L}} \right) \left(.05 \text{ L} \right) = .0025 \text{ mol}$$

$$\left(\frac{393.83 \text{ g}}{\text{mol}} \right) \left(.0025 \text{ mol} \right)$$

= .984 g

$\text{HAuCl}_4 \cdot 3\text{H}_2\text{O}$

~~.005~~ mM NaBH_4

.0189 g NaBH_4
50 mL

~~.0025~~ mM $\text{HAuCl}_4 \cdot 3\text{H}_2\text{O}$

.0984 g
50 mL

- Set up u-tube w/

and .01 M and ~~.005 M~~ NaBH_4 solution

and .005 M $\text{HAuCl}_4 \cdot 3\text{H}_2\text{O}$ membrane

* part 1 (2000) *

O2 DNA

- 1: ~~017~~ .036
- 2: ~~023~~ .020
- 3: ~~024~~ .043
- 4: ~~016~~ .040
- 5: ~~034~~ .030
- 6: .015

25 μL H ₂ O	60
30 μL buffer 2	
1 μL Rst+1	
1 μL EcoRI	

- 7: .020
- 8: .012
- 9: .032
- 10: .020

10 → 25 μL = 250 μL → 12
 12 x 25 = 300 μL total
 - 12 x 20 = 240

60 μL total
30 μL buffer
30 μL
1 μL enzyme
1 μL enzyme
28 μL - 120

- Extraction DNA from cells from plate 9.
using the RapidPURE Plasmid Mini Kit
Protocol

1. Prepare Wash Solution by adding 100% ethanol
2. Spin 1.5 mL of overnight culture, discard supernatant.
3. Resuspend cell pellet in 50 μ L Pre-Lysis Buffer
4. Add 100 μ L Alkaline Lysis Solution, invert to mix.
5. Incubate @ room temp. 3 mins.
6. Add 100 μ L Neutralizing Sol'n + incubate 1 min. @ room temp. mixing or vortexing gently several times.
7. Microfuge 2 mins. @ room temp.
8. Transfer supernatant (cleared lysate) to new tube.
9. Add 250 μ L RapidPURE Mini Salt Sol'n to cleared lysate.
10. Transfer mixture to assembled RapidPURE Turbo cartridge.

11. Spin 30 sec. @ $\leq 14000 \times g$.
12. Add 350 μL prepared Wash Sol'n.
13. Spin 30 sec. @ $\leq 14000 \times g$.
14. Discard a flow through.
15. Spin 2 mins.
16. Insert Turbo Filter into a clean recovery tube.
17. Add 100 μL water or TE to filter membrane.
18. Incubate 5 mins. @ room temp.
19. Spin 1 min. to collect DNA.

$$\frac{\text{Igem 2}}{\text{WT}} = .94$$

$$\text{Igem 7} = .82$$

- Cells for PCDF-igem-7 and PCDF-igem-2
(W6) were grown in LB/spec @
30°C and samples were taken every
hour.

Making gold nanowires

$\text{HAuCl}_4 \cdot 3\text{H}_2\text{O}$ 0.005 M (100 mL)

NaBH_4

0.0189 g NaBH_4 in 100 mL H_2O

Igen 9

#8: 0.103

#6: 0.136

#3: 0.098

#1: 0.075

#2: 0.130

#5: 0.097

#7: 0.155

#4: 0.137

Igen - 2

#3: 0.095

black dot ✓
 iGem gel = left → top
 Brian gel = right → bottom
 8/2/17

STS - PAGE

- took igem-plate ~~2 w/ IPTG~~ and iGem-7 samples w/o IPTG and ho
 and added the blue buffer (100 μ L) S20
-w.
- heated for 20 min @ 92°C
- Spin them down.

Gel ladder

Gel order (Expression)
 for STS - PAGE

ladder	iGem-7	iGem-7 w/ IPTG 1 hr	iGem-7 w/ IPTG 2 hr	iGem-7 w/ IPTG 3 hr	iGem-7 w/ IPTG final	iGem-2 w/	iGem-2 w/ IPTG 1 hr	iGem-2 w/ IPTG 2 hr	iGem-2 w/ IPTG 3 hr

① Western Blot.

② Coomassie Stain

100 V
 Time: 1:30
 Amp: 1.00A
 WAT: 300W

ladder

IPTG and hourly samples w/ IPTG

Digest for DNA extractions from plate - 9

Master Mix

2.5 mL buffer x 9 = 22.5 mL	→ 22.5 mL buff
6 mL DNA x 9 = 54 mL DNA	1 mL P ₅ FI
16.5 mL H ₂ O x 9 = 148.5 mL H ₂ O	1 mL EcoRI
	146.5 mL H ₂ O
	171 mL

9
= 19 mL of master mix in each

Gem-9
w/ IPTG
3 hrs
2
1

- Incubate @ 37°C for 1 hr.

- make agarose gel (17.)

↳ (50 mL TBE)

- add red dye (0.5g agarose DNA) (Suber safe) to gel. (5 mL)

- 5 mL of dye in each tube.

Gel electrophoresis (digest)

ladder

iGem-9 #1	iGem-9 #2	iGem-9 #3	iGem-9 #4	iGem-9 #5	iGem-9 #6	iGem-9 #7	iGem-9 #8	iGem-9 #9
-----------	-----------	-----------	-----------	-----------	-----------	-----------	-----------	-----------

08/03/17

Run samples for Coomassie on 12% gel

- dilute .005 M HAuCl_4 by half

- not do @ equal concentrations.

if it still contamination dilute NaBH_4 by half

NaBH_4
.005 M.

MW: 37.83 g/mol

$$\left(\frac{.005 \text{ mol}}{K} \right) \left(\frac{1 \text{ L}}{K} \right) = .0005 \text{ mol}$$

$$\frac{37.83 \text{ g}}{\text{mol}} \left(.0005 \text{ mol} \right)$$

$$= .0189$$

$\text{HAuCl}_4 \cdot 3\text{H}_2\text{O}$
.005 M.

MW: 393.83 g/mol

$$\left(\frac{.005 \text{ mol}}{K} \right) \left(\frac{1 \text{ L}}{K} \right) = .0005 \text{ mol}$$

$$\frac{393.83 \text{ g}}{\text{mol}} \left(.0005 \text{ mol} \right)$$

$$= .196$$

ladder

Try: NaBH_4 .005 M

and HAuCl_4 .0025 M


- List of things to do 08/15/17

- ✓ ① Running SDS Page for verification of test expression for iGem-7 / iGem-2
- ② Same gel for expression # 7
- ✓ ③ Pick ^{cells} from plate # 9, and let grow overnight.
- ④ Re-do # 8, # 9.
- ✓ ⑤ Glycerol stock of # 7, # 2

15m SDS PAGE


adder	iGem-2 wt	iGem-2 w/ IPTG 1hr	iGem-2 wt 2hr	iGem-2 3hr	iGem-2 final	iGem-7	iGem-7 w/ IPTG 1hr	iGem-7 IPTG 2hr	iGem-7 IPTG 3hr	iGem-7 final
-------	--------------	--------------------------	---------------------	---------------	-----------------	--------	--------------------------	-----------------------	-----------------------	-----------------

Making of Au nanowires :



• 0.0025 M

• 0.00945 g in
100 mL
H₂O



• 0.0025 M

Dic

08/16/17

- Isolating DNA from plates heat 9, - Cut 8,
Cut ⊖.

Ran out of pre-lysis buffer:

Following has pre-lysis buffer and followed
RapidPURE-Plasmid Mini.lit protocol:

Heat 9: 2-6

Cut 8: 1-6

Used different neutralization solution for heat 9-6

Control Digest: PCDF iGem-2 # 3

Digest for DNA extraction from Plates heat 9, cut 8,
cut ⊖

2.5 mL buffer × 19 = 47.5 mL buffer

6 mL DNA × 19 = 114 mL DNA.

16.5 mL H₂O × 19 = 313.5 mL H₂O

47.5 mL buffer

1 mL Pst I

1 mL EcoR I

313.5 mL H₂O

361 mL

19

- Incubate @ 37°C for 1 hour.

19 μL of
master mix
in each

08/17/17
- Gel electrophoresis for DNA from
heat 9, cut 8, cut \ominus plates

00544A Gel 1

<u>ladder</u>	<u>control</u>	<u>control</u>	<u>cut 8</u> 1	<u>cut 8</u> 2	<u>cut 8</u> 3	<u>cut 8</u> 4	<u>cut 8</u> 5	<u>cut 8</u> 6	<u>cut 8</u> 7
<u>cut 8</u> 5	<u>cut 8</u> 6	<u>Heat 9</u> 1	<u>Heat 9</u> 2	<u>Heat 9</u> 3					

ladder control Gel 2
cut \ominus 6 Heat 9 4 Heat 9 5 Heat 9 6 cut \ominus 1 cut \ominus 2 cut \ominus 3 cut \ominus 4 cut \ominus 5

Heat 9 - 4 looks like it worked.

Au nanowires

.005 M HAuCl₄ w/ .0025 M NaBH₄
200 nm membrane

- 8 and 9

cut 8
5

at 8
5

08/17

- Ligation of PCDF - Duet naked -

Digests

PCDF = 15 μ L total

①

last

- 10 μ L DNA
- 1.5 μ L Buffer 2
- 1 μ L Pst I
- 1 μ L Eco RI
- 1.5 μ L H₂O

iGem 8

Same

iGem 9 - 25 μ L total

- 20 μ L DNA
- 2.5 μ L Buffer 2
- 1 μ L Pst I
- 1 μ L Eco RI
- .5 μ L H₂O

iGem 8/PCDF - 35 μ L total

- 10 PCDF
- 20 iGem 8
- 30 μ L DNA
- 3.5 μ L buffer
- 1 μ L Pst I
- 1 μ L Eco RI

- Incubate for 2hr

PCDF vector: iGem-# insert.

Ligation:

Total = (10 μ L)

1: 3

- ✓ 1 μ L PCDF
 - ✓ 7 μ L insert iGem-#
 - ✓ 1 μ L buffer
 - ✓ 1 μ L ligase
- last

- neg. control
- ✓ 7 μ L H₂O
 - ✓ 1 μ L PCDF
 - ✓ 1 μ L buffer
 - ✓ 1 μ L ligase → last

✓ Reminder

- Run Gel for expressions

Ligation for PCDF/i

Control	total
✓ - 3.5 mL PCDF	10 mL.
✓ - 1 mL buffer	
✓ - 1 mL ligase	
✓ - <u>4.5</u> H ₂ O	

✓ 8 mL of tube w/ PCDF & Gem-8
✓ 1 mL buffer
✓ 1 mL ligase

for

Digest

09/05/17

PCDF - 20 μ L total

10 μ L DNA

2 μ L Buffer 2

lost } 1 μ L Pst I

1 μ L Eco RI

6 μ L H₂O

@ 37°C for 1 hr.

then 80°C for 20 min.

agrose gel

Control - 20 μ L

rec. tube
6 μ L DNA

2 μ L Buffer 2

lost } 1 μ L Pst I

1 μ L Eco RI

10 μ L H₂O

0.126 A F6EM2

0.102 A I6EM7

- 3 times → 3 tubes

09/21/17

- combine @ step 5

DNA extraction from gel using centrifuge
for PCDF/EcoRI/
PstI

Digest of 8 & 9: (total = 15 μ L)

10 μ L DNA

1.5 μ L Buffer 2

last { 1 μ L PstI

{ 1 μ L EcoRI

1.5 μ L H₂O

borate

Gold nanowires

09/27/17

• 0.0025 M HAuCl_4

• 0.005 M NaBH_4

↳ good results

• 0.0025 M HAuCl_4

• 0.05 M NaBH_4

↳ too high [NaBH_4]

• 0.0025 M HAuCl_4

• 0.0025 M NaBH_4

9/28/17

• 0.0025 M HAuCl_4

• 0.005 M NaBH_4

↳ 0.189 g

10/2/17

1) EcoRV XBTB Digest of Tgen⁹
And KCl₃
↳ Successful (presence of gel)
SAVED

2) B121 transformation using - KCl₃
Tgen #9 sample 1

Protocol followed same as

7/13/17

Except 1ul of DNA was
used (higher concentration)
SAMPLE

10/2/17

1) Eosin YSTB Digest of Igem 9
And K13

↳ Successful (picture of gel)
SAVED

↳ last two

2) B121 Handwritten using K13

Igem #9 sample 1

Protocol followed same as

7/13/17

Except 1ul of DNA was
used (higher concentration)
SAMPLE

10/4/17 ^{50mL}
- Cell culture : 500 mL of (LB & cells)
450 μ L of spec

Absorbance = .301

- Waited 45 min then saw absorbance.

Absorbance = .905

- Save 1 mL of it

- record time of putting in 500 μ L IPTG

↓
12:55 p.m.

tube

✓
✓
✓
✓

BioBricks

* Digest linearized pSB1C3 as outlined on iGEM page:

Master Mix:
5 μ L NEB 2
5 μ L BSA - (2 mg/mL stock)
0.5 μ L EcoRI - HF (NEB)
0.5 μ L PstI (NEB)
0.5 μ L PpnI (Thermo)
13.5 μ L H₂O

25 μ L total.

Digest:
4 μ L pSB1C3 (linearized)
(25 μ g/ μ L, 100 μ g total)
+ 4 μ L enzyme mix
8 μ L total (12.5 μ g/ μ L)

37° 45 min 1 DT # 2, 8, 9, 7 digested w/ PstI / EcoRI
80° 20 min # 2 @ 8 μ g/ μ L.
7, 8, 9 @ 6.6 μ g/ μ L

Ligation: (Transform 2 μ L)
Red colonies no good! 25 ng of # 2 # 7, 8, 9
3.1 μ L 3.8 μ L

Ligation #2

✓ 2 μ L vector
✓ 3.1 μ L insert
✓ 1 μ L Buffer
✓ 0.5 μ L Ligase
✓ 3.4 μ L water

10 μ L

2 μ L pSB1C3 (25 μ g total)
25 ng total insert
1 μ L T4 Buffer
0.5 μ L T4 ligase
H₂O to 10 μ L
16° C 30 min, 80° C 20 min

for Mtr-CAB

10/9/17.

DNA extraction ↑ from GEL following protocol

A from THERMO SCIENTIFIC GeneJET ~~is~~

a few steps were changed. we followed the
exact procedure from 07/25/17.

00260 Absorbance: 0.642

45 μ L H₂O