

## DNA Electrophoresis

|                  |  |
|------------------|--|
| <b>Goal</b>      | To separate DNA fragments according to their size  |
| <b>Materials</b> | <ul style="list-style-type: none"><li>• DNA samples</li><li>• Appropriate DNA ladder</li><li>• 6X loading dye</li><li>• Pipette tips</li><li>• Agarose</li><li>• Microwave</li><li>• Erlenmeyer flask</li><li>• Graduated Cylinder</li><li>• 1X TAE buffer</li><li>• Gel box</li><li>• Gel casting tray</li><li>• Ethidium bromide</li><li>• UV gel imager</li><li>• Staining container</li></ul>  |
| <b>Procedure</b> | <ol style="list-style-type: none"><li>1. Measure 30 mL of 1X TAE buffer and put in Erlenmeyer flask</li><li>2. Weight 0.24g of agarose</li><li>3. Add to TAE buffer in Erlenmeyer flask</li><li>4. Microwave flask to dissolve agarose. Swirl around to ensure that all agarose has been dissolved</li><li>5. Pour melted agarose into gel casting tray and carefully insert the comb. Let the gel sit for 20 minutes. After, carefully remove the comb from gel</li><li>6. Place gel into gel box with wells nearest the black electrode</li><li>7. Add 1X TAE buffer to gel box. Pour enough buffer so there is approximately 1cm of buffer overtop of gel</li><li>8. Prepare desired samples for loading. Add 6X loading dye to each sample</li><li>9. Load the ladder and desired samples into each well</li></ol> |

|  |  |
|--|--|
|  | <ol style="list-style-type: none"><li>10. Close lid for gel box and make sure to match electrodes properly</li><li>11. Turn on power and run at 80V. Run until DNA samples have run to an appropriate location on gel</li><li>12. Remove gel, put in staining container and fill with 1X TAE buffer until staining is needed</li><li>13. To stain, add 2 <math>\mu</math>L of ethidium bromide to buffer in staining container. Gently swirl container to mix ethidium bromide. Let gel stain sit for around 20 minutes</li><li>14. Carefully remove gel from staining container and place onto an appropriate platform</li><li>15. Turn on UV light and visualize/analyze gel</li></ol> |
|--|--|