

Lab Notebook - Week 14 (9/11/2017-9/28/2017)

Project: NU iGEM 2017 Shared Project

Authors: Katerina Cheronis

Dates: 2017-09-11 to 2017-09-17

MONDAY, 9/11/17

things are going pBAD today

Sequencing Sent:

	A	B	C
1	Name	Plasmid	Primer
2	iPY1	pC86	P2
3	iPY2	pC87	P2
4	iPY3	pC88	P2

Nanodrop results

PK90: pet28a

37.8 ng/uL

260/280: 1.91

260/230: 1.99

PC91: pBAD

9.9 ng/uL or 40 ng/uL

260/280: 1.90

260/230: 1.48

TUESDAY, 9/12/17

NEW DAY GOOD DAY PRODUCTIVE DAY

cjCas9 gBlock into a random vector

Selected random vector: YcdB_saCas9_His6 (pSB1C3)

Phosphorylation:

- Treat the DNA with T4 polynucleotide kinase (PNK) to phosphorylate the 5'-ends
 - 10x T4 DNA ligase buffer (1 μ L)
 - linearized DNA (2 μ L)
 - T4 PNK (1 μ L)
 - ddH₂O (6 μ L)
- Incubate at 37 °C for 30 min, then cool to room temperature.

Single digest:

1. Add components to a clean tube in the order shown:
 - 1 μ L DNA (concentration 1 μ g/ μ L)

- 2 μL 10x buffer
 - 1 μL restriction enzyme
 - 16 μL sterile water
2. Incubate the reaction at digestion temperature (usually 37°C) for 1 hour.
 3. Stop the digestion by heat inactivation (65°C for 15 minutes) or addition of 10mM final concentration EDTA.
 4. The digested DNA is ready for use in research applications.

Ligation:

- Add the following to the reaction mix:
 - (Total Volume: 20 μL)
 - ddH₂O: 8 μL
 - 10x T4 DNA ligase buffer: 1 μL
 - T4 DNA ligase: 1 μL
- Incubate at room temperature for 60 min.

Gel electrophoresis:

- Run gel for linearized plasmid and ligated product with a 2-log ladder
-

Double Restriction Digest - cjCas9 into pET28a Expression Vector

Plasmid goal:



Restriction sites: BamHI and NotI

Digestion Setup		
	A	B
1	Insert (cjCas9)	Vector (pET28a)
2	100 ng DNA (10 μL)	600 ng DNA (15 μL)
3	10x Buffer (3 μL)	10x Buffer (3 μL)
4	NF Water (15 μL)	NF Water (15 μL)
5	BamHI (1 μL)	BamHI (1 μL)
6	NotI (1 μL)	NotI (1 μL)

- Mix above components into separate 1.5 mL microcentrifuge tubes
- Incubate at 37 C for 1 hour
- Heat inactivate at 65 C for 15 mins

WEDNESDAY, 9/13/17

ClyA-saCas9-His6 PCR

INSERT PCR:

Plasmid: pBAD (plasmid pC37 and pC38)

Primers

PCR Primers - Insert			
	A	B	C
1	Primer Name	Sequence	Melting Temperature
2	P77 (FW)	gagcacgaaagaggggacaaTTAAAGAGGAGAAAGGTCATGA	51.7, 64.5
3	P78 (REV)	agaatgtaattctctcatGACTTCAGGTACCTCAAAGA	56.4, 65.9

PCR Materials:

- nf water 22 uL
- 1 uL 10uM P77
- 1 uL 10uM P78
- 1 uL plasmid (~1 ng/uL)
- 25 uL Phusion HS Flex 2xmm
- PCR tube

Total volume = 50 uL

PCR protocol (3 tubes used)

98 deg C for 30s

Repeat 10x:

98 deg C for 15s

54.7 deg C for 30s

72 deg C for 30s

Repeat 25x:

98 deg C for 15s

67.5 deg C for 30s

72 deg C for 30s

72 deg C for 5 minutes

4 deg C for inf

Gel of pC37 and pC38:

pC38 worked the best

BACKBONE PCR:

Plasmid: try pC59 and pC74 --> dilute to 1ng/uL

Primers

Table2			
	A	B	C
1	Primer Name	Sequence	Melting Temperature
2	P79 (FW)	TCTTTGAGGTACCTGAAGTCatgaagagaaattacattctggga	50.6, 66.6
3	P80 (REV)	ATGACCTTTCTCCTCTTTAAttgtcccctcttctgctc	50.4, 61.7

PCR Materials:

- nf water 22 uL
- 1 uL 10uM P79
- 1 uL 10uM P80
- 1 uL plasmid (~1 ng/uL)

- 25 uL Phusion HS Flex 2xmm
- PCR tube

Total volume = 50 uL

PCR protocol (3 tubes used)

98 deg C for 30s

Repeat 10x:

98 deg C for 15s

53.4deg C for 30s

72 deg C for 2:00 mins

Repeat 25x:

98 deg C for 15s

64.7 deg C for 30s

72 deg C for 2:00 mins

72 deg C for 5 minutes

4 deg C for inf

Gibson Assembly:

DPNI Digest (1uL added) of backbone and insert at 37deg C for 1 hour

PCR wash: Backbone, Insert

- Add 250uL of PB and mix well with each PCR reaction
- Pipette everything into a column from the miniprep kit and spin for 1 minute
 - Dump out PB into miniprep waste
- Wash with CWC (500 uL) and spin for 1 minute
 - Dump out CWC into miniprep waste
- Spin dry column for 1 minute
 - Dump out any elution
- Add 30uL of nf water into sterile 2mL tube, spin for 1 minute to elute

Gibson concentrations					
	A	B	C	D	E
1		ng/uL	bp length	nmol/uL	nM
2	Backbone	63.8	6103	6.37123e-5	63.71231483
3	Insert	36.0	3159	1.9743e-5	19.7429951

Gibson Volumes			
	A	B	C
1		Assembly	Negative control
2	BB	1 uL	1 uL
3	Insert	1 uL	1 uL
4	nf water	8 uL	18 uL
5	Gibson MM	10 uL	---

Place Gibson mix into PCR machine, 15 minutes @ 50°C

Transformation using gibson protocol attached. Transformation successful.

Table3				
	A	B	C	
1	SS	Primers	Sequence	
2	NapA (1 a,b)	REV	TGCCACACCGGGGACCGACAGA	With pC
3		FWD	aagagaaattacattctgggattagatataggcattactagcgtgggtacgga	
4	YcbK (2 a,b)	REV	AGCGAAGGCAGGAGTAGGCAGGAT	With pC (pSB1C)
5		FWD	aagagaaattacattctgggattagatataggcattactagcgt	
6	YcdB (3 a,b)	REV	GGCGTGTGCCACGGGGCATGA	With pC
7		FWD	aagagaaattacattctgggattagatataggcattactagcgtgggtacggaata	

- Materials (Per RXN)

- o 1 uL 10 uM FW
- o 1 uL 10 uM REV
- o 1.5 uL DMSO
- o 1 uL plasmid DNA (1ng/uL)
- o 25 uL of Phusion HS Flex 2x MM
- o PCR Tube

Thermocycler conditions:

98 C for 30s

98 C for 15s (Repeat Red 35x)

64 C for 30s

72 C for 2.2 minutes

72 C for 10 min


4 C for inf. Time

Gel results

Table4

	Lane number	Description
1	1	NapA a
2	2	NapA b
3	3	YcbK a
4	4	YcbK b
5	5	YcdB a
6	6	YcdB b
7	7	pBAD A
8	8	pBAD B

Gel results



Digestion with DpnI

- Add DpnI (1uL) to the purified PCR product and incubate for 1-4 hours (37 degC)
 - Incubated for 60 minutes

PCR wash

- Add 250uL of PB and mix well with each PCR reaction
 - Used the one in the Qiagen kit
- Pipette everything into a column and spin for 1 minute
- Wash with CWC (500 uL) and spin for 1 minute
- Elute using 30uL of ddH2O

Nanodrop results:

Nanodrop				
	A	B	C	D
1		Concentration	260/280	260/230
2	NapA	88 ng/uL		
3	YcbK - pSB1C3	27 ng/uL		
4	YcdB	55 ng/uL		
5	YcbK - pBAD18cm	84ng/ul		

FRIDAY, 9/15/17

Ligation

Procedure:

1. Calculate the volume of the DpnI-treated PCR product containing 50-150 ng and fill in the table below.
2. Calculate the volume of water to make up to 10 μ L.
3. Mix the components below:
4. Incubate for 120 min at room temperature (~22°C).

Reaction Volumes for Ligat...					
	A	B	C	D	E
1	Reagent	Nap A	YcbK - pSB1C3	YcdB	Ycbk - pBad18cm
2	ddH2O	7 uL	6 uL	7 uL	7 uL
3	Purified PCR product (50 ng)	1 uL	2 uL	1 uL	1 uL
4	10x reaction buffer	1 uL	1 uL	1 uL	1 uL
5	T4 DNA ligase, 5U/uL	1 uL	1 uL	1 uL	1 uL
6	Total:	10 uL	10 uL	10 uL	10 uL

Perform Transformations:

Volumes for iPCR Transformations			
	A	B	C
1	30 mL C-cells	5 uL iPCR ligated product	35 uL rescue media

Transformations performed and cells incubated at 37 deg C.

SUNDAY, 9/17/17

Overnight cultures for: ClyA-saCas9-His6, NapA-saCas9-His6, YcdB-saCas9-His6, YcbK-saCas9-His6 (pBad and pSB1C3) started (2 of each). Cm antibiotic used.

For TEM experiment:

3 cultures (5 mL) of Top10 cells and 3 cultures of JC8031 were started

Transformation with Gibson Product

Introduction

Transformation into competent cells after Gibson assembly

Protocol given by Chelsea

Materials

- › 50 uL competent cells per Gibson assembly DNA
- › 15 uL of Gibson assembly DNA
- › 65 uL SOC media
- › Agar plates
- ›

Procedure

Transformation protocol

- ✓ 1. Place 50 uL of competent cells into a microcentrifuge tube
(Or use original vial, if it already has 50 uL)
Keep on ice
- ✓ 2. Pipette 15 uL of Gibson assembly DNA into each tube
- ✓ 3. Incubate on ice for 20 minutes
- ✓ 4. Heat shock tubes in water bath at 42°C for 60s
- ✓ 5. Incubate on ice for 5 minutes
- ✓ 6. Add 65 uL of SOC media, rescue for 1 hour & 15 minutes at 37°C in shaker
- ✓ 7. Plate 80 uL onto plate with correct antibiotics