AHK4 Assay Protocol

Materials

E. coli KMI002 strain was used in the experiments. The KMI002 (ΔrcsC, *cps::lacZ*) strain lacks *rcsC* gene and habors *cps::lacZ* fusion gene.

A high-copy plasmid, pSB1C3, was used in the experiments and the *ahk4* gene was inserted downstream the BAD/*araC* promoter, an L-arabinose inducible promoter, followed by a ribosome binding site

Qualitative experiment

- 1. LB agar plates containing chloramphenicol (34 μg/mL) were prepared.
- 2. 50 μ l of X-Gal (50 mg/ml), 10 μ l of 100 mM iP or DMSO as a control, and 40 μ l of LB medium was mixed in microtubes. Then the solutions were applied to the agar plates.
- 3. Samples were inoculated and incubated at room temperature.
- 4. Photographs were taken after sufficient blue color was developed.

Quantitative experiment

- 1. Overnight culture of samples were prepared in 2 ml of LB medium containing chloramphenicol (34 μg/mL) at 25°C.
- 2. Samples were diluted for 2000-fold in 1ml of fresh LB medium containing chloramphenicol (34 μ g/mL) and various concentration of IP (10 nM-100 μ M). Cells were also inoculated into medium containing DMSO instead of iP.
- 3. Samples were cultured overnight at 900 rpm at 25°C.
- 4. Cells were collected by centrifugation at $10,000 \times g$ for 10min.
- All of supernatant was discarded and then cells were resuspended in 500 μL of PBS buffer containing 1 mM MgSO₄ and 1 mM dithiothreitol (DTT). Also 500 μL of the same buffer in was prepared as a control for spontaneously splitting of ONPG.
- 20 μL of each suspension was added into 180μL of the buffer used above and Abs600 was measured and recorded by a microplate reader.
- 7. 10μL of 0.1% SDS and 10 μL of chloroform was added into each tube including the control and vortexed for 15sec.
- 8. Tubes were heated at 28°C for 5min.
- 9. 100 μL of ONPG (4 mg/mL) was added to each tube and incubated at 37°C for 30min. ONPG was dissolved in the buffer used above.
- 10. After 30min incubation, tubes were heated at 65°C for 10min to inactivate β -galactosidase.

- 11. All samples were centrifuged at 15,000 rpm for 10min.
- 12. Abs420 of supernatant was measured and recorded by a microplate reader. The control was used as a blank.
- 13. Relative β -galactosidase activity was calculated by following formula:

Relative
$$\beta$$
 -galactosidase activity =
$$\frac{Abs420}{Abs600 \cdot 10 \cdot 30min}$$