

Transformation – CaCl Method

(from Christopher)

- inoculate 50 μ l from overnight culture (XL1 blue or DH5 α) in 10 ml LB
(for 10 transformations)
- let grow to $OD_{600} = 0,4 - 0,6$ at 37°C with agitation (at least 200 rpm)
- prepare 0,2 M CaCl solution and **put it on ice!**
(everything should be on ice all the time)
- centrifuge the cells at fullspeed for 1 min
- resuspend cells with 1 ml of 0,2 M CaCl solution
- use 100 μ l cells for transformation
- add DNA (ca. 100 ng of plasmid)
- put it on ice for 30 min
- heat shock for 90 s at 42°C
- put it on ice again for 2 min
- add 400 μ l of fresh LB
- let grow for 30 – 60 min at 37°C with agitation (at least 200 rpm)
- centrifuge the cells for 1 min and resuspend each tube with 100 μ l LB
- plate on selective media