

## Protein purification protocol

### A. Note :

Before entering the laboratory: Put on facemasks and gloves to prevent enzymes dissolving the proteins.

### B. Theory :

1. The experiment focuses on Immobilized Metal Affinity Chromatography (IMAC)
  - There are some special sequences in many proteins and they can help proteins to specifically bind some metal ions.
  - However, the specific binding character needs a kind of molecule called "supports".
  - "Supports" can combine with groups that can covalently bond to metals, such as nitrilotriacetic acid, these groups would combine with metal ions, becoming affinity absorptions.
  - We usually modify our proteins and add a special sequence called "His-Tag", which consists of 6 histidines. It appears in the end of the recombinant protein.
  - After recombinant protein attaches onto the supports with nickel, they can be washed out with imidazole.
2. Histidine tag (His-tag)
  - Histidine has high affinity towards metal ions.
  - His-tag could be added into N or C terminal of a recombinant protein.
  - It can be expressed using different protein expression systems: bacteria, yeast, insect cells, and mammalian cells etc.
  - There are several advantages :
    - ✓ Small molecular size (4 to 10 amino acids sequence)
    - ✓ High purification efficiency
    - ✓ It would not change the function and activity of protein (in most cases).
    - ✓ Broad applications

### C. Materials and equipment :

- Column-purified
- Equivalent buffer, Wash buffer, Elution buffer
- Nickel Resin
- A few eppendorfs, 15 ml and 50 ml centrifuge tubes

### D. Laboratory procedures :

1. Add NaCl and Imidazole into the supernatant of bacterial cell lysate and make the final concentration to the chart below and shake it for 30 mins at 4°C

NaCl	150 mM
Imidazole	10 mM

2. Nickel-agarose preparation :
  - Because nickel resin would form layers of clear supernatant liquid and precipitation, so it would be better to mix well.
  - Add 2 ml resin into the column
  - Open the valve of the column and discard the effusive liquid.
  - Wash the column with equivalent buffer buffer and discard the effusive liquid
  - Close the valve
3. Batch binding :
  - Add the cell lysate mixture into the column and mix it with resin well
  - Transfer all the liquid into a tube and shake it for 30 mins at 4°C
4. Flow Through :
  - After shaken, move the mixture back into the column
  - Open the valve and collect the liquid flowing out of the column.
  - Label the liquid as "Flow Through (FT)" and put it on ice.

**Note** : All the samples should be put on ice from this step.
5. Collection-Wash buffer :
  - Wash the column with wash buffer
  - Use Eppendorfs to collect the liquid(1 ml for each Eppendorf)
  - Label the first 3 Eppendorfs as "W<sub>1</sub>", "W<sub>2</sub>", "W<sub>3</sub>"
  - Label the last Eppendorf as "WF"
6. Collection-Elution buffer :
  - Elute the column with elution buffer
  - Use Eppendorfs to collect the liquid(1 ml for each Eppendorf)
  - Collect around 10 ml and label them as "E<sub>1</sub>", "E<sub>2</sub>", ... "E<sub>10</sub>"
7. Store all the sample at -20°C
8. Wash resin :
  - Wash the column with solution I (concentration : 500 mM NaCl, 500 mM Imidazole)
  - Wash the column with ddH<sub>2</sub>O
  - Wash the column with 30% alcohol
  - Close the valve
  - Mix the resin with 30% alcohol and store at 4°C.