

Legend:

aiiA: Autoinducer inactivation enzyme from Bacillus; hydrolyzes acetyl homoserine lactone (BBa_C0060)

yhjH: yhjH Gene From E.coli str. K12 (BBa_K861090)

epsE: EpsE Molecular Clutch Gene of B. subtilis (BBa_K143032)

P+RBS: Promoter T7 and RBS (BBa_K525998)

RBS: RBS (Elowitz 1999) -- defines RBS efficiency (BBa_B0034)

T: T1 from E. coli rrnB (BBa_B0010)

DT: double terminator (B0010-B0012)

May

May/29/2017

Protocol: Sterilization and antibiotic stock

Lab materials were sterilized at 121°C for 15 minutes.

- 10 µL tips box (3)
- 200 µL tips box (3)
- 1 mL tips box (3)
- 1.5 mL Eppendorf tubes bottle (2)
- L shaped glass spreader (2)
- Pipettes (5)

First, calculations were done using as a reference the webpage openwetware.com. Each antibiotic was weighed and added to a 1.5 mL amber eppendorf tube, then, 1.5 mL of absolute alcohol was added and mixed until homogeneous. Stored at -20°C.

Chloramphenicol: 0.051g, labeled as "Cl".

Ampicillin: 0.075g, labeled as "Am".

May/30/2017

Protocol: Antibiotic stock and competent cells

First, calculations were done in order to get a concentration of 100 µg/mL. 15g of ampicillin was weighed and added to a 1.5 mL amber eppendorf tube, then, 1.5 mL of absolute alcohol was added and mixed until homogeneous. Stored at -20°C.

TOP10 competent cells

May/31/2017

Protocol: Parts resuspension

From the distribution kit, transformed *E. coli* competent cells and plated on CAM and AMP plates.

DNA	Plate	Well	Antibiotic
P+RBS	1 (2015)	3M	CAM
aiiA	2	4H	CAM
T1	4 (2015)	16G	AMP
DT	3	3F	CAM
RFP	4	4B	CAM

June

June/01/2017

Protocol: 100 mM and 50 mM CaCl₂ preparation and overnight cultures

CaCl₂ (5.59 g) added to 500 mL and 1 L of distilled water.

LB medium (10 mL) inoculated with TOP 10 (200 µL). 37°C, no agitation.

June/02/2017

Protocol: Competent cells

TOP 10 competent cells.

June/03/2017

Protocol: Transformation

Transformation P+RBS (CAM), aiiA (CAM), RFP (CAM), T (AMP) and RFP (AMP) into TOP 10 competent cells.

June/04/2017

Results: Inoculated plates after 22 incubation hours showed:

Negative Control (NC) - LB (w/o antibiotic) → Countless extended growth.

NC - LB (CAM) → Extended growth with 8 punctuated colonies, translucent and yellowish.

P+RBS (CAM) → Two colonies, medium and big size, both translucent and yellowish; entire edges and slightly elevated.

aiiA (CAM) → Two colonies, one punctuated colony and the other one small and round, both translucent and yellowish; entire edges and slightly elevated.

RFP (CAM) → No growth.

T (AMP) → One whitish, opaque, slightly dry, big size colony.

RFP (AMP) → Same morphology as T (AMP), one whitish, opaque, slightly dry, big size colony.

Protocol: Overnight cultures

8 LB medium tubes (10 mL) inoculated by duplicate with: P+RBS, aiiA, T and RFP. It was taken only half colony from those plates with just one colony, also, an extra tube was inoculated with a colony from the NC - LB (CAM) in order to see why this occurred.

June/05/2017

Results: From the overnight cultures, both aiiA and T, showed turbidity. There was no growth at all in NC - LB (CAM) tube which leads to conclude that contamination occurred during plate extension.


Protocol: Miniprep, Nanodrop and Electrophoresis

Plasmid extraction of: aiiA, P+RBS and T.

Nanodrop

Part	First Measurement	Second Measurement	Third Measurement	Average
P+RBS	36.3/1.98	36.6/1.94	35.9/1.93	36.26/1.95
aiiA	116.7/2.03	119/2.02	118.7/2.02	118.13/2.02
T	35.3/1.82	34.1/1.88	33.9/1.89	34.43/1.86

Measurements: Concentration / 260/280



June/06/2017

Protocol: Electrophoresis, psB1A3 resuspension, Nanodrop and competent cells



Plasmid resuspension

DNA	Plate	Well	Antibiotic
psB1A3	4	2H	AMP

Nanodrop

Part	First Measurement	Second Measurement	Average
psB1A3	2.6/1.63	1.5/1.74	2.05/1.68

Measurements: Concentration / 260/280

TOP 10 competent cells.

June/07/2017

Protocol: DNA resuspension and transformation

Parts taken from the 2017 kit:

DNA	Plate	Well
aiiA	2	4H
yhjH	1	13D
P+RBS	1	3M
GFP	3	6I
RFP	3	6G
CFP	3	19C

Transformation of yhjH, aiiA, P+RBS, RFP, CFP and GFP into TOP 10 competent cells.

June/08/2017

Results: From RFP, CFP and GFP plates, only RFP and CFP showed fluorescence.

Protocol: Overnight cultures

LB (CAM) plates inoculated with aiiA, yhjH, P+RBS, RFP, GFP and CFP due to some irregularities seen on the first plates.

June/09/2017

Results: RFP and CFP showed fluorescence. LB medium (CAM) tubes were inoculated with RFP and CFP (1 μ L/mL), 37°C for 24 hours.

Protocol: Overnight cultures

LB medium inoculated with DH5 α , 37°C and 225 rpm agitation.

June/10/2017

Protocol: Antibiotic stock and transformation

1.5 mL of CAM prepared. Transformation of yhjH, aiiA, P+RBS and RFP into DH5 α competent cells.

June/11/2017

Protocol: Overnight cultures

10 mL LB medium tube inoculated with TOP 10 cells. Incubated for 18 hours, 37°C and 225 rpm agitation.

June/12/2017

Protocol: Alkaline Lysis Miniprep solutions preparation and overnight cultures

Solution 1 (100 mL): Sucrose (15 g), Tris (0.6057 g), EDTA (1.8612 g) and distilled water (80 mL). pH adjusted to 7.9 and sterilized for 15 min.

Solution 2 (100 mL): NaOH (8 g), SDS (1 g) and distilled water (80 mL).

Solution 3 (100 mL): Potassium acetate (29.442 g), Glacial acetic acid (11.5 mL) and distilled water (88.5 mL).

40 mL LB medium tube inoculated with TOP 10 cells. 37°C and 225 rpm agitation.

June/13/2017

Protocol: Competent cells

TOP 10 competent cells.

June/14/2017

Protocol: Transformation

Transformation of yhjH, aiiA, P+RBS, RFP and NC into TOP 10 competent cells.

June/15/2017

Results: RFP, yhjH and aiiA showed growth of a single colony on each one of them, RFP's colony with a reddish color.

Protocol: Overnight cultures and sodium acetate solution preparation

3 LB medium tubes (10 mL) inoculated with 300 µL of: RFP (2) and CFP. 37°C and 225 rpm agitation.

Solution 3M (50 mL): Sodium acetate (12.3045 g) and distilled water (30 mL).

June/16/2017

Protocol: Miniprep

Alkaline lysis plasmid extraction by duplicate of RFP1 and RFP2.


June/17/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Third Measurement	Average
RFP1	2592.3/2.01	2594.8/1.97	2454.5/1.97	2547.2/2
RFP2	3750.5/2.07	3746.1/2.07	3572/2.08	3689.5/2.07

Measurements: Concentration / 260/280



Results: The weight expected was of 3139 pb, the sixth strip has a weight around 3000 pb, slightly below than expected but contaminated by proteins.

June/18/2017

Protocol: Overnight cultures

6 LB+CAM medium (60 mL) tubes inoculated with 300 μ L of aiiA (2), yhjH (2), RFP and NC. 37°C and 225 rpm agitation.

June/19/2017

Protocol: Competent cells and transformation

TOP 10

Transformation of P+RBS, epsE, RFP1, RFP2, RFP3 and NC into TOP 10 competent cells.

June/20/2017

Results: From the transformation, RFP1, RFP2, RFP3 and P+RBS showed growth, it was not the case for epsE because it was cultured in CAM instead of AMP.

Protocol: Miniprep and Nanodrop

Plasmid extraction (alkaline lysis) by duplicated of aiiA, yhjH and RFP.

Nanodrop

Part	First Measurement	Second Measurement	Average
yhjH1	5362/2.05/2.36	5383/2.05/2.35	5372.5/2.05/2.355
aiiA1	2932.8/1.98/2	3026.7/1.98/2	2979.75/1.98/2
RFP1	992.1/1.95/2.21	1156.5/1.91/1.92	1074.3/1.93/2.06
yhjH2	288.4/1.83/1.81	286/1.84/1.82	287.2/1.835/1.815
aiiA2	1903.8/2.03/2.28	1936.6/2.01/2.29	1920.2/2.02/2.285

RFP2	179/1.93/2.27	185.5/1.92/2.28	182.25/1.925/2.275
------	---------------	-----------------	--------------------

Measurements: Concentration / 260/280 / 260/230

June/21/2017

Protocol: Miniprep, overnight cultures and electrophoresis

Plasmid extraction (alkaline lysis) by duplicated of P+RBS.

2 plates inoculated with 200 μ L of epsE and NC each one, and one last plate as a NC. 37°C and 225 rpm agitation.



June/22/2017

Protocol: Antibiotic stock, Miniprep, Nanodrop and electrophoresis


1.5 mL of AMP prepared.

Plasmid extraction (alkaline lysis) by duplicated of P+RBS.

Nanodrop

Part	First Measurement	Second Measurement	Third Measurement	Average
P+RBS1	38.3/1.5/1.04	47.3/1.47/0.99	50.9/1.52/1.03	45.5/1.496/1.02
P+RBS2	1385.4/1.68/1.2	733.5/1.63/1.28	303.4/1.7/1.26	807.4/1.67/1.24

Measurements: Concentration / 260/280 / 260/230



June/23/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Third Measurement	Average
P+RBS1	240.6/1.98/2.24	244.5/1.96/2.23	251.8/1.94/2.18	245.6/1.96/2.21
P+RBS2	30.3/1.69/1.89	16.2/1.59/1.78	15.4/1.7/1.92	19.76/1.78/1.86

Measurements: Concentration / 260/280 / 260/230

June/24/2017

Protocol: Overnight cultures

Erwinia Amylovora cultured on a nutrient agar plate by the quadrant streak plate method, 28°C. epsE and NC cultured on LB agar plates added with AMP, 37°C.

June/25/2017

Protocol: Digestion and Ligation

aiiA digestion (Upstream Mix: 2.0358 µL of P+RBS, 1 µL of EcoRI, 1 µL of SpeI, 5 µL of 10x NE Buffer and 40.9642 µL of nuclease free water. Downstream Mix: 2.0833 µL of aiiA, 1 µL of XbaI, 1 µL of PstI, 5 µL of 10x NE Buffer and 40.9167 µL of nuclease free water. Vector Mix: 30.39 µL of psB1A3, 1 µL of EcoRI, 1 µL of PstI, 5 µL of 10x NE Buffer and 12.61 µL of nuclease free water).

June/26/2017

Protocol: TAE 1x preparation, overnight cultures ,Nanodrop and electrophoresis

TAE 1x (500 mL): 10 mL of TAE 50x and 490 mL of distilled water, sterilized (121.5 °C) for 15 min.

Cultures of P+RBS, aiiA, yhjH and *Erwinia Amylovora*.

Nanodrop

Part	First Measurement	Second Measurement	Average
Downstream Mix	194.3/1.88/1.22	220/1.88/1.23	207.15/1.88/1.225
Upstream Mix	107.3/1.42/0.49	103.9/1.44/0.52	105.6/1.43/0.505
Vector Mix	38.9/1.03/0.26	74.3/1.45/0.66	56.6/1.24/0.46
Ligation	728.1/3.95/2.12	702.9/4.31/2.6	71505/4.13/2.36

Measurements: Concentration / 260/280 / 260/230



June/27/2017

Protocol: Miniprep, Buffer TE preparation and overnight cultures

Plasmid extraction of: P+RBS, aiiA and yhjH.

Plasmid extraction by alkaline lysis of: P+RBS, aiiA and yhjH.

Buffer TE: 0.0605 g of Tris and 0.0186 g of EDTA, pH 7.9 and 40 mL of distilled water.

3 LB medium (15 mL) tubes inoculated with 200 µL of aiiA, P+RBS and yhjH. 37°C and 225 rpm agitation.


June/28/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Average
P+RBS	5.5/1.62/4.67	6.9/1.69/3.17	6.2/1.655/3.92
yhjH	18.5/1.84/2.86	13.2/1.67/3	15.85/1.755/2.43
aiiA	10.4/1.62/6.05	7.5/1.6/3.74	8.95/1.61/4.895

Measurements: Concentration / 260/280 / 260/230



June/29/2017

Protocol: Digestion, Ligation and electrophoresis

aiiA digestion (Upstream Mix: 10.266 μ L of P+RBS, 1.5 μ L of EcoRI, 1.5 μ L of SpeI, 5 μ L of Cutsmart and 31.734 μ L of nuclease free water. Downstream Mix: 6.311 μ L of aiiA, 1.5 μ L of XbaI, 1.5 μ L of PstI, 5 μ L of Buffer 3.1 and 35.686 μ L of nuclease free water. Vector Mix: 30.39 μ L of psB1A3, 1.5 μ L of EcoRI, 1.5 μ L of PstI, 5 μ L of Buffer 3.1 and 11.6048 μ L of nuclease free water.)

Ligation Mix: 2 μ L of Upstream Mix, 2 μ L of Downstream Mix, 2 μ L of Vector Mix, 2 μ L of Ligase Buffer, 1.5 μ L of Ligase and 10.5 μ L of nuclease free water.



June/30/2017

Protocol: Digestion, Ligation and overnight cultures

aiiA digestion (Upstream Mix: 10.266 μ L of P+RBS, 1.5 μ L of EcoRI, 1.5 μ L of SpeI, 5 μ L of Cutsmart and 31.734 μ L of nuclease free water. Downstream Mix: 6.311 μ L of aiiA, 1.5 μ L of XbaI, 1.5 μ L of PstI, 5 μ L of Buffer 3.1 and 35.7 μ L of nuclease free water. Vector Mix: 30.39 μ L of psB1A3, 1.5 μ L of EcoRI, 1.5 μ L of PstI, 5 μ L of Buffer 3.1 and 11.6048 μ L of nuclease free water.)

Ligation Mix: 2 μ L of Upstream Mix, 2 μ L of Downstream Mix, 2 μ L of Vector Mix, 2 μ L of Ligase Buffer, 1.5 μ L of Ligase enzyme and 10.5 μ L of nuclease free water.

Nutrient broth (10 mL) tube inoculated with TOP 10, 37°C.

July

July/01/2017

Protocol: Competent cells and transformation

TOP10 competent cells

Transformation of epsE, DT and RFP into TOP 10 competent cells.

July/02/2017

Protocol: Overnight cultures

Two nutrient broth tubes (50 mL) inoculated with TOP 10 cells.

July/03/2017

Protocol: Competent cells and transformation

TOP10 competent cells

Transformation of epsE, DT, RFP and NC into TOP 10 competent cells.

July/04/2017

Protocol: Overnight cultures

P+RBS, yjhH and aiiA cultured in nutrient broth 2x, 37°C and 225 rpm agitation.

July/05/2017

Protocol: Digestion, Ligation and competent cells

This protocol was done twice, a four hours aiiA digestion and twelve hours aiiA digestion (Upstream Mix: 10.266 µL of P+RBS, 1.5 µL of EcoRI, 1.5 µL of SpeI, 5 µL of Cutsmart and 31.734 µL of nuclease free water. Downstream Mix: 6.311 µL of aiiA, 1.5 µL of XbaI, 1.5 µL of PstI, 5 µL of Buffer 3.1 and 35.686 µL of nuclease free water. Vector Mix: 30.39 µL of psB1A3, 1.5 µL of EcoRI, 1.5 µL of PstI, 5 µL of Buffer 3.1 and 11.6048 µL of nuclease free water.)

Ligation Mix: 2 µL of Upstream Mix, 2 µL of Downstream Mix, 2 µL of Vector Mix, 2 µL of Ligase Buffer, 1.5 µL of Ligase enzyme and 10.5 µL of nuclease free water.

TOP 10 cells.


July/06/2017

Protocol: Transformation

Transformation of epsE, DT, RFP and NC into TOP 10 competent cells.

July/07/2017

Protocol: Electrophoresis



July/10/2017

Protocol: Competent cells

TOP 10 competent cells.

July/11/2017

Protocol: Competent cells, transformation

TOP 10 competent cells.

Transformation of epsE, DT and RFP into TOP 10 competent cells.

July/12/2017

Protocol: Overnight cultures

Four nutrient broth tubes (10 mL) inoculated with aiiA, yhjH and P+RBS, 37°C and 225 rpm agitation.

July/13/2017

Protocol: Overnight cultures

Four LB broth tubes (10 mL) inoculated with epsE, DT and NC, 37°C and 225 rpm agitation.

July/14/2017

Protocol: Miniprep

Plasmid extraction by duplicate (alkaline lysis) of DT, epsE and NC.

July/15/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Average
epsE1	1897.5/1.77/1.49	1234.5/1.97/1.93	1566/1.87/1.71
epsE2	970.9/2.02/2.16	1024.3/1.97/1.93	997.7/1.995/2.045
DT1	1160.4/1.96/1.83	1164.7/1.97/1.85	1162.55/1.965/1.84
DT2	836.6/1.98/1.91	836.3/1.98/1.91	836.45/1.98/1.91
NC1	771.2/1.98/2.03	805.1/1.96/1.93	788.15/1.97/1.98
NC2	539.1/1.98/2.16	571.3/1.99/2.21	555.2/1.985/2.185

Measurements: Concentration / 260/280 / 260/230

July/16/2017

Protocol: Miniprep

Alkaline lysis plasmid extraction by duplicate of DT, epsE and NC.

July/17/2017
 Protocol: Nanodrop
 Nanodrop

Part	First Measurement	Second Measurement	Average
epsE1	977.9/1.87/1.2	948.1/1.84/1.19	963/1.855/1.195
epsE2	7773.8/1.93/1.49	6721.3/1.92/1.47	7247.55/1.925/1.48
DT1	1704.3/1.97/1.82	1704.6/1.97/1.82	1704.45/1.97/1.82
DT2	5590.7/1.99/1.9	5548.9/1.99/1.9	5569.8/1.99/1.9
NC1	1991/2.01/1.89	1981.3/2.01/1.89	1986.15/2.01/1.89
NC2	6203.6/2.01/1.87	6150.8/2.01/1.88	6177.2/2.01/1.875


Measurements: Concentration / 260/280 / 260/230

July/18/2017
 Protocol: Miniprep
 Alkaline lysis plasmid extraction by duplicate of DT, epsE and NC.

July/19/2017
 Protocol: Overnight cultures, Nanodrop and electrophoresis
 Nutrient agar added with sucrose inoculated with *Erwinia Amylovora*. Two nutrient broth tubes (10 mL) inoculated with TOP10 cells.
 Nanodrop

Part	First Measurement	Second Measurement	Average
epsE1	134.6/1.93/2.06	137/1.92/2.05	135.8/1.925/2.055
epsE2	89.4/1.76/1.02	99.7/1.73/0.95	94.55/1.745/0.985
DT1	594.2/1.95/2.16	581.3/1.95/2.17	587.75/1.95/2.165
DT2	1082.3/1.98/2	1080.5/1.97/2.09	1081.4/1.975/2.045
NC1	297.6/1.95/2.14	346.7/1.97/2.17	322.15/1.96/2.155
NC2	1172.8/1.99/2.17	1194.5/1.98/2.21	1183.65/1.985/2.19

Measurements: Concentration / 260/280 / 260/230



July/20/2017

Protocol: Competent cells, transformation

TOP10 competent cells

Transformation of NC, epsE, DT, Ligation 4h and Ligation 12h into TOP10 competent cells.

July/21/2017

Protocol: Overnight cultures

LB and nutrient broth added with sucrose, inoculated with DT, epsE, Ligation 4h and Ligation 12h.

July/22/2017

Protocol: Miniprep

Alkaline lysis plasmid extraction by duplicate of epsE, Ligation 4h and Ligation 12h.

July/23/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Average
epsE1	1591.4/2.05/2.03	1569.2/2.06/2.07	1580.3/2.055/2.05
epsE2	1285.9/2.09/2.36	1278.2/2.1/2.37	1282.05/2.095/2.365
L4h1	1780.5/2.1/2.34	1736.9/2.11/2.34	1758.7/2.105/2.34

L4h2	116.5/2.02/2.11	143.1/2.02/2.18	129.8/2.02/2.145
L12h1	2449.3/2.12/2.29	2325.2/2.11/2.29	2387.25/2.115/2.29
L12h2	1603.7/2.1/2.29	1604.7/2.08/2.28	1604.2/2.09/2.285

Measurements: Concentration / 260/280 / 260/230

July/24/2017

Protocol: Miniprep

Plasmid extraction by alkaline lysis of epsE, Ligation 4h and Ligation 12h.


July/25/2017

Protocol: Nanodrop and electrophoresis

Nanodrop

Part	First Measurement	Second Measurement	Average
epsE1	3355.5/2.12/2.28	3307.2/2.11/2.3	3331.35/2.115/2.29
epsE2	1483.9/2.03/1.91	1504.5/2.03/1.93	1494.2/2.03/1.92
L4h1	2655.7/2.11/2.23	2615.1/2.11/2.25	2635.4/2.11/2.24
L4h2	1359.9/2.1/2.36	1391.6/2.11/2.38	1375.75/2.105/2.37
L12h1	3515.6/2.14/2.25	3585.2/2.14/2.24	3550.4/2.14/2.245
L12h2	3857.3/2.14/2.23	3839.9/2.15/2.24	3845.6/2.145/2.235

Measurements: Concentration / 260/280 / 260/230



July/27/2017

Protocol: Miniprep

Plasmid extraction by alkaline lysis of epsE, Ligation 4h (colony 1,2,3) and Ligation 12h (colony 1,2,3).

July/28/2017

Protocol: Nanodrop, parts resuspension, competent cells, transformation and digestion
Nanodrop

Part	First Measurement	Second Measurement	Average
L4hC1	604.2/2.06/2.29	591.8/2.06/2.32	598/2.06/2.305
L4hC2	1973.9/2.13/2.23	1973.5/2.13/2.23	1973.7/2.13/2.23
L4hC3	2482.9/2.13/2.23	2496.8/2.13/2.23	2489.85/2.13/2.23
L12hC1	744.1/2.08/2.47	710.5/2.08/2.46	727.3/2.08/2.465
L12hC2	780.3/2.07/2.32	769/2.1/2.29	774.65/2.085/2.305
L12hC3	994.5/2.09/2.41	767.5/2.1/2.42	881/2.095/2.415
epsE	1880.7/2.12/2.21	1916.8/2.13/2.22	1898.75/2.125/2.214

Measurements: Concentration / 260/280 / 260/230

Parts resuspension

Part	Plate	Well
(T) BBa_B0011	4 (2016)	1B
(DT) BBa_B0015	3	3F

TOP10 competent cells

Transformation of T and DT into TOP10 competent cells.

yhjH digestion (Digestion Mix: 6.2 µL of yhjH, 1.5 µL of PstI, 1.5 µL of XbaI, 5 µL of 3.1 NE buffer and 35.8 µL of nuclease free water).

July/29/2017

Protocol: Miniprep, Nanodrop, ligation and electrophoresis


Plasmid extraction of P+RBS+ aiiA.

Nanodrop

Part	First Measurement	Second Measurement	Third Measurement	Average
P+RBS+aiiA	24.1/1.84/2.21	24.4/1.74/1.77	24/1.83/1.21	24.16/1.8/1.73

Measurements: Concentration / 260/280 / 260/230

yhjH ligation (Ligation Mix: 10.5 μ L of nuclease free water, 2 μ L of T4 ligase buffer, 2 μ L of digested Vector Mix, 2 μ L of Upstream Mix (P+RBS), 2 μ L of Downstream Mix and 1.5 μ L of T4 ligase).




August

August/01/2017

Protocol: Miniprep

Plasmid extraction of DT.

Plasmid extraction by alkaline lysis of yhjH+RBS.



August/02/2017

Protocol: Partial digestion, Miniprep, digestion, Nanodrop and electrophoresis

Partial digestion of T (Mix: 5 μ L of DNA, 2.5 μ L Buffer cutsmart, 0.5 μ L XbaI enzyme and 17 μ L of nuclease free water) and double terminator (Mix: 10 μ L of DNA, 2.5 μ L of buffer cutsmart, 0.5 μ L of XbaI enzyme and 12 μ L of nuclease free water).


Plasmid extraction of T.

Nanodrop

Part	First Measurement	Second Measurement	Average
yhjH+RBS1	3822.2/2.14	3855.3/2.14	3838.75/2.14
yhjH+RBS2	1124,3/2.25	1598/2.11	1361.15/2.18
yhjH+RBS3	3776.9/2.12	3718.4/2.01	3747.65/2.065

Measurements: Concentration / 260 / 280

12-hour digestion of T (Mix: 1.5 μ L XbaI enzyme, 1.5 μ L PstI enzyme, 5 μ L buffer, 33.3 μ L nuclease free water and 8.7 μ L of DNA) and CAM (PSB1C3) vector (Mix: 1.5 μ L EcoRI enzyme, 1.5 μ L PstI enzyme, 5 μ L buffer, 21.7 μ L nuclease free water and 20.3 μ L of the vector)



August/03/2017

Protocol: Miniprep, partial digestion, Nanodrop and Ligation

Plasmid extraction of P+RBS+yhjH


Partial digestion P+RBS+yhjH (Mix: 5 μ L of DNA, 2.5 μ L Buffer cutsmart, 0.5 μ L XbaI enzyme and 17 μ L of nuclease free water) and RBS (Mix: 5 μ L of DNA, 2.5 μ L Buffer cutsmart, 0.5 μ L XbaI enzyme and 17 μ L of nuclease free water), by an hour.

Nanodrop

Part	First Measurement	Second Measurement	Average
RBSC1	3066.1/2.1	3071.7/2.09	3068.5/2.95
RBSC2	86.1/1.91	86.6/1.94	86.35/1.925

Measurements: Concentration / 260/280

12-hour ligation for biobrick of aiiA (Mix: 2 μ L (P+RBS+aiiA) + 2 μ L T + 2 μ L CAM vector + 2 μ L Ligation buffer + 2 μ L ligase enzyme and 10.5 μ L nuclease free water).



August/04/2017


Protocol: Miniprep, partial digestion, competent cells and transformation

Plasmid extraction of RBS

Partial digestion of RBS (Mix: 5 10 μ L of RBS, 2.5 μ L Buffer cutsmart. 0.5 μ L XbaI enzyme and 12 μ L of nuclease free water),

TOP10 competent cells

Transformation of aiiA final biobrick into TOP10 competent cells.



August/08/2017

Protocol: Ligation

12-hour ligation of the yhjH final biobrick (Mix:2 μ L of P+RBS+yhjH + 2 μ L T+ 2 μ L CAM backbone + 2 μ L Ligation buffer + ligase enzyme + 10.5 μ L nuclease free water).

August/09/2017

Protocol: Antibiotic stock, chemically competent cells, transformation and miniprep

TetR stock with a concentration of 35mg/mL in a volume of 1.5mL

Plasmid extraction of aiiA biobrick (6 colonies).

BL21(DE3) Chemically competent cells.

Transformation of P+RBS+epsE (TetR), biobrick of aiiA (CAM) and yhjH biobrick (CAM).

August/10/2017

Protocol: Overnight cultures, Nanodrop, partial digestion, miniprep and electrophoresis
5 LB broth tubes inoculated with 5 colonies of yhjH biobrick, 37 °C, 225 r.p.m.

Nanodrop


Part	First Measurement	Second Measurement	Average
Biobrick yhjH C1	2559.8/2.09	2559.3/2.08	2559.55/2.085
Biobrick yhjH C2	1611.2/2.08	1336/2.08	1473.6/2.08
Biobrick yhjH C3	26.8/1.91	27.4/1.88	27.1/1.895
Biobrick yhjH C4	4284/2.07	4326.9/2.09	4305.45/2.08
Biobrick yhjH C5	1681.9/2.13	1775.1/ 2.14	1728.5/2.135

Measurements: Concentration / 260/280

Plasmid extraction by alkaline lysis of yhjH+RBS

Partial digestion of the aiiA biobrick (5 colonies), all colonies followed this mix: (5 μ L of DNA, 2.5 μ L Buffer cutsmart. 0.5 μ L XbaI enzyme and 17 μ L of nuclease free water).

Partial digestion of aiiA biobrick at the standard conditions (5 colonies).



August/11/2017

Protocol: Digestion, Ligation, overnight cultures and miniprep

12-hour digestion of yhjH+RBS (Mix: 5.74 μ L of (P+RBS+yhjH) + 5 μ L Cutsmart buffer + 1.5 μ L SpeI + 1.5 μ L EcoRI and 10.5 μ L nuclease free water) and 12-hour digestion of RBS (Mix: 5.74 μ L of RBS + 5 μ L Cutsmart buffer + 1.5 μ L SpeI + 1.5 μ L EcoRI and 10.5 μ L nuclease free water). 12-hour ligation of epsE+T (Mix: 2 μ L of epsE + 2 μ L T+ 2 μ L CAM backbone + 2 μ L Ligation buffer T4 + 10.5 μ L nuclease free water + T4 Ligase enzyme).

LB broth tubes inoculated with aiiA biobrick.

Plasmid extraction by alkaline lysis of 5 colonies of the yhjH biobrick.

August/12/2017

Protocol: Miniprep, 12-hour digestion, partial digestion, Nanodrop and electrophoresis

Plasmid extraction by alkaline lysis of yhjH.

Plasmid extractions of yhjH biobrick and aiiA biobrick (2 colonies).



12-hour digestion of kanR backbone (7.2 μ L vector, 5 μ L of buffer 3.1, 1.5 μ L PstI, 1.5 μ L EcoRI, 35 μ L of nuclease free water).

Partial digestion of yhjH (5 colonies), aiiA biobrick (2 colonies).

Nanodrop

Part	First Measurement	Second Measurement	Average
yhjH C1	781.2/2.08	699.2/1.94	740.2/2.01
yhjH C2	1807.1/ 2.09	1773.8/2.09	1790.45/2.09
yhjH C3	1472.1 / 2.10	1339/2.08	1405.55/2.09
yhjH C4	1946.8/2.12	1420.4/2.08	1683.6/2.1
kanR Backbone	10.7/2.3	11.9/ 2.27	11.3/2.285

Measurements: Concentration / 260/280



August/13/2017

Protocol: Ligation and overnight cultures

12-hour ligation of yhjH in kanR vector.

LB broth tubes inoculated with 4 colonies of the yhjH biobrick.

August/14/2017

Protocol: Miniprep, competent cells, transformation and overnight cultures

Plasmid extraction of aiiA biobrick.

TOP10 competent cells

Transformation of RFP, RBS + yhjH, epsE+T and RBS+epsE into TOP10 competent cells.

4 LB agar plates inoculated with epsE+T (CAM), RFP (CAM), RBS+yhjH (TetR) and RBS+epsE (TetR).

August/15/2017

Protocol: Overnight cultures, partial digestion, miniprep, chemically competent cells and transformation

LB agar plate supplemented with sucrose inoculated with Erwinia Amylovora.

Partial digestion by alkaline lysis of yhjH biobrick (5 colonies).

Plasmid extraction of aiiA biobrick.

TOP10 competent cells

Transformations of epsE+T (CAM), RBS+yhjH (TetR), RBS+epsE (TetR) and RFP (TetR).

August/16/2017

Protocol: Overnight cultures and electrophoresis

LB agar plate (CAM) inoculated with RBS + epsE.



August/17/2017

Protocol: Overnight cultures and calibration curve of *Erwinia Amylovora*

12 LB agar plates inoculated with yhjH (CAM), epsE+T (CAM), RBS+epsE (5 colonies) (CAM) and RBS + yhjH (2 colonies) (TetR).

Calibration curve of *Erwinia Amylovora* was prepared to predict the growth rate, to know the exact time when *Erwinia* reaches its optimal density to begin with the competent cells protocol.

August/18/2017

Protocol: Miniprep

Plasmid extraction of yhjH biobrick.

Plasmid extraction by alkaline lysis of epsE+T, yhjH + RBS (2 colonies) and epsE+RBS (5 colonies).

August/19/2017

Protocol: Transformation

Transformation of *Erwinia Amylovora* with aiiA into BL21(DE3) electrocompetent cells.

August/20/2017

Protocol: Chemically competent cells and overnight cultures

BL21(DE3) Chemically competent cells.

LB agar (CAM) inoculated with P+RBS+epsE, LB agar (TetR) inoculated with RBS+epsE, RBS+yhjH, epsE+T and RFP.

August/25/2017

Protocol: Competent cells and transformation

TOP10 competent cells.

Transformation of yhjH biobrick and RBS+yhjH.

August/26/2017

Protocol: Overnight cultures, competent cells and transformation

5 LB broths supplemented with sucrose tubes (CAM) inoculated with *Erwinia Amylovora* previously transformed with biobrick of aiiA. 3 LB broth tubes (CAM) inoculated with yhjH biobrick. LB broth tube (TetR) inoculated with yhjH+RBS, 28 °C, 255 rpm.

TOP10 competent cells.

Transformations of RBS+yhjH and P+RBS+epsE (TetR).


August/27/2017

Protocol: Overnight cultures, miniprep and electrophoresis

LB agar supplemented with sucrose inoculated with *Erwinia Amylovora*.

LB broth tubes inoculated with yhjH biobrick and RBS+yhjH.

Plasmid extraction by alkaline lysis of yhjH, aiiA and RBS+epsE.



August/28/2017

Protocol: Miniprep, partial digestion and electrophoresis

Plasmid extraction by alkaline lysis of biobrick yhjH (3 extractions) and RBS+yhjH.

Partial digestion of the extraction at the standard conditions.



August/29/2017

Protocol: Partial digestion

Partial digestion of yhjH biobrick (3 colonies) and RBS + yhjH.

August/30/2017

Protocol: Partial digestion, miniprep and Nanodrop

Partial digestion of the plasmid extractions of yhjH biobrick and RBS+ yhjH both of them with the standard mix of partial digestion at the standard conditions.

Plasmid extraction of yhjH biobrick and RBS+epsE.

Nanodrop

Part	Concentration	260/280	260/230
RBS+epsE_1	36.6	1.92	6.25
RBS+epsE_2	38.5	1.93	5.09
yhjHBL21_1	85.3	1.87	2.63
yhjHBL21_2	84.1	1.82	2.33



August/31/2017

Protocol: Ligation and electrophoresis

12-hour ligation of RBS+yhjH) in TetR vector and (P+RBS+epsE) in TetR vector.



September

September/02/2017

Protocol: Miniprep

Plasmid extraction of biobrick aiiA and biobrick of yhjH.

September/03/2017

Protocol: Miniprep

Plasmid extraction by alkaline lysis of RBS+yhjH (5 colonies) and P+RBS+epsE (5 colonies).

September/05/2017

Protocol: Overnight cultures

8 LB broth tubes (5 mL) 4 tubes added with CAM and the other 4 tubes added with TetR, inoculated with *Erwinia Amylovora*, RBS+yhjH, P+RBS+epsE (5 colonies) and RBS+aiiA.

September/06/2017

Protocol: Miniprep and characterization

Plasmid extraction by alkaline lysis of *Erwinia* wild and transformed aiiA *Erwinia* (4 colonies).
Trials of characterization: Inoculation of *Erwinia Amylovora* in apple slices, in a total of 6 plates, *Erwinia* was inoculated in different ways: plates 1 through 4 containing transformed *Erwinia* with different aiiA colonies, plate 5 containing *Erwinia* wild and plate 6 containing an apple without inoculation used as a negative control. Each plate divided in two, with each part inoculated by different techniques.

September/07/2017

Protocol: Miniprep and overnight cultures

Plasmid extraction by alkaline lysis of P+RBS+epsE and *Erwinia Amylovora*.

LB broth (TetR) tubes inoculated with RBS+yhjH and RBS+epsE (CAM).

September/08/2017

Protocol: Partial digestion and miniprep.

Standard mix in standard conditions to perform partial digestion of aiiA *Erwinia*, *Erwinia* wild and P+RBS+epsE.

Plasmid extraction of RBS+yhjH and RBS+epsE.

September/09/2017

Protocol: Nanodrop, Miniprep, partial digestion, overnight cultures and 12-hour digestion

Plasmid extraction of RBS+yhjH and RBS+epsE.

Partial digestion of RBS+yhjH and RBS+epsE.

LB broth tubes inoculated with RBS.

12-hour digestion of epsE+RBS.

Nanodrop

Part	First Measurement	Second Measurement
P+RBS+epsE	43.4/1.82	44.5/1.81
P+RBS+yhjH	64.4/1.81	68.4/1.82

Measurements: Concentration / 260/280

September/10/2017

Protocol: Miniprep, partial digestion, transformation and Ligation.

Plasmid extraction of (Two colonies)

Partial digestion of P+RBS+epsE.

Transformation of *Erwinia Amylovora* into electrocompetent aiiA biobrick and yhjH biobrick.

12-hour ligation of of RBS+epsE+T in TetR vector and RBS + yhjH in TetR vector.

September/11/2017

Protocol: Partial digestion, digestion, chemically competent cells and transformation and overnight cultures

Partial digestion with standard mix of RBS+epsE and RBS+yhjH.

12-hour digestion of P+RBS+epsE (Mix: 7.3 μ L DNA, 1.5 μ L XbaI enzyme, 1.5 μ L PstI enzyme, buffer 5 μ L and 34.7 μ L nuclease free water).

LB agar plate(TetR) inoculated with ligated transformations.

September/12/2017

Protocol: 12-hour ligation and overnight cultures

12-hour ligation of epsE biobrick (P+RBS+epsE+T)

LB broth tubes inoculated with transformations, RBS+epsE+T and RBS+yhjH inoculated in LB agar plates.

September/13/2017:

Protocol: Overnight cultures, miniprep, transformation

Erwinia Amylovora and BL21(DE3) inoculated in LB broth tubes supplemented with sucrose.

Plasmid extraction by alkaline lysis of RBS+epsE+T and RBS+yhjH.

September/14/2017

Protocol: Solutions preparation for colorimetry tube test, digestion and overnight cultures

Solutions preparation for the colorimetry microtiter-plate test.

12-hour digestion of RBS.

LB broth tubes inoculated with BL21(DE3) and transformed epsE biobrick.

September/15/2017

Protocol: Miniprep, digestion, overnight cultures

Plasmid extraction by alkaline lysis of BL21(DE3) and *Erwinia*

12-hour digestion of RBS+yhjH.

LB broth tubes inoculated with RBS.

SDS polyacrylamide gel electrophoresis with SDS page of *aiiA Erwinia Amylovora*.

September/16/2017

Protocol: Miniprep, partial digestion, digestion and Ligation

Plasmid extraction of RBS and RBS+yhjH (two colonies)

Plasmid extraction by alkaline lysis of epsE biobrick (5 colonies) and RBS+epsE+T (3 colonies).

Partial digestion of RBS and RBS+yhjH.

12-hour digestion and ligation of RBS+ yhjH in TetR vector.

September/17/2017

Protocol: Overnight cultures and Ligation

LB broth tubes inoculated with plasmid extractions.


12-hour ligation of RBS+yhjH in TetR vector, RBS+epsE in CAM vector and RBS+epsE+T in TetR vector.

September/18/2017

Protocol: Miniprep and Partial digestion

Plasmid extraction of epsE biobrick.

Partial digestion of epsE biobrick.



September/19/2017

Protocol: Transformation

Transformation of epsE biobrick into BL21(DE3).

September/20/2017

Protocol: Electroporation

Electroporation of *Erwinia Amylovora* with epsE biobrick.

September/21/2017

Protocol: Overnight cultures

LB broth tubes inoculated with epsE biobrick transformations.

September/22/2017

Protocol: Miniprep

Plasmid extraction by alkaline lysis of RBS+yhjH (two colonies) and epsE biobrick.

September/23/2017

Protocol: Overnight cultures

LB agar plates supplemented with sucrose inoculated with *Erwinia Amylovora* and LB broth tubes inoculated with RBS+yhjH.

September/24/2017

Protocol: Miniprep and electrophoresis

Plasmid extraction by alkaline lysis of RBS+yhjH.



October

October/01/2017

Protocol: Overnight cultures

LB broth tubes supplemented with sucrose inoculated with *Erwinia Amylovora* and LB broth tubes inoculated with *aiiA* and *yhjH* biobricks.

October/02/2017

Protocol: Overnight cultures and miniprep

LB broth tubes inoculated with TOP10

Plasmid extraction by alkaline lysis of the biobricks (15 colonies).

October/03/2017

Protocol: Partial digestion

Partial digestion of the 15 plasmid extractions (biobricks and *Erwinia Amylovora*).

October/04/2017

Protocol: Transformation

Transformation by electroporation of *epsE* into *Erwinia Amylovora*.

October/05/2017

Protocol: Overnight cultures

Inoculation of 24-hour Elisa plates, BL21(DE3) with *epsE* biobrick transformation.

October/06/2017

Protocol: Miniprep, partial digestion and colorimetry

Plasmid extraction by alkaline lysis of *epsE* biobrick.

Partial digestion of *epsE* biobrick.

Biofilm quantification using 24-hour tubes with safranin and microtiter plate test.

October/09/2017

Protocol: DNA resuspension

Resuspension of IDT aiiA and IDT epsE.

October/10/2017

Protocol: Overnight cultures

LB agar plates inoculated with BL21(DE3). Two plates inoculated with transformations of IDT aiiA (AMP) and IDT epsE (AMP). Another plate inoculated with epsE (CAM).

October/11/2017

Protocol: Digestion and overnight cultures

Digestion of IDT aiiA and IDT epsE.

Inoculation of electroporations and BL21(DE3).


October/12/2017

Protocol: Miniprep, partial digestion, overnight cultures and electrophoresis

Plasmid extraction of epsE biobrick in *Erwinia Amylovora*.

Partial digestion of epsE biobrick in *Erwinia Amylovora*.

Inoculation of *Erwinia Amylovora* and transformations of yhjH, aiiA and epsE.



October/13/2017

Protocol: Chemically competent cells, transformation and characterization

BL21(DE3) Chemically competent cells.

Transformation of IDT aiiA and IDT epsE parts.

Inoculation on the apple leaves for characterization trials.

October/14/2017

Protocol: Transformation, digestion and overnight cultures

Transformation of IDT sequences (aiiA and epsE).

12-hour digestion of IDT sequences (aiiA and epsE).

LB agar plates inoculated with RFP and BL21(DE3). Inoculation of *Erwinia Amylovora* for the collaboration with TEC CEM team.

October/15/2017

Protocol: Miniprep, digestion, overnight cultures and electroporation

Plasmid extraction of RFP.

12-hour digestion of RFP.

Inoculation of BL21(DE3) and *Erwinia Amylovora*.

Electroporation of IDT aiiA.

October/16/2017

Protocol: Chemically competent cells, overnight cultures and transformation

BL21(DE3) Chemically competent cells.

LB agar plates supplemented with sucrose inoculated with *Erwinia Amylovora*.

Transformation of IDT dilutions sets of 48 hours *Erwinia* and 48 hours of BL21(DE3).

October/17/2017

Protocol: Overnight cultures and colorimetry solutions preparation

Inoculation of aiiA and epsE transformations and aiiA *Erwinia Amylovora*.

Solutions preparation for the colorimetry microtiter-plate test for quantification of biofilm formation.


October/19/2017

Protocol: Overnight cultures, digestion and electrophoresis

Inoculation of yhjH, epsE and pseudomonas as a positive control in *Erwinia Amylovora*.

Separation of Elisa plates.

12-hour digestion of IDT aiiA and IDT epsE.



October/20/2017

Protocol: Miniprep

Plasmid extraction of IDT parts *Erwinia Amylovora*.

October/21/2017

Protocol: Digestion, transformation and characterization

12-hour digestion of RFP.

Transformation by electroporation of epsE IDT.

Inoculation in apple slices as a part of the characterization trials.

•ERWINIONS•