

# GREEN YOUR LAB



## 1. Freezer


Defrost and clean your freezers regularly!

Become a sample minimalist – store only things you still need.




## 2. Fume Hood

Close the sash of your fume hood each time you stop working there and remind your lab mates to do the same!



## 3. Autoclave

Do not run the autoclave only half filled but do not overload it either.



## 4. Electrical devices

PCR Cylers, Centrifuges, Shakers, ... etc.

Remember to turn them off every time you finish your task!



## 5. Consolidate purchases

Reduce packaging waste and delivery trips by combining your orders with other labs or by planning ahead. Additionally, if sterility is not an issue, try to buy items that are packed in bulk.



## 6. Dry orders

Order your oligos dry.

It saves shipping weight, package material and is more stable.



## 7. Paper and printing products

Change the default settings of your printer to double-sided and use recycled paper. Better: avoid printing if possible.



## 8. Computers

Shut down when not in use or share a computer with lab mates



→ For more information read our guideline / Follow us

- [Guideline](#) →
- [Instagram: igem\\_goes\\_green](#)

