

7/6/17


Thursday, July 6, 2017 10:55 AM

Run Gel Electrophoresis to confirm our PCR amplification products

Gel

- 50 mL TAE buffer
- .500 g Agarose
- 5 uL Midori Green
- Set for 45 minutes

Loaded



PCR products tested from 7/5/17:

Lane 1: 2-log DNA ladder

Lane 3: JOE PCR Product

Lane 5: YCP PCR Product

JOE is much shorter than its expected 2kb amplification. YCP appears to be on target at 700bp