

7/5/17

Wednesday, July 5, 2017 1:44 PM

Who's in lab today: FULL TEAM !
UIUC Sequencing requires 10uL in 1.7mL tubes

LC1539	T7 Pro
LC1539	T7 Term
LC1853	T7 Pro
LC1853	T7 Term

Made Stock Solution of Primers: JoeF, JOER, YCPF, YCPR

Added the appropriate amount of sterile water to the primer tubes to make it 100uM (the recommended volume is on the sheets that come with the primers)

Made working solution of primers:

Mixed 10uL of stock solution of primers with 90uL of sterile water to create the working solution of primers