

Technical Data

Jensen's Medium M710

Jensens Medium is recommended for detection and cultivation of nitrogen fixing bacteria.

Composition**

Ingredients	Gms / Litre
Sucrose	20.000
Dipotassium phosphate	1.000
Magnesium sulphate	0.500
Sodium chloride	0.500
Ferrous sulphate	0.100
Sodium molybdate	0.005
Calcium carbonate	2.000
Agar	15.000

^{**}Formula adjusted, standardized to suit performance parameters

Directions

Suspend 39.1 grams in 1000 ml distilled water. Heat just to boiling . Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

Note: Due to presence of calcium carbonate, the medium forms opalescent solution with white precipitate .

Principle And Interpretation

Nitrogen-fixing organisms are free-living bacteria, which grow well on a nitrogen-free medium. These bacteria utilize atmospheric nitrogen gas for their cell protein synthesis. This cell protein is then mineralized in soil after the death of the cells thereby contributing towards the nitrogen availability of the crop plants (1). Nitrogen fixing bacteria enter into symbiosis only with leguminous plants, by infecting their roots and forming nodules on them. Jensens Medium is formulated according to Jensen and is recommended for detection and cultivation of nitrogen fixing bacteria (2).

Sucrose acts as the energy source. Sodium molybdate in the media increases the fixation of nitrogen (3). Sodium chloride maintains osmotic equilibrium of the media. Calcium stimulates nodulation when present as chloride or sulphate.

Quality Control

Appearance

White to cream homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Cream coloured, slightly opalescent gel with a slight precipitate forms in Petri plates.

Cultural Response

M710: Cultural characteristics observed after incubation at 30°C for 8 days.

Organism	Growth
Rhizobium leguminosarum	luxuriant
ATCC 10004	
Rhizobium meliloti ATCC	luxuriant
9930	
Rhizobium oryzae ATCC	luxuriant
9363	

Storage and Shelf Life

Store below 30°C in tightly closed container and use freshly prepared medium. Use before expiry period on the label.

HiMedia Laboratories Technical Data

Reference

1.Subba Rao N. S., 1977, In: Soil Microorganisms and Plant Growth, Oxford and IBH Publishing Co., New Delhi, Pages 254-255.

2.Jensen. H. L., 1942, Pro Line Soc. N.S.W., 57,205-212.

3.Ranganayaki S., Mohan C., Ally Z., 1981; 21 (8): 607-10.

Revision: 2 / 2015

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.